

R e m o t e – C o n t r o l C o m m a n d s f o r R a d E y e

Version check

[illegible]

*) Category K: editorial correction
 V: explanatory improvement
 S: substantial change

An explanation is required at least for category S.

Contents :

1. Connection to a PC	1-1
2. Data transmission method.....	2-1
3. General.....	3-1
3.1 HEX-Format	3-1
3.2 Date and time as a decimal value.....	3-1
4. Common remote-control commands	4-1
4.1 History	4-1
4.2 Event log	4-1
4.3 Date and time	4-1
4.4 EEPROM	4-2
4.5 Configuration	4-2
4.6 Serial interface	4-3
4.7 Calibration	4-3
4.8 RadEye Type.....	4-4
4.9 Automatic sending	4-4
4.10 Device description.....	4-4
4.11 Measurement values.....	4-4
4.11.1 High voltage	4-5
4.12 Battery test	4-5
4.12.1 During start up	4-5
4.12.2 During operation.....	4-6
4.12.3 Power down	4-6
4.12.4 Bluetooth configuration.....	4-7

5. RadEye specific remote-control commands5-1

5.1	RadEye PRD / PRD-ER (<V3.00).....	5-1
5.1.1	Limit values	5-1
5.1.2	Measurement values	5-3
5.1.3	Configuration flags.....	5-4
5.1.4	High voltage	5-7
5.1.5	Dead time correction	5-7
5.1.6	History output	5-7
5.1.7	Event log.....	5-9
5.1.8	Automatic sending.....	5-9
5.1.9	Status information.....	5-10
5.2	RadEye PRD-S / PRD-ER-S (<V3.00).....	5-12
5.2.1	Limit values	5-12
5.2.2	Measurement values	5-14
5.2.3	Configuration flags.....	5-15
5.2.4	High voltage	5-20
5.2.5	Dead time correction	5-20
5.2.6	History output	5-20
5.2.7	Event log.....	5-22
5.2.8	Automatic sending.....	5-22
5.2.9	Status information.....	5-23
5.3	RadEye PRD / PRD-ER / PRD-S / PRD-ER-S (>V3.00).....	5-25
5.3.1	Limit values	5-25
5.3.2	Measurement values	5-27
5.3.3	Configuration flags.....	5-28
5.3.4	High voltage	5-33
5.3.5	Dead time correction	5-33
5.3.6	History output	5-34
5.3.7	Event log.....	5-35
5.3.8	Nuclide calibration data (from V3.04).....	5-35
5.3.9	Automatic sending.....	5-37
5.3.10	Status information.....	5-38
5.4	RadEye G/G-10.....	5-39
5.4.1	Limit values	5-39
5.4.2	Measurement values	5-39
5.4.3	Configuration flags.....	5-40
5.4.4	High voltage	5-42
5.4.5	Dead time correction	5-42
5.4.6	History output	5-43
5.4.7	Event log.....	5-44
5.4.8	Automatic sending.....	5-44
5.4.9	Status information.....	5-45

5.5	RadEye N	5-47
5.5.1	Limit values	5-47
5.5.2	Measurement values	5-48
5.5.3	Configuration flags	5-49
5.5.4	High voltage	5-54
5.5.5	Dead time correction	5-54
5.5.6	History output	5-54
5.5.7	Event log	5-56
5.5.8	Automatic sending	5-57
5.5.9	Status information	5-58
5.6	RadEye AB100 (V1.52)	5-59
5.6.1	Limit values	5-59
5.6.2	Measurement values	5-60
5.6.3	Configuration flags	5-61
5.6.4	High voltage	5-65
5.6.5	Dead time correction	5-65
5.6.6	History output	5-65
5.6.7	Event log	5-68
5.6.8	Automatic sending	5-69
5.6.9	Status information	5-70
5.7	RadEye B20/G20/G20-10	5-71
5.7.1	Limit values	5-71
5.7.2	Measurement values	5-72
5.7.3	Configuration flags	5-73
5.7.4	High voltage	5-77
5.7.5	Dead time correction	5-77
5.7.6	History output	5-78
5.7.7	Nuclide table	5-80
5.7.8	Event log	5-82
5.7.9	Automatic sending	5-83
5.7.10	Status information	5-84
5.8	RadEye DW	5-85
5.8.1	Limit values	5-85
5.8.2	Measurement values	5-86
5.8.3	Configuration flags	5-87
5.8.4	High voltage	5-89
5.8.5	Dead time correction	5-90
5.8.6	History output	5-90
5.8.7	Event log	5-91
5.8.8	Automatic sending	5-91
5.8.9	Status information	5-93
5.8.10	Miscellaneous	5-93

5.9	RadEye DLW	5-95
5.9.1	Limit values	5-95
5.9.2	Measurement values	5-95
5.9.3	Configuration flags	5-96
5.9.4	High voltage	5-98
5.9.5	Dead time correction	5-99
5.9.6	History output	5-99
5.9.7	Event log	5-100
5.9.8	Automatic sending	5-100
5.9.9	Status information	5-101
5.10	RadEye AB100 (from V1.67)	5-103
5.10.1	Limit values	5-103
5.10.2	Measurement values	5-105
5.10.3	Scaler remote control (from V1.73)	5-105
5.10.4	Configuration flags	5-106
5.10.5	High voltage	5-110
5.10.6	History output	5-111
5.10.7	Nuclide table	5-113
5.10.8	Probe data	5-114
5.10.9	Event log	5-116
5.10.10	Automatic sending	5-117
5.10.11	Status information	5-118
5.10.12	Miscellaneous	5-118
5.11	RadEye GX	5-120
5.11.1	Limit values	5-120
5.11.2	Measurement values	5-122
5.11.3	Scaler remote control (from V1.66)	5-122
5.11.4	Configuration flags	5-122
5.11.5	High voltage	5-128
5.11.6	Dead time correction	5-128
5.11.7	History output	5-128
5.11.8	Nuclide table	5-130
5.11.9	Probe data	5-131
5.11.10	Event log	5-133
5.11.11	Automatic sending	5-134
5.11.12	Status information	5-135
5.12	RadEye SX	5-136
5.12.1	Limit values	5-136
5.12.2	Measurement values	5-138
5.12.3	Scaler remote control (from V1.73)	5-139
5.12.4	Configuration flags	5-139
5.12.5	High voltage	5-143
5.12.6	History output	5-143
5.12.7	Nuclide table	5-146

5.12.8	Probe data	5-147
5.12.9	Event log.....	5-149
5.12.10	Automatic sending	5-150
5.12.11	Status information.....	5-152
5.12.12	Miscellaneous	5-153
5.13	RadEye PX.....	5-154
5.13.1	Limit values	5-154
5.13.2	Measurement values	5-156
5.13.3	Configuration flags.....	5-157
5.13.4	High voltage	5-161
5.13.5	History output.....	5-162
5.13.6	Nuclide table.....	5-164
5.13.7	Probe data	5-165
5.13.8	Event log.....	5-167
5.13.9	Automatic sending.....	5-168
5.13.10	Status information.....	5-169
5.13.11	Miscellaneous	5-169
5.14	RadEye HEC	5-171
5.14.1	Limit values	5-171
5.14.2	Measurement values	5-172
5.14.3	Scaler remote control (from V1.69).....	5-172
5.14.4	Configuration flags.....	5-173
5.14.5	High voltage	5-177
5.14.6	Remote control (from V1.68)	5-177
5.14.7	History output.....	5-178
5.14.8	Nuclide table.....	5-180
5.14.9	Parameter data	5-181
5.14.10	Event log	5-183
5.14.11	Automatic sending	5-184
5.14.12	Status information.....	5-185
5.14.13	Checksource.....	5-186
5.15	RadEye GFex/GF-10ex/Gex/G-10ex	5-187
5.15.1	Limit values	5-187
5.15.2	Measurement values	5-187
5.15.3	Configuration flags.....	5-188
5.15.4	High voltage	5-191
5.15.5	Dead time correction	5-191
5.15.6	History output.....	5-191
5.15.7	Event log.....	5-193
5.15.8	Automatic sending.....	5-193
5.15.9	Status information.....	5-194

5.16	RadEye G/G-10/B20/G20/G20-10/GF/GF-10 (from V2.00)	5-196
5.16.1	Limit values	5-196
5.16.2	Measurement values	5-197
5.16.3	Configuration flags	5-199
5.16.4	High voltage	5-204
5.16.5	Dead time correction	5-204
5.16.6	History output	5-204
5.16.7	Nuclide table	5-206
5.16.8	Event log	5-208
5.16.9	Automatic sending	5-209
5.16.10	Status information	5-211
5.17	RadEye G/G-10/B20/G20/G20-10/GF/GF-10 (from V3.00)	5-211
5.17.1	Limit values	5-212
5.17.2	Measurement values	5-213
5.17.3	Configuration flags	5-214
5.17.4	High voltage	5-219
5.17.5	Dead time correction	5-219
5.17.6	History output	5-219
5.17.7	Nuclide table	5-221
5.17.8	Event log	5-223
5.17.9	Automatic sending	5-224
5.17.10	Status information	5-225
5.18	RadEye GN/GN+	5-226
5.18.1	Limit values	5-226
5.18.2	Measurement values	5-227
5.18.3	Configuration flags	5-228
5.18.4	High voltage	5-232
5.18.5	History output	5-232
5.18.6	Nuclide table	5-234
5.18.7	Event log	5-234
5.18.8	Automatic sending	5-236
5.18.9	Status information	5-237
5.18.10	Miscellaneous	5-237
5.19	RadEye PRD-CD	5-238
5.19.1	Limit values	5-238
5.19.2	Measurement values	5-241
5.19.3	Configuration flags	5-242
5.19.4	Spill-over	5-246
5.19.5	High voltage	5-247
5.19.6	Dead time correction	5-247
5.19.7	History output	5-248
5.19.8	Event log	5-249

5.19.9	Calibration log	5-250
5.19.10	Automatic sending	5-250
5.19.11	Status information.....	5-251

1.Connection to a PC

The RadEye features an infrared interface allowing for direct connection to a PC with RS232 interface when a special adapter cable (42540/29 or 42540/26) is used.

In this manner, transfer of stored measured values and the configuration of the unit are possible.

For these tasks special Windows PC programs are available.

2.Data transmission method

Data transmission is performed in ASCII code using the following parameters:

- 9600 Baud
- Start bit
- 7 data bits
- Parity bit: even
- 2 stop bits

To allow for remote control of the RadEye unit via a PC, the following data transmission procedure has to be observed:

- PC sends <@> (40) without <CR> (0D)/<LF> (0A). to the RadEye
- RadEye replies with the character ">".
- PC waits 500µs and send a remote-control command where the last character must be a line feed character <LF> (0A). Before the line feed character, there may be a carriage return character <CR> (0D):
- The RadEye answers with „#“ as positive acknowledgement, or with „?“ in case of an unknown command
- If an output was prompted, the output data will follow now.
- The RadEye terminates each transmission with <CR> <LF>.

Following control lines are needed for the power supply of the interface adapter cable with RS232 connector (42506/29) and must be set by the PC:

- RTS (Ready To Send) must always have at positive voltage level! No handshake!
- DTR (Data Terminal Ready) must always have a negative voltage level!

3. General

3.1 HEX-Format

The flags are coded bit-by-bit and are put out hexadecimally.

E.g.: 84

means:

Hexadezimally	8	4
Binary:	1 0 0 0	0 1 0 0
Bit number	7 6 5 4	3 2 1 0

E.g.: 84DA

means:

Hexadezimally	8	4	D	A
Binary	1 0 0 0	0 1 0 0	1 1 0 1	1 0 1 0
Bit number	15 14 13 12	11 10 9 8	7 6 5 4	3 2 1 0

E.g.: 84DA1F50

means:

Hexadezimally	8	4	D	A	1	F	5	0
Binary	1 0 0 0	0 1 0 0	1 1 0 1	1 0 1 0	0 0 1 1	1 1 1 0	1 0 1 0	0 0 0 0
Bit number	31 30 29 28	27 26 25 24	23 22 21 20	19 18 17 16	15 14 13 12	11 10 9 8	7 6 5 4	3 2 1 0

3.2 Date and time as a decimal value

Date and time for history and eventlog comes as a decimal value:

value converted to binary:

010101	0101	01010	10101	010101	010101	
-----	----	-----	-----	-----	-----	
						6 Bit seconds
						6 Bit minutes
						5 Bit hour
						5 Bit day
						4 Bit month
						6 Bit year

4.Common remote-control commands

4.1 History

TR	Read history cycle time Response: Number in seconds e.g. 300 means 300s
TW <i>Number</i>	Set history cycle time <i>Number</i> : Time in second from 0 to 43200
HI	Initialize history readout Response: --
+	Read next history data set Response: see RadEye specific remote control command
-	Read last history data set Response: see RadEye specific remote control command
ph	Clear history Response: --

4.2 Event log

EI	Initialize readout of event log Response: --
E+	Read next data set Response: see RadEye specific remote control command
E-	Read last data set Response: see RadEye specific remote control command
EC	Clear event log Response: --

4.3 Date and time

ZR	Reading date and time Response in the format: year, month, day, hour, minute, second (JJMMTTHHMMSS).
----	--

e.g.: 100927172845

ZWNumber Setting date and time
Format: as described above

4.4 EEPROM

EW Store configuration to EEPROM
Response:--

4.5 Configuration

mR Read menu configuration
Response: hexadecimal value
Meaning: see RadEye specific commands

mWHex Write menu configuration
Meaning: see RadEye specific commands

fR Read configuration flag 1
Response: hexadecimal value
Meaning: see RadEye specific commands

fWHex Write configuration flag 1
Meaning: see RadEye specific commands

kR Read configuration flag 2
Response: hexadecimal value
Meaning: see RadEye specific commands

kWHex Write configuration flag 2
Meaning: see RadEye specific commands

KR Read configuration flag 3
Response: hexadecimal value
Meaning: see RadEye specific commands

KWHex Write configuration flag 3
Meaning: see RadEye specific commands

sR Read menu language
Response: Number
0: English
1: German
2: French

<i>sWNumber</i>	Write menu language Number: see above
ART	Read Timeout of alarm latching Response: Number from 0..255 in seconds
AWT <i>Number</i>	Setting Timeout of alarm latching <i>Number</i> : Number from 0..255 in seconds

4.6 Serial interface

ARS	Read serial time out Response: Number in seconds. Default 2s
AWS	Write serial time out Description: see above
t	Reading and resetting transfer error counter. <u>Response:</u> 0 – no error n – error number

4.7 Calibration

WR	Read calibration date Response: Date as JJMMDD
\$	Read calibration factor Response: see RadEye specific command
#R	Read device serial number Response: Number from 0 to 65535

4.8 RadEye Type

Vx Read RadEye type
Response: RadEye type, firmware version and firmware checksum
e.g. RadEye PRD V1.52 AB48

4.9 Automatic sending

X0 Deactivate cyclic sending of the measurement value.

X1 Activate cyclic sending of the measurement value. The dose rate measurement value is sent every second. Format: see RadEye specific remote control command

4.10 Device description

DR Reading a text stored in the device. This text cannot be displayed
Response: Text, up to 200 characters.

DW*Number*Text Writing a text
Number: 0..9
Text: 20 characters

dR*Number* Reading text info. This text is displayed via menu item “Text info”
Number: Representing line number. Range 0...3. 0 means bottom line, 3 means top line
Response: Text

dW*Number*Text Writing text info. This text is displayed via menu item “Text info”
Each line contains up to 16 characters
Number see above

4.11 Measurement values

Z Read raw count rates with dead time correction
Format: see RadEye specific command

z Read filtered count rate with dead time correction
Format: see RadEye specific command

R Read measured dose rate
see RadEye specific command


D	Read accumulated dose see RadEye specific command
clr	Clear accumulated dose and reset overload flag
Ux	Read battery voltage Response: Number in 0.1V units
F	Read status Response: see RadEye specific command
tR	Read temperature Response: Temperature in °C

4.11.1 High voltage


HR	Reading nominal value of high voltage Response: see RadEye specific command
----	--

4.12 Battery test


4.12.1 During start up



4.12.2 During operation



4.12.3 Power down



4.12.4 Bluetooth configuration

XRa Get bluetooth remote address. Return: string with up to 11 characters
XRp Get bluetooth pin. Return: string with up to 16 characters
XRf Get bluetooth parameter. Return Hex-value (see below)

XWa *String* Set bluetooth remote address. *String* with up to 11 characters
XWp *String* Get bluetooth pin. *String* with up to 16 characters
XWf *hex* Get bluetooth parameter. *hex*-value see below

Bit number			
0	Mode	0: Remote	1: Master
1	Btcom cover	0: off	1: on
2	BTcom cover off at low battery	0: no	1: yes
3	Use bluetooth pin	0: no	1: yes
4	Use this settings	0: no	1: yes
5	Firewall	0: no	1: yes
6	LED	0: no	1: yes
7	Bluetooth active	0: no	1: yes

BT0 Bluetooth module not available (low Batt)
BT1 Bluetooth adapter is on, but not connected
BT2 Bluetooth adapter is connected
BT3 Bluetooth adapter was connected, now disconnected
BTE Bluetooth module error

5.RadEye specific remote-control commands

5.1 RadEye PRD / PRD-ER (<V3.00)

Used firmware version: 1.52

5.1.1 Limit values

AR0	Reading computed sigma alarm threshold 1 . Response: number in cps.
AR1	Reading the threshold 1 for level alarm. Response: number from 2..9. Number 2 means 32 cps. Each subsequent level has twice as much cps as the preceding level.
AR2	Reading the threshold 2 for level alarm. Response: see command AR1
AR3	Reading the threshold 1 for count rate alarm. Response: number in cps.
AR4	Reading the threshold 2 for count rate alarm. Response: number in cps.
AR5	Reading the threshold 1 for dose rate alarm. Response: number in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ unit. Depending on used measuring unit e.g. 1 2 3 means 1 2 3 $\mu\text{R/h}$ resp. 1.23 $\mu\text{Sv/h}$
AR6	Reading the threshold 2 for dose rate alarm. Response: see command AR5
AR7	Reading the threshold 1 for the dose alarm Response: number in μSv , $100\mu\text{R}$ or $100\mu\text{rem}$ units
AR8	Reading the threshold 2 for the dose alarm see command AR7
AR9	Reading sigma value Response: Number from 2..9

ARM	Reading minimum count rate for sigma alarm Response: Number in cps
ARN	NBR alarm threshold low energy Response: Number in 0.01% units
ARH	NBR alarm threshold high energy Response: Number in 0.01% units
SR3	Reading minimum count rate for NBR alarm Response: Number in cps
AW1Number	Setting the threshold 1 for level alarm. <i>Number</i> = 2...9. Number 2 means 32 cps. Each subsequent level has twice as much cps as the preceding level
AW2Number	Setting the threshold 2 for level alarm. <i>Number</i> = 2...9. Number 2 means 32 cps. Each subsequent level has twice as much cps as the preceding level
AW3Number	Setting the threshold 1 for count rate alarm. <i>Number</i> = Number in cps.
AW4Number	Setting the threshold 2 for count rate alarm. <i>Number</i> = see command AW3.
AW5Number	Setting the threshold 1 for dose rate alarm. <i>Number</i> = in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units. e. g. AW3123 means $123\mu\text{R/h}$.
AW6Number	Setting the threshold 2 for dose rate alarm. see command AW3
AW7Number	Setting the threshold 1 for dose alarm. <i>Number</i> = in $100\mu\text{R}$, $100\mu\text{rem}$ or $1\mu\text{Sv}$ units. e. g. AW7123 means 12.3mR .
AW8Number	Setting the threshold 2 for dose alarm. <i>Number</i> = see command AW7
AW9Number	Setting the sigma value. <i>Number</i> = form 2 to 9
AWMNumber	Setting minimum count rate for sigma alarm Number= value from 0 to 255

<i>AWNNumber</i>	Setting NBR alarm threshold level low energy Number=value in 0.01% units. From 1.00% to 2.55%
<i>AWHNumber</i>	Setting NBR alarm threshold level high energy Number=value in 0.01% units. From 0% to 1.00%
<i>SW3 Number</i>	Setting minimum count rate for NBR alarm Number= value from 0 to 255

5.1.2 Measurement values

<i>Z</i>	Read raw count rates with dead time correction Response: <ul style="list-style-type: none"> - Counter 1 in cps - Counter 2 in cps - Counter 3 in cps - HV power in cps
<i>z</i>	Read filtered count rate 1 Response: Number in cps

5.1.3 Configuration flags

5.1.3.1 Configuration flags 1 with kR / kW

Bit number	
0	0x01: Measuring unit: Level
1	0x02: Measuring unit: Count rate
2	0x04: Measuring unit: Dose rate
3	Alarming Sound 0: off 1: on
4	Alarming LED 0: off 1: on
5	Alarming Vibration 0: off 1: on
6	NBR 0: off 1: on
7	Single Pulse 0: off 1: on

5.1.3.2 Configuration flags 2 with fR / fW

Bit number	
0	10: Dose rate Sievert ($H^*(10)$)
1	01: Dose rate Roentgen (PRD: H_x , PRD-ER: $H^*(10)$) 11: Dose rate rem ($H^*(10)$)
2	not used, write "0"
3	Alarm threshold read-only 0: off 1: on
4	Flag for overload (readonly)
5	Temp.display 0: off 1: on
6	Temperature unit 0: °C 1: °F
7	In cps-mode display of dose rate 0: no 1: yes

5.1.3.3 *Configuration flags 3 with KR / KW*

Bit number	
0	Key tone 0:off 1:on
1	Autosend 0:off 1:on
2	Finder 0:off 1:on
3	Mean-Max 0:off 1:on
4	not used, write "0"
5	display of dose 0:off 1:on
6	not used, write "0"
7	not used, write "0"

5.1.3.4 *Configuration flags 4 with jR / jW*

Bit number	
0	Safety-Alarm 0:off 1:on
1	Battery type 0:Alkaline 1:NiMh
2	Display rotation 0:No 1:yes
3	not used, write "0"
4	not used, write "0"
5	not used, write "0"
6	not used, write "0"
7	not used, write "0"
8	not used, write "0"
9	not used, write "0"
10	not used, write "0"
11	not used, write "0"
12	not used, write "0"
13	not used, write "0"
14	not used, write "0"
15	not used, write "0"

5.1.3.5 *Menu configuration*

Bit number			
0	Switch off	0:hidden	1:visible
1	Sound	0: hidden	1:visible
2	LED	0: hidden	1:visible
3	Vibrator	0: hidden	1:visible
4	Level	0: hidden	1:visible
5	Count Rate	0: hidden	1:visible
6	Dose Rate	0: hidden	1:visible
7	Alarm Counter	0: hidden	1:visible
8	Alarm Dose Rate	0: hidden	1:visible
9	Alarm Dose	0: hidden	1:visible
10	Alarm Level	0: hidden	1:visible
11	Autosend	0: hidden	1:visible
12	NBR	0: hidden	1:visible
13	Clear Dose	0: hidden	1:visible
14	Finder	0: hidden	1:visible
15	Single Pulse	0: hidden	1:visible
16	Backlight	0: hidden	1:visible
17	Show Alarm	0: hidden	1:visible
18	Settings	0: hidden	1:visible
19	Text Info	0: hidden	1:visible
20	not used, write "0"		
21	not used, write "0"		
22	not used, write "0"		
23	not used, write "0"		
24	not used, write "0"		
25	not used, write "0"		
26	not used, write "0"		
27	not used, write "0"		
28	not used, write "0"		
29	not used, write "0"		
30	not used, write "0"		
31	not used, write "0"		

5.1.4 High voltage

HR	Reading high voltage bit value Response: Bit value 0...255.
hR	Reading high voltage correction bit date and time of last successful Luthetium check Response: Number 0...255 with offset 128 and date and time as YYMMDDhhmm
hWNumber	Setting of high voltage correction bit Number from 0...255

5.1.5 Dead time correction

x	Read dead time Response: dead time in ns for - Counter 1 (Rate 1) - Counter 2 (Rate 2) - Counter 3 (Rate 3)
---	---

5.1.6 History output

5.1.6.1 History readout

History status for 1. readout and change of history cycle time

1073741944 520534610 0 0

```
-----
|           |           |
|           |           | \-- Always 0
|           | \----- Date and Time as a decimal value
| \----- History status as a decimal value
```

following readout:

15 11 166 105 +19

```
-----
| | | | |
| | | | | \--- Temperature in °C
| | | | | \----- Rate mean in µR/h, µrem/h or 0.01µSv/h units
| | | | | \----- Rate max in µR/h, µrem/h or 0.01µSv/h units
| | | | | \----- Counter 1 mean
| | | | | \----- Counter 1 max
```

End of History:

End

5.1.6.2 *History status*

decimal value converted in HEX:

0x40000078

```
|      ----  
|      |  
|      \-- History cycle time (0x78 -> 120s)  
|  
\----- 4: identifier for status information  
          6: identifier for status information +  
            last power on without power off
```


5.1.7 Event log

6656 520549251

```

-----
|          |
|          | \----- Date and time as a decimal value (see 3.2)
|          | \----- Event log as a decimal value
\-----

```

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	0x01 Level display
9	0x02: Display count rate
10	0x04: Display dose rate
11	Sound 0: off 1:on
12	LED 0: off 1:on
13	Vibration alarm 0: off 1:on
14	1: Dose cleared
15	1: Alarm threshold changed
16	1: Count rate, dose rate or level-alarm
17	1: Dose alarm
18	1: Safety alarm
19	Not used
20	Value > alarm threshold 1 (depending on count rate, dose rate or level display)
21	Value > alarm threshold 2 (depending on count rate, dose rate or level display)
22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24	Low energy alarm
25	High energy alarm
26	Power off
27	Power on
28	NBR-alarm
29	Not used
30	Not used
31	Not used

5.1.8 Automatic sending

<STX>7 2 0 0 14 FH41PR 123 AB<ETX><CR><LF>

16	1: Count rate, dose rate or level-alarm
17	1: Dose alarm
18	1: Safety alarm
19	Not used
20	Value > alarm threshold 1 (depending on count rate, dose rate or level display)
21	Value > alarm threshold 2 (depending on count rate, dose rate or level display)
22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24	Low energy alarm
25	High energy alarm
26	Not used
27	Not used
28	Not used
29	Not used
30	Not used
31	Fixed to "0"

5.2 RadEye PRD-S / PRD-ER-S (<V3.00)

Used firmware version: V3.05

5.2.1 Limit values

AR0	Reading computed sigma alarm threshold 1. Response: number in cps.
AR3	Reading the threshold 1 for count rate alarm. Response: number in cps.
AR4	Reading the threshold 2 for count rate alarm. Response: number in cps.
AR5	Reading the threshold 1 for dose rate alarm. Response: number in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ unit. Depending on used measuring unit e.g. 123 means 123 $\mu\text{R/h}$ resp. 1.23 $\mu\text{Sv/h}$
AR6	Reading the threshold 2 for dose rate alarm. Response: see command AR5
AR7	Reading the threshold 1 for the dose alarm Response: number in μSv , $100\mu\text{R}$ or $100\mu\text{rem}$ units
AR8	Reading the threshold 2 for the dose alarm see command AR7
AR9	Reading sigma value Response: Number from 2..9
ARM	Reading minimum count rate for sigma alarm Response: Number in cps
ARN	Reading: NBR alarm threshold low energy Response: Number in 0.01% units Background preset count Response: Number in counts Background preset time Response: Number in seconds
ARP	Reading scaler parameter: Preset count. Response: Number in counts Preset time. Response: Number in seconds Scaler wait. time Response: Number in seconds
ARB	Reading scaler background values: Background value Response: Number in 0.01 cps units

	Used time for background value Response: Number in seconds
ARH	NBR alarm threshold high energy Response: Number in 0.01% units
SR3	Reading minimum count rate for NBR alarm Response: Number in cps
AW3Number	Setting the threshold 1 for count rate alarm. Number = Number in cps.
AW4Number	Setting the threshold 2 for count rate alarm. Number = see command AW3.
AW5Number	Setting the threshold 1 for dose rate alarm. Number = in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units. e. g. AW3123 means $123\mu\text{R/h}$.
AW6Number	Setting the threshold 2 for dose rate alarm. see command AW3
AW7Number	Setting the threshold 1 for dose alarm. Number = in $100\mu\text{R}$, $100\mu\text{rem}$ or $1\mu\text{Sv}$ units. e. g. AW7123 means 12.3mR .
AW8Number	Setting the threshold 2 for dose alarm. Number = see command AW7
AW9Number	Setting the sigma value. Number = form 2 to 9
AWMNumber	Setting minimum count rate for sigma alarm Number= value from 0 to 255
AWNNumber Number Number	Setting NBR alarm threshold level low energy Number=value in 0.01% units. From 1.00% to 2.55% Setting background preset count Number=counts. From 0 to 9999 Setting background preset time Number=seconds. From 0 to 9999s
AWHNumber	Setting NBR alarm threshold level high energy Number=value in 0.01% units. From 0% to 1.00%

AWBNumber	Number	Setting background value Number=0.01 cps units. From 0 to 10000 (100cps) Setting used time for background value Number=seconds. From 0 to 9999s
AWPNumber	Number	Setting scaler preset count Number=counts. From 0 to 9999 Setting scaler preset time Number=seconds. From 0 to 9999s
SW3	Number	Setting minimum count rate for NBR alarm Number= value from 0 to 255

5.2.2 Measurement values

Z	Read count rates with dead time correction Response: <ul style="list-style-type: none"> - Counter 1 in cps - Counter 2 in cps - Counter 3 in cps - Counter 4 in cps - PMT current index
z	Read filtered count rate 1 Response: Number in cps

Configuration flags

Configuration flags 1 with kR / kW

Bit number			
0	not used, write “0”		
1			
2			
3	Alarming Sound	0: off	1: on
4	Alarming LED	0: off	1: on
5	Alarming Vibration	0: off	1: on
6	NBR	0: off	1: on
7	Single Pulse	0: off	1: on

Configuration flags 2 with fR / fW

Bit number	
0	Disable keylock: 0: no 1:yes
1	not used, write “0”
2	not used, write “0”
3	Alarm threshold read-only 0: off 1:on
4	Flag for overload (readonly)
5	Temp.display 0: off 1:on
6	Temperature unit 0: °C 1: °F
7	In cps-mode display of dose rate 0: no 1: yes

5.2.3.3 *Configuration flags 3 with KR / KW*

Bit number	
0	Key tone 0:off 1:on
1	Autosend 0:off 1:on
2	Finder 0:off 1:on
3	Mean-Max 0:off 1:on
4	not used, write "0"
5	display of dose 0:off 1:on
6	Use click divider for single pulse 0:off 1:on
7	not used, write "0"

5.2.3.4 *Configuration flags 4 with jR / jW*

Bit number	
0	Safety-Alarm 0:off 1:on
1	Battery type 0:Alkaline 1:NiMh
2	Display rotation 0:No 1:yes
3	Scaler mode 0: preset count 1: preset time
4	Scaler "Auto restart" 0:off 1:on
5	Net Scaler 0:off 1:on
6	Net Ratemeter 0:off 1:on
7	not used, write "0"
8	not used, write "0"
9	not used, write "0"
10	not used, write "0"
11	not used, write "0"
12	not used, write "0"
13	not used, write "0"
14	not used, write "0"
15	not used, write "0"

5.2.3.5 *Measuring unit with uR /uW*

uR Read measuring unit and operation mode
 Response: Hex-value. See below
uRHex Write measuring unit and operation mode
 See below

Bit number	
0	0x00 Display unit cps 0x05 Display unit Sv/h
1	0x08 Display unit Bq/g 0x06 Display unit R/h
2	0x07 Display unit rem/h
3	
4	0: Ratemeter, 1: Scaler
5	0x01 Display unit Sv/h
6	0x02 Display unit R/h 0x03 Display unit rem/h
7	Fixed to "0"

5.2.3.6 *Menu configuration*

mR Read menu configuration for

- Main menu
- Submenu “Settings”
- Submenu “Alarm indication”

Response: Hex-values. See below

mWHex Hex Hex Write menu configuration for

- Main menu
- Submenu “Settings”
- Submenu “Alarm indication”

See below

5.2.3.6.1 Main menu

Bit number			
0	Switch off	0:hidden	1:visible
1	Background	0: hidden	1:visible
2	Not used, write “0”		
3	Backlight	0: hidden	1:visible
4	Measuring unit	0: hidden	1:visible
5	Operation mode	0: hidden	1:visible
6	Scaler parameter	0: hidden	1:visible
7	Nuclide table	0: hidden	1:visible
8	Alarm count rate	0: hidden	1:visible
9	Alarm dose rate	0: hidden	1:visible
10	Alarm dose	0: hidden	1:visible
11	Settings	0: hidden	1:visible
12	Alarm indication	0: hidden	1:visible
13	Show alarm	0: hidden	1:visible
14	Text info	0: hidden	1:visible
15..31	not used, write “0”		

5.2.3.6.2 Submenu “Settings”

Bit number			
0	Batt. type	0:hidden	1:visible
1	Autosend	0: hidden	1:visible
2	Single Pulse	0: hidden	1:visible
3	Finder	0: hidden	1:visible
4	Set Date/Time	0: hidden	1:visible
5	Lu-Test	0: hidden	1:visible
6..15	not used, write “0”		

5.2.3.6.3 Submenu “Alarm indication”

Bit number			
0	Sound	0:hidden	1:visible
1	LED	0: hidden	1:visible
2	Vibrator	0: hidden	1:visible
3...8	not used, write “0”		

5.2.4 High voltage

HR	Reading high voltage bit value Response: Bit value 0...255.
hR	Reading high voltage correction bit date and time of last successful Luthetium check Response: Number 0...255 with offset 128 and date and time as YYMMDDhhmm
hWNumber	Setting of high voltage correction bit Number from 0...255

5.2.5 Dead time correction

x	Read dead time Response: dead time in ns for - Counter 1 (Rate 1) - Counter 2 (Rate 2) - Counter 3 (Rate 3) - Counter 4 (Rate 4)
---	---

5.2.6 History output

5.2.6.1 History readout

Ratemeter:

256 716612088 721 12 120 23

```
---  ---  ---  ---  ---  ---
|      |      |      |      |      |
|      |      |      |      | \--- Temperature in °C
|      |      |      |      | \----- Measuring time in seconds
|      |      |      | \----- Max value
|      |      |      |      Count rate in 0.01 cps, cpm
|      |      |      |      Activity in 0.01 Bq, dps, dpm, Bq/cm²
|      |      |      |      Dose rate µR/h, rem/h or 0.01µSv/h
|      |      | \----- Mean value
|      |      |      Count rate in 0.01 cps, cpm
|      |      |      Activity in 0.01 Bq, dps, dpm, Bq/cm²
|      |      |      Dose rate µR/h, rem/h or 0.01µSv/h
|      | \----- Date and Time as a decimal value
| \----- History status as a decimal value
| see below
```

Scaler:

15 11 166 105 +19

```

-- -- -- -- --
| | | | |
| | | | \--- Temperature in °C
| | | \----- Rate mean in µR/h, µrem/h or 0.01µSv/h units
| | \----- Rate max in µR/h, µrem/h or 0.01µSv/h units
| \----- Counter 1 mean
\----- Counter 1 max

```

End of History:

End

5.2.6.2 *History status*

decimal value converted in HEX:

0x40000078

```

| ----
| |
| \-- History cycle time (0x78 -> 120s)
|
\----- 4: identifier for status information
          6: identifier for status information +
            last power on without power off

```

5.2.7 Event log

6656 520549251

```

|      |
|      \-----
\-----

```

Date and time as a decimal value (see 3.2)

Event log as a decimal value

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	0x01 Level display
9	0x02: Display count rate
10	0x04: Display dose rate
11	Sound 0: off 1:on
12	LED 0: off 1:on
13	Vibration alarm 0: off 1:on
14	1: Dose cleared
15	1: Alarm threshold changed
16	1: Count rate, dose rate or level-alarm
17	1: Dose alarm
18	1: Safety alarm
19	Not used
20	Value > alarm threshold 1 (depending on count rate, dose rate or level display)
21	Value > alarm threshold 2 (depending on count rate, dose rate or level display)
22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24	Low energy alarm
25	High energy alarm
26	Power off
27	Power on
28	NBR-alarm
29	Not used
30	Not used
31	Not used

5.2.8 Automatic sending

<STX>7 2 0 0 14 FH41PR 123 AB<ETX><CR><LF>

```

| | | | | | | |
| | | | | | | \- BCC in HEX
| | | | | | | \----- Dose in µR/h
| | | | | | | \----- information RadEye PRD
| | | | | | | \----- Status as Hex-value
| | | | | | | \----- „0“ not used
| | | | | | | \----- „0“ not used
| | | | | | | \----- „2“ not used
| | | | | | | \----- Dose rate in µR/h

```

Formation of a BCC (block check character):

Modulo 256 sum of <STX> up to the last character before the BCC (including), coded as hexadecimal ASCII-number (e.g. 1F).

Status

Bit number	
0	Not used
1	1: Overload
2	1: Count rate, dose rate or level-alarm
3	1: Dose alarm
4	1: NBR alarm
5	1: Battery voltage low
6	not used
7	not used

5.2.9 Status information

F

Reading status information

Response Number with status information

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	0x01 Level display
9	0x02: Display count rate
10	0x04: Display dose rate
11	not used,
12	Alarm threshold read-only 0: off 1:on
13	Flag for overload
14	Temperature display 0: off 1:on
15	In cps-mode display of dose rate 0: no 1: yes

16	1: Count rate, dose rate or level-alarm
17	1: Dose alarm
18	1: Safety alarm
19	Not used
20	Value > alarm threshold 1 (depending on count rate, dose rate or level display)
21	Value > alarm threshold 2 (depending on count rate, dose rate or level display)
22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24	Low energy alarm
25	High energy alarm
26	Not used
27	Not used
28	Not used
29	Not used
30	Not used
31	Fixed to "0"

5.3 RadEye PRD / PRD-ER / PRD-S / PRD-ER-S (>V3.00)

Used firmware version: 3.05

5.3.1 Limit values

AR0	Reading computed sigma alarm threshold 1 . Response: number in cps.
AR1	Reading the threshold 1 for level alarm. Response: number from 2..9. Number 2 means 32 cps. Each subsequent level has twice as much cps as the preceding level.
AR2	Reading the threshold 2 for level alarm. Response: see command AR1
AR3	Reading the threshold 1 for count rate alarm. Response: number in cps.
AR4	Reading the threshold 2 for count rate alarm. Response: number in cps.
AR5	Reading the threshold 1 for dose rate alarm. Response: number in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ unit. Depending on used measuring unit e.g. 1 2 3 means 1 2 3 $\mu\text{R/h}$ resp. 1.23 $\mu\text{Sv/h}$
AR6	Reading the threshold 2 for dose rate alarm. Response: see command AR5
AR7	Reading the threshold 1 for the dose alarm Response: number in μSv , $100\mu\text{R}$ or $100\mu\text{rem}$ units
AR8	Reading the threshold 2 for the dose alarm see command AR7
AR9	Reading sigma value Response: Number from 2..9
ARM	Reading minimum count rate for sigma alarm Response: Number in cps
ARN	NBR alarm threshold low energy Response: Number in 0.01% units
ARH	NBR alarm threshold high energy Response: Number in 0.01% units

ARk	Reading user display contrast Return: Number from -15...+15
SR3	Reading minimum count rate for NBR alarm Response: Number in cps
AW1Number	Setting the threshold 1 for level alarm. Number = 2...9. Number 2 means 32 cps. Each subsequent level has twice as much cps as the preceding level
AW2Number	Setting the threshold 2 for level alarm. Number = 2...9. Number 2 means 32 cps. Each subsequent level has twice as much cps as the preceding level
AW3Number	Setting the threshold 1 for count rate alarm. Number = Number in cps.
AW4Number	Setting the threshold 2 for count rate alarm. Number = see command AW3.
AW5Number	Setting the threshold 1 for dose rate alarm. Number = in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units. e. g. AW3123 means $123\mu\text{R/h}$.
AW6Number	Setting the threshold 2 for dose rate alarm. see command AW3
AW7Number	Setting the threshold 1 for dose alarm. Number = in $100\mu\text{R}$, $100\mu\text{rem}$ or $1\mu\text{Sv}$ units. e. g. AW7123 means 12.3mR .
AW8Number	Setting the threshold 2 for dose alarm. Number = see command AW7
AW9Number	Setting the sigma value. Number = form 2 to 9
AWMNumber	Setting minimum count rate for sigma alarm Number= value from 0 to 255
AWNNumber	Setting NBR alarm threshold level low energy Number=value in 0.01% units. From 1.00% to 2.55%
AWHNumber	Setting NBR alarm threshold level high energy Number=value in 0.01% units. From 0% to 1.00%

<i>AWkNumber</i>	Set user display contrast Number=value from -15 to +15
<i>SW3 Number</i>	Setting minimum count rate for NBR alarm Number= value from 0 to 255

5.3.2 Measurement values

<i>Z</i>	Read raw count rates with dead time correction Response: <ul style="list-style-type: none"> - Counter 1 in cps - Counter 2 in cps - Counter 3 in cps - Counter 4 in cps - HV power in cps - PMT current in 0.1uA units
<i>z</i>	Read filtered count rate 1 Response: Number in cps

5.3.3 Configuration flags

5.3.3.1 Configuration flags 1 with *kR* / *kW*

Bit number	
0	Disable key lock 0: no 1: yes
1	Show live graphic 0: no 1: yes
2	
3	Alarming Sound 0: off 1: on
4	Alarming LED 0: off 1: on
5	Alarming Vibration 0: off 1: on
6	NBR 0: off 1: on
7	Single Pulse 0: off 1: on

5.3.3.2 Configuration flags 2 with *fR* / *fW*

Bit number	
0	not used, write "0"
1	
2	
3	Alarm threshold read-only 0: off 1:on
4	Flag for overload (read only)
5	Temp.display 0: off 1:on
6	Temperature unit 0: °C 1: °F
7	In cps-mode display of dose rate 0: no 1: yes

5.3.3.3 *Configuration flags 3 with KR / KW*

Bit number	
0	Key tone 0:off 1:on
1	Autosend 0:off 1:on
2	Finder 0:off 1:on
3	Mean-Max 0:off 1:on
4	not used, write "0"
5	display of dose 0:off 1:on
6	Single pulse click divider 0:off 1:on
7	not used, write "0"

5.3.3.4 *Configuration flags 4 with jR / jW*

Bit number	
0	Safety-Alarm 0:off 1:on
1	Battery type 0:Alkaline 1:NiMh
2	Display rotation 0:No 1:yes
3	Scaler mode 0: preset count 1: preset time
4	Scaler "Auto restart" 0:off 1:on
5	Net Scaler 0:off 1:on
6	not used, write "0"
7	Net Ratemeter 0:off 1:on
8	not used, write "0"
9	Rated alarm 0:off 1:on
10	not used, write "0"
11	Alarm LCD-LED 0:off 1:on
12	not used, write "0"
13	not used, write "0"
14	not used, write "0"
15	not used, write "0"

5.3.3.5 *Measuring unit with uR /uW*

uR Read measuring unit and operation mode
 Response: Hex-value. See below
uRHex Write measuring unit and operation mode
 See below

Bit number	
0	0x00 Display unit cps 0x05 Display unit Sv/h
1	0x02 Display unit Bq 0x06 Display unit R/h
2	0x07 Display unit rem/h
3	
4	0: Ratemeter, 1: Scaler
5	0x01 Dose rate unit Sv/h
6	0x02 Dose rate unit R/h 0x03 Dose rate unit rem/h
7	Fixed to "0"

5.3.3.6 *Menu configuration*

5.3.3.6.1 Main menu RadEye PRD/PRD-ER

Bit number			
0	Switch off	0:hidden	1:visible
1	Sound	0: hidden	1:visible
2	LED	0: hidden	1:visible
3	Vibrator	0: hidden	1:visible
4	Level	0: hidden	1:visible
5	Count Rate	0: hidden	1:visible
6	Dose Rate	0: hidden	1:visible
7	Alarm Counter	0: hidden	1:visible
8	Alarm Dose Rate	0: hidden	1:visible
9	Alarm Dose	0: hidden	1:visible
10	Alarm Level	0: hidden	1:visible
11	Autosend	0: hidden	1:visible
12	Alarm-NBR	0: hidden	1:visible
13	Clear Dose	0: hidden	1:visible
14	Finder	0: hidden	1:visible
15	Single Pulse	0: hidden	1:visible
16	Backlight	0: hidden	1:visible
17	Show Alarm	0: hidden	1:visible
18	Settings	0: hidden	1:visible
19	Text Info	0: hidden	1:visible
20	Bluetooth	0: hidden	1:visible
21..31	not used, write "0"		

5.3.3.6.2 Submenu "Settings" RadEye PRD / PRD-ER

Bit number			
0	Batt. type	0:hidden	1:visible
1	Set Date/Time	0: hidden	1:visible
2	Lu-Test	0: hidden	1:visible
3	Language	0: hidden	1:visible
4	Rated alarm (from V3.05)	0: hidden	1:visible
5	Contrast (from V3.05)	0: hidden	1:visible
7..15	not used, write "0"		

5.3.3.6.3 Main menu RadEye PRD-S/PRD-ER-S

Bit number			
0	Switch off	0:hidden	1:visible
1	Background	0: hidden	1:visible
2	Not used, write "0"		
3	Backlight	0: hidden	1:visible
4	Measuring unit	0: hidden	1:visible
5	Operation mode	0: hidden	1:visible
6	Scaler parameter	0: hidden	1:visible
7	Nuclide table	0: hidden	1:visible
8	Alarm count rate	0: hidden	1:visible
9	Alarm dose rate	0: hidden	1:visible
10	Alarm-NBR	0: hidden	1:visible
11	Alarm dose	0: hidden	1:visible
12	Clear dose	0: hidden	1:visible
13	Settings	0: hidden	1:visible
14	Alarm indication	0: hidden	1:visible
15	Show alarm	0: hidden	1:visible
16	Text info	0: hidden	1:visible
17	Bluetooth	0: hidden	1:visible
18..29	not used, write "0"		
30	not used, write "1"		
31	not used, write "0"		

5.3.3.6.4 Submenu "Settings" RadEye PRD-S / PRD-ER-S

Bit number			
0	Batt. type	0:hidden	1:visible
1	Autosend	0: hidden	1:visible
2	Single Pulse	0: hidden	1:visible
3	Finder	0: hidden	1:visible
4	Set Date/Time	0: hidden	1:visible
5	Lu-Test	0: hidden	1:visible
6	Language	0: hidden	1:visible
7	Rated alarm (from V3.05)	0: hidden	1:visible
8	Contrast (from V3.05)	0: hidden	1:visible
9..15	not used, write "0"		

5.3.3.6.5 Submenu “Alarm indication”

Bit number			
0	Sound	0:hidden	1:visible
1	LED	0: hidden	1:visible
2	Vibrator	0: hidden	1:visible
3	LCD-LED	0: hidden	1:visible
4...8	not used, write “0”		

5.3.4 High voltage

HR Reading high voltage bit value
 Response: value in Volt.

hR Reading high voltage correction bit date and time of last successful Lutetium check
 Response: Number 0...255 with offset 128
 and date and time as YYMMDDhhmm

hW*Number* Set of high voltage correction
 Number from 0...255

5.3.5 Dead time correction

x Read dead time
 Response: dead time in ns for
 - Counter 1 (Rate 1)
 - Counter 2 (Rate 2)
 - Counter 3 (Rate 3)

5.3.6 History output

5.3.6.1 History readout

Following readout:

1536 716612088 1239 1600 30 5 120 23 4

```

-----
|      |      |      |      |      |      |      |
|      |      |      |      |      |      | \---- Temperature in °C
|      |      |      |      |      | \----- Measuring time in seconds
|      |      |      |      | \----- Max value dose rate
|      |      |      |      |      | μR/h, μrem/h, 0.01μSv/h
|      |      |      |      | \----- Mean value dose rate
|      |      |      |      |      | 0.1μR/h, 0.1μrem/h, 0.001μSv/h
|      |      |      | \----- Max value
|      |      |      |      | Count rate in 0.01 cps
|      |      |      |      | Activity in 0.01 Bq
|      |      | \----- Mean value
|      |      |      | Count rate in 0.01 cps
|      |      |      | Activity in 0.01 Bq
|      | \----- Date and Time as a decimal value
\----- History status as a decimal value

```

see below

End of History:

End

5.3.6.2 History status

decimal value converted in HEX:

Bit number	
0	Ratemeter/Scaler net value 0: no 1: yes
1	Operation mode Scaler 0: no 1: yes
2	Not used, read as '0'
3	Not used, read as '0'
4	Background measurement 0: no 1: yes
5	Preset time 0: no 1: yes
6	Not used, read as '0'
7	Not used, read as '0'
8..10	5:Sv/h, 6:R/h, 7:rem/h
11	Contamination (Bq)
12..15	Not used, read as '0'

5.3.7 Event log

6656 520549251

| |
| \-----
\-----

Date and time as a decimal value (see 3.2)

Event log as a decimal value

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	0x01 Level display
9	0x02: Display count rate
10	0x04: Display dose rate
11	Sound 0: off 1:on
12	LED 0: off 1:on
13	Vibration alarm 0: off 1:on
14	1: Dose cleared
15	1: Alarm threshold changed
16	1: Count rate, dose rate or level-alarm
17	1: Dose alarm
18	1: Safety alarm
19	Not used
20	Value > alarm threshold 1 (depending on count rate, dose rate or level display)
21	Value > alarm threshold 2 (depending on count rate, dose rate or level display)
22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24	Low energy alarm
25	High energy alarm
26	Power off
27	Power on
28	NBR-alarm
29	Not used
30	Not used
31	Not used

5.3.8 Nuclide calibration data (from V3.04)

nRNumber	Reading calibration data. Number: consecutive number Response: Data (see below)
nWNumberString	Write calibration data Number: consecutive number String: nuclide data e.g. nW02Co-60 50
nRA	Reading number of stored nuclides. Response: value from 0...15
nWANumber	Write the number of stored nuclides. Number: value from 0...15
nRG	Reading active nuclide. Response: value from 0...15
nWANumber	Write the number of active nuclides. Number: value from 0...15
yR	Read custom measuring unit (e.g. "g" for Bq/g)
yWText	Write custom measuring unit (e.g. "g" for Bq/g)

Nuclide data:

Co-60 20

```

-----
|      |
|      \----- Factor for activity calculation in (Bq)/cps
\----- Nuclide name. Up to 14 characters. Do not use space (0x20)

```

5.3.9 Automatic sending

```

<STX>7 2 9 5 14 FH41PR 123 AB<ETX><CR><LF>
      | | | | |      |      |      |
      | | | | |      |      |      \- BCC in HEX
      | | | | |      |      \----- Dose in µR/h
      | | | | |      \----- information RadEye
      | | | | \----- Status as Hex-value
      | | | \----- „5“ not used
      | | \----- Count rate (cps)
      | \----- „2“ not used
      \----- Dose rate in µR/h

```

Formation of a BCC (block check character):

Modulo 256 sum of <STX> up to the last character before the BCC (including), coded as hexadecimal ASCII-number (e.g. 1F).

Information RadEye

Reading	
FH41PR	RadEye PRD
PRDER	RadEye PRD-ER
PRDS	RadEye PRD-S
PRDERS	RadEye PRD-ER-S

Status

Bit number	
0	Not used
1	1: Overload
2	1: Count rate, dose rate or level-alarm
3	1: Dose alarm
4	1: NBR alarm
5	1: Battery voltage low
6	not used
7	not used

5.3.10 Status information

F Reading status information
 Response Number with status information

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	0x01 Level display
9	0x02: Display count rate
10	0x04: Display dose rate
11	not used,
12	Alarm threshold read-only 0: off 1: on
13	Flag for overload
14	Temperature display 0: off 1: on
15	In cps-mode display of dose rate 0: no 1: yes
16	1: Count rate, dose rate or level-alarm
17	1: Dose alarm
18	1: Safety alarm
19	Not used
20	Value > alarm threshold 1 (depending on count rate, dose rate or level display)
21	Value > alarm threshold 2 (depending on count rate, dose rate or level display)
22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24	Low energy alarm
25	High energy alarm
26	Not used
27	Not used
28	Not used
29	Not used
30	Not used
31	Fixed to "0"

5.4 RadEye G/G-10

Used firmware version: V1.52

5.4.1 Limit values

AR5	Reading the threshold 1 for dose rate alarm. Response: number in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ unit. Depending on used measuring unit e.g. 123 means 123 $\mu\text{R/h}$ resp. 1.23 $\mu\text{Sv/h}$
AR6	Reading the threshold 2 for dose rate alarm. Response: see command AR5
AR7	Reading the threshold 1 for the dose alarm Response: number in μSv , 100 μR or 100 μrem units
AR8	Reading the threshold 2 for the dose alarm see command AR7
AW5Number	Setting the threshold 1 for dose rate alarm. Number = in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units. e. g. AW3123 means 123 $\mu\text{R/h}$.
AW6Number	Setting the threshold 2 for dose rate alarm. see command AW3
AW7Number	Setting the threshold 1 for dose alarm. Number = in 100 μR , 100 μrem or 1 μSv units. e. g. AW7123 means 12.3mR.
AW8Number	Setting the threshold 2 for dose alarm. Number = see command AW7

5.4.2 Measurement values

Z	Read raw count rates with dead time correction Response: <ul style="list-style-type: none">- Counter 1 in cps- TTP value- HV power in cps
---	--

5.4.3 Configuration flags

5.4.3.1 Configuration flags 1 with *kR* / *kW*

Bit number	
0	Fixed to “0”
1	Fixed to “0”
2	Fixed to “1”
3	Alarming Sound 0: off 1: on
4	Alarming LED 0: off 1: on
5	Alarming Vibration 0: off 1: on
6	not used, write “0”
7	Single Pulse 0: off 1: on

5.4.3.2 Configuration flags 2 with *fR* / *fW*

Bit number	
0	0x01: RadEye G, measuring unit R/h
1	0x10: RadEye G-10, measuring unit Sv/h 0x11: RadEye G-10, measuring unit rem/h
2	Fixed to “1”
3	Alarm threshold read-only 0: off 1: on
4	Flag for overload (readonly)
5	Temp.display 0: off 1: on
6	Temperature unit 0: °C 1: °F
7	not used, write “0”

5.4.3.3 *Configuration flags 3 with KR / KW*

Bit number	
0	Key tone 0:off 1:on
1	Autosend 0:off 1:on
2	Finder 0:off 1:on
3	Mean-Max 0:off 1:on
4	not used, write "0"
5	display of dose 0:off 1:on
6	Fixed to "1"
7	not used, write "0"

5.4.3.4 *Configuration flags 4 with jR / jW*

Bit number	
0	not used, write "0"
1	Battery type 0:Alkaline 1:NiMh
2	Display rotation 0:No 1:yes
3	not used, write "0"
4	not used, write "0"
5	not used, write "0"
6	not used, write "0"
7	not used, write "0"
8	not used, write "0"
9	not used, write "0"
10	not used, write "0"
11	not used, write "0"
12	not used, write "0"
13	not used, write "0"
14	not used, write "0"
15	not used, write "0"

5.4.3.5 *Menu configuration*

Bit number			
0	Switch off	0:hidden	1:visible
1	Sound	0: hidden	1:visible
2	LED	0: hidden	1:visible
3	Vibrator	0: hidden	1:visible
4	not used, write "0"		
5	not used, write "0"		
6	not used, write "0"		
7	not used, write "0"		
8	Alarm Dose Rate	0: hidden	1:visible
9	Alarm Dose	0: hidden	1:visible
10	not used, write "0"		
11	Autosend	0: hidden	1:visible
12	not used, write "0"		
13	Clear Dose	0: hidden	1:visible
14	Finder	0: hidden	1:visible
15	Single Pulse	0: hidden	1:visible
16	Backlight	0: hidden	1:visible
17	Show Alarm	0: hidden	1:visible
18	Settings	0: hidden	1:visible
19	Text Info	0: hidden	1:visible
20	not used, write "0"		
21	not used, write "0"		
22	not used, write "0"		
23	not used, write "0"		
24	not used, write "0"		
25	not used, write "0"		
26	not used, write "0"		
27	not used, write "0"		
28	not used, write "0"		
29	not used, write "0"		
30	not used, write "0"		
31	not used, write "0"		

5.4.4 **High voltage**

HR Reading high voltage bit value
Response: Bit value 0...255.

5.4.5 **Dead time correction**

x Read dead time
Response: dead time in ns

5.4.6 History output

5.4.6.1 History readout

History status for 1. readout and change of history cycle time

1073741944 520534610

```
-----
|           |
|           |
|           | \----- Date and Time as a decimal value
| \----- History status as a decimal value
```

following readout:

15 11 +19

```
-- -- --
| | |
| | | \--- Temperature in °C
| \----- Rate mean in µR/h, µrem/h or 0.01µSv/h units
| \----- Rate max in µR/h, µrem/h or 0.01µSv/h units
```

End of History:

End

5.4.6.2 History status

decimal value converted in HEX:

0x40000078

```
|      ----
|      |
|      | \-- History cycle time (0x78 -> 120s)
|
| \----- 4: identifier for status information
|           6: identifier for status information +
|           last power on without power off
```

5.4.7 Event log

6656 520549251

```

|      |
|      \-----
\-----

```

Date and time as a decimal value (see 3.2)

Event log as a decimal value

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	Fixed to "0"
9	Fixed to "0"
10	Fixed to "1"
11	Sound 0: off 1:on
12	LED 0: off 1:on
13	Vibration alarm 0: off 1:on
14	1: Dose cleared
15	1: Alarm threshold changed
16	1: Alarm dose rate
17	1: Alarm dose
18	Not used
19	Not used
20	Dose rate > alarm threshold 1
21	Dose rate > alarm threshold 2
22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24	Not used
25	Not used
26	Power off
27	Power on
28	Not used
29	Not used
30	Not used
31	Fixed to "0"

5.4.8 Automatic sending

<STX>7 2 0 0 14 FH41B2 123 AB<ETX><CR><LF>

```

| | | | | | | |
| | | | | | | \- BCC in HEX
| | | | | | | \----- Dose in µR/h
| | | | | | | \----- information RadEye G
| | | | | | | \----- Status as Hex-value
| | | | | | | \----- „0“ not used
| | | | | | | \----- „0“ not used
| | | | | | | \----- „2“ not used
| | | | | | | \----- Dose rate in µR/h

```

Formation of a BCC (block check character):

Modulo 256 sum of <STX> up to the last character before the BCC (including), coded as hexadecimal ASCII-number (e.g. 1F).

Status

Bit number	
0	Not used
1	1: Overload
2	1: Alarm dose rate
3	1: Alarm dose
4	not used
5	1: Battery voltage low
6	not used
7	not used

5.4.9 Status information

F

Reading status information

Response Number with status information

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	0x04: Display dose rate
9	
10	
11	not used
12	Alarm threshold read-only 0: off 1:on
13	Flag for overload
14	Temperature display 0: off 1:on
15	Fixed to “0”

16	1: Alarm dose rate
17	1: Alarm dose
18	Fixed to "0"
19	Not used
20	Dose rate > alarm threshold 1
21	Dose rate > alarm threshold 2
22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24	Low energy alarm
25	High energy alarm
26	Not used
27	Not used
28	Not used
29	Not used
30	Not used
31	Not used

5.5 RadEye N

Used firmware version: 1.53

5.5.1 Limit values

AR3	Reading the threshold 1 for count rate alarm. Response: number in cps.
AR4	Reading the threshold 2 for count rate alarm. Response: number in cps.
AR5	Reading the threshold 1 for dose rate alarm. Response: number in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ unit. Depending on used measuring unit e.g. 123 means 123 $\mu\text{R/h}$ resp. 1.23 $\mu\text{Sv/h}$
AR6	Reading the threshold 2 for dose rate alarm. Response: see command AR5
AR7	Reading the threshold 1 for the dose alarm Response: number in μSv , 100 μR or 100 μrem units
AR8	Reading the threshold 2 for the dose alarm see command AR7
AW3Number	Setting the threshold 1 for count rate alarm. Number = Number in cps.
AW4Number	Setting the threshold 2 for count rate alarm. Number = see command AW3.
AW5Number	Setting the threshold 1 for dose rate alarm. Number = in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units. e. g. AW3123 means 123 $\mu\text{R/h}$.
AW6Number	Setting the threshold 2 for dose rate alarm. see command AW3
AW7Number	Setting the threshold 1 for dose alarm. Number = in 100 μR , 100 μrem or 1 μSv units. e. g. AW7123 means 12.3mR.
AW8Number	Setting the threshold 2 for dose alarm. Number = see command AW7

5.5.2 Measurement values

Z	Read raw count rates with dead time correction Response: <ul style="list-style-type: none">- Counter 1 in cps- Counter 2 in cps- Counter 3 in cps- HV power in cps
z	Read filtered count rate 1 Response: Number in 0.01 cps units
A	Read display value and status Response: <ul style="list-style-type: none">- display value in 0.1 cps or $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units- Status (see 5.5.6.2)

5.5.3 Configuration flags

5.5.3.1 Configuration flags 1 with *kR* / *kW*

Bit number	
0	Fixed to “0”
1	Fixed to “0”
2	Fixed to “0”
3	Alarming Sound 0: off 1: on
4	Alarming LED 0: off 1: on
5	Alarming Vibration 0: off 1: on
6	Fixed to “0”
7	Single Pulse 0: off 1: on

5.5.3.2 Configuration flags 2 with *fR* / *fW*

Bit number	
0	Fixed to “0”
1	Fixed to “0”
2	not used, write “0”
3	Alarm threshold read-only 0: off 1: on
4	Flag for overload (readonly)
5	Temp.display 0: off 1: on
6	Temperature unit 0: °C 1: °F
7	Fixed to “0”

5.5.3.3 *Configuration flags 3 with KR / KW*

Bit number	
0	Key tone 0:off 1:on
1	Autosend 0:off 1:on
2	Finder 0:off 1:on
3	Mean-Max 0:off 1:on
4	not used, write "0"
5	display of dose 0:off 1:on
6	not used, write "0"
7	not used, write "0"

5.5.3.4 *Configuration flags 4 with jR / jW*

Bit number	
0	Fixed to "0"
1	Battery type 0:Alkaline 1:NiMh
2	Display rotation 0:No 1:yes
3	Scaler mode 0: Preset count 1: Preset time
4	Scaler display 0: mean value 1: accumulated counts
5	Scaler, after measurement 0: Stop 1: automatic restart
6	Scaler net 0:No 1:yes
7	Ratemeter net 0:No 1:yes
8	not used, write "0"
9	not used, write "0"
10	not used, write "0"
11	not used, write "0"
12	not used, write "0"
13	not used, write "0"
14	not used, write "0"
15	not used, write "0"

5.5.3.5 *Measuring unit with uR /uW*

uR Read measuring unit and operation mode
Response: Hex-value. See below
uRHex Write measuring unit and operation mode
See below

Bit number	
0	0x00 Display unit cps 0x05 Display unit Sv/h
1	0x01 Display unit cpm 0x06 Display unit R/h
2	0x07 Display unit rem/h
3	
4	0: Ratemeter, 1: Scaler
5	Fixed to "0"
6	Fixed to "0"
7	Fixed to "0"

5.5.3.6 *Menu configuration*

mR	Read configuration for main menu and submenu “Settings” Response: Hex-values. See below
mRHex Hex	Write configuration for main menu and submenu “Settings” See below

5.5.3.6.1 Main menu

Bit number			
0	Switch off	0:hidden	1:visible
1	Background	0: hidden	1:visible
2	not used, write “0”		
3	Backlight	0: hidden	1:visible
4	Measuring unit	0: hidden	1:visible
5	Operation mode	0: hidden	1:visible
6	Scaler parameter	0: hidden	1:visible
7	Select DR factor	0: hidden	1:visible
8	Alarm Count Rate	0: hidden	1:visible
9	Alarm Dose Rate	0: hidden	1:visible
10	not used, write “0”		
11	Alarm Dose	0: hidden	1:visible
12	Clear Dose	0: hidden	1:visible
13	Settings	0: hidden	1:visible
14	Alarm indication	0: hidden	1:visible
15	Show alarm	0: hidden	1:visible
16	Text Info	0: hidden	1:visible
17..31	not used, write “0”		

5.5.3.6.2 Submenu “Settings”

Bit number			
0	Batt. type	0:hidden	1:visible
1	Autosend	0: hidden	1:visible
2	Single Pulse	0: hidden	1:visible
3	Finder	0: hidden	1:visible
4	Set Date/Time	0: hidden	1:visible
5..15	not used, write “0”		

5.5.3.6.3 Submenu “Alarm indication”

Bit number	
------------	--

0	Sound	0:hidden	1:visible
1	LED	0: hidden	1:visible
2	Vibrator	0: hidden	1:visible
3...8	not used, write "0"		

5.5.4 High voltage

HR Reading high voltage bit value
Response: Bit value 0...255.

hR Reading high voltage correction bit
Response: Number 0...255 with offset 128
and date and time as YYMMDDhhmm

hWNumber Setting of high voltage correction bit
Number from 0...255

5.5.5 Dead time correction

x Read dead time
Response: dead time in ns for
- Counter 1 (Rate 1)
- Counter 2 (Rate 2)
- Counter 3 (Rate 3)

5.5.6 History output

5.5.6.1 History readout

256 716612088 1 1 120 23

```
-----
|      |      | | | |
|      |      | | | \----- Temperature in °C
|      |      | | | \----- Measuring time in seconds
|      |      | | \----- Rate max in µR/h, µrem/h, 0.01µSv/h or
|      |      | | or 0.1 cps units
|      |      | \----- mean in µR/h, µrem/h or 0.01µSv/h units
|      |      | or accumulated counts in cps-mode
|      |      | \----- Date and Time as a decimal value
|      | \----- History status as a decimal value see below
\-----
```

End of History:
End

5.5.6.2 History status

decimal value converted in HEX:

```
0101 0101 0101 0101
|||| |||| |||| ||||
```

```

|||| |  ||| |  ||| |  ||| \-- Nulleffektabzug 0: Nein 1:Ja
|||| |  ||| |  ||| |  ||| \--- 1: Scaler
|||| |  ||| |  ||| \----- Nulleffektmessung 0: Nein    1:Ja
|||| |  ||| |  ||| \----- 0: Preset Counts  1: Preset Time
|||| |  ||| |  ||| \----- --
|||| |  ||| |  ||| \----- 1: Accumulated counts
|||| |  ++++----- gewähltes Nuklid max.15.
++++----- 0:cps, 1:cpm, 5:Sv/h, 6:R/h, 7:rem/h

```

Bit number	
0	Net value 0:No 1:Yes
1	Operation mode 0: Ratemeter 1:Scaler
2	not used
3	not used
4	Background measurement 0:No 1:Yes
5	Scaler with 0: Preset Counts 1: Preset Time
6	not used
7	Accumulated counts 0:No 1:Yes
8	Number of used calibration factor
9	
10	
11	
12	0x00 Display unit cps 0x05 Display unit Sv/h 0x01 Display unit cpm 0x06 Display unit R/h 0x07 Display unit rem/h
13	
14	
15	

5.5.7 Event log

6656 520549251

```

-----
|      |
|      | \----- Date and time as a decimal value (see 3.2)
|      | \----- Event log as a decimal value
\-----

```

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	Not used
9	
10	
11	Sound 0: off 1:on
12	LED 0: off 1:on
13	Vibration alarm 0: off 1:on
14	1: Dose cleared
15	1: Alarm threshold changed
16	1: Count rate or dose rate -alarm
17	1: Dose alarm
18	Not used
19	Not used
20	Value > alarm threshold 1 (depending on count rate or dose rate display)
21	Value > alarm threshold 2 (depending on count rate or dose rate display)
22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24	1: Scaler or Background parameter changed
25	Not used
26	Power off
27	Power on
28	Not used
29	Not used
30	Not used
31	Not used

5.5.8 Automatic sending

```

<STX>7 2 0 0 14 FH41PR 123 AB<ETX><CR><LF>
| | | | | | | |
| | | | | | | \- BCC in HEX
| | | | | | \----- Dose in µR/h
| | | | | \----- information RadEye PRD
| | | | \----- Status as Hex-value
| | | \----- Fixed to "0"
| | \----- Fixed to "0"
| \----- Fixed to "2"
\----- Dose rate in µR/h, µrem/h or 0.01µSv/h

```

Formation of a BCC (block check character):

Modulo 256 sum of <STX> up to the last character before the BCC (including), coded as hexadecimal ASCII-number (e.g. 1F).

Status

Bit number	
0	Not used
1	1: Overload
2	1: Count rate, dose rate or level-alarm
3	1: Dose alarm
4	Not used
5	1: Battery voltage low
6	not used
7	not used

5.5.9 Status information

F Reading status information
 Response Number with status information

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	Not used
9	
10	
11	not used,
12	Alarm threshold read-only 0: off 1:on
13	Flag for overload
14	Temperature display 0: off 1:on
15	Not used
16	1: Count rate or dose rate-alarm
17	1: Dose alarm
18	Not used
19	Not used
20	Value > alarm threshold 1 (depending on count rate or dose rate display)
21	Value > alarm threshold 2 (depending on count rate or dose rate display)
22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24	Not used
25	Not used
26	Not used
27	Not used
28	Not used
29	Not used
30	Not used
31	Not used

5.6 RadEye AB100 (V1.52)

Used firmware version: 1.52

5.6.1 Limit values

AR1	Reading the threshold 1 for activity alarm. Alpha/Beta channel Response: number in 0.01 Bq units.
AR2	Reading the threshold 1 for activity alarm. Alpha/Beta channel Response: number in 0.01 Bq units
AR3	Reading the threshold 1 for activity alarm. Alpha channel. Response: number in 0.01 Bq units.
AR4	Reading the threshold 2 for activity alarm. Alpha channel. Response: number in 0.01 Bq units.
AR5	Reading the threshold 1 for count rate alarm. Alpha/Beta channel. Response: number in cps
AR6	Reading the threshold 2 for count rate alarm. Alpha/Beta channel. Response number in cps
AR7	Reading the threshold 1 for count rate alarm. Alpha channel. Response: number in cps
AR8	Reading the threshold 2 for count rate alarm. Alpha channel. Response number in cps
AW1Number	Setting the threshold 1 for activity alarm. Alpha/Beta channel Number: value in 0.01 Bq units.
AW2Number	Setting the threshold 1 for activity alarm. Alpha/Beta channel Number: value in 0.01 Bq units
AW3Number	Setting the threshold 1 for activity alarm. Alpha channel. Number: value in 0.01 Bq units.
AW4Number	Setting the threshold 2 for activity alarm. Alpha channel. Number: value in 0.01 Bq units.
AW5Number	Setting the threshold 1 for count rate alarm. Alpha/Beta channel. Number: value in cps

<i>AW6Number</i>	Setting the threshold 2 for count rate alarm. Alpha/Beta channel. <i>Number</i> : value in cps
<i>AW7Number</i>	Setting the threshold 1 for count rate alarm. Alpha channel. <i>Number</i> : value in cps
<i>AW8Number</i>	Setting the threshold 2 for count rate alarm. Alpha channel. <i>Number</i> : value in cps

5.6.2 Measurement values

Z	Read raw count rates with dead time correction Response: <ul style="list-style-type: none"> - Counter 1 in cps - Counter 2 in cps - Counter 3 in cps - HV power in cps - filtered HV power in cps
z	Read filtered count rate alpha/beta and alpha channel Response: <ul style="list-style-type: none"> - Number in 0.01 cps units, alpha/beta channel - Number in 0.01 cps units, alpha channel
A	Read display value and status Response: <ul style="list-style-type: none"> - display value in 0.01 cps, 1μR/h, 1μrem/h, 0.01μSv/h or 0.01 Bq units - Status (see 5.5.6.2)

5.6.3 Configuration flags

5.6.3.1 Configuration flags 1 with *kR* / *kW*

Bit number	
0	Fixed to “0”
1	Fixed to “0”
2	Fixed to “0”
3	Alarming Sound 0: off 1: on
4	Alarming LED 0: off 1: on
5	Alarming Vibration 0: off 1: on
6	Fixed to “0”
7	Single Pulse 0: off 1: on

5.6.3.2 Configuration flags 2 with *fR* / *fW*

Bit number	
0	Fixed to “0”
1	Fixed to “0”
2	Fixed to “1”
3	Alarm threshold read-only 0: off 1: on
4	Flag for overload (readonly)
5	Display of temperature 0: off 1: on
6	Temperature unit 0: °C 1: °F
7	Fixed to “0”

5.6.3.3 *Configuration flags 3 with KR / KW*

Bit number	
0	Key tone 0:off 1:on
1	Autosend 0:off 1:on
2	Finder 0:off 1:on
3	Mean-Max 0:off 1:on
4	not used
5	not used
6	not used
7	not used

5.6.3.4 *Configuration flags 4 with jR / jW*

Bit number	
0	Not used
1	Battery type 0:Alkaline 1:NiMh
2	Display rotation 0:No 1:yes
3	Scaler mode 0: Preset count 1: Preset time
4	Scaler display 0: mean value 1: accumulated counts
5	Scaler, after measurement 0: Stop 1: automatic restart
6	Scaler net 0:No 1:yes
7	Ratemeter net 0:No 1:yes
8	Alpha LED 0:No 1:yes
9	Alpha pulse 0:No 1:yes
10	Beta only 0:No 1:yes
11	not used, write "0"
12	not used, write "0"
13	not used, write "0"
14	not used, write "0"
15	not used, write "0"

5.6.3.5 *Measuring unit with uR /uW*

uR Read measuring unit and operation mode
Response: Hex-value. See below
uRHex Write measuring unit and operation mode
See below

Bit number	
0	0x00 Display unit cps 0x03 Display unit dps
1	0x01 Display unit cpm 0x04 Display unit dpm
2	0x02 Display unit Bq 0x08 Display unit Bq/cm ²
3	
4	0: Ratemeter, 1: Scaler
5	0x00 Display alpha/beta channel
6	0x01 Display alpha channel 0x02 Dual display
7	Fixed to "0"

5.6.3.6 *Menu configuration*

mR Read configuration for main menu and submenu “Settings”
 Response: Hex-values. See below

mRHex Hex Write configuration for main menu and submenu “Settings”
 See below

Bit number			
0	Switch off	0:hidden	1:visible
1	Background	0: hidden	1:visible
2	not used, write “0”		
3	Backlight	0: hidden	1:visible
4	Measuring unit	0: hidden	1:visible
5	Operation mode	0: hidden	1:visible
6	Scaler parameter	0: hidden	1:visible
7	Nuclide table	0: hidden	1:visible
8	Alarm alpha/beta	0: hidden	1:visible
9	Alarm alpha/beta	0: hidden	1:visible
10	Alarm activity alpha/beta	0: hidden	1:visible
11	Alarm activity alpha	0: hidden	1:visible
12	not used, write “0”		
13	Settings	0: hidden	1:visible
14	Alarm indication	0: hidden	1:visible
15	Show alarm	0: hidden	1:visible
16	Text Info	0: hidden	1:visible
17	not used, write “0”		
18	not used, write “0”		
19	not used, write “0”		
20	not used, write “0”		
21	not used, write “0”		
22	not used, write “0”		
23	not used, write “0”		
24	not used, write “0”		
25	not used, write “0”		
26	not used, write “0”		
27	not used, write “0”		
28	not used, write “0”		
29	not used, write “0”		
30	not used, write “0”		
31	not used, write “0”		

Bit number			
0	Batt. type	0:hidden	1:visible
1	Autosend	0: hidden	1:visible

2	Single Pulse	0: hidden	1:visible
3	Finder	0: hidden	1:visible
4	Alpha LED	0: hidden	1:visible
5	Alpha Sound	0: hidden	1:visible
6	Set Date/Time	0: hidden	1:visible
7	Display alpha/beta or beta	0: hidden	1:visible
8	not used, write "0"		
9	not used, write "0"		
10	not used, write "0"		
11	not used, write "0"		
12	not used, write "0"		
13	not used, write "0"		
14	not used, write "0"		
15	not used, write "0"		

5.6.4 High voltage

HR Reading high voltage bit value
Response: Bit value 0...255.

hR Reading high voltage correction bit
Response: Number 0...255 with offset 128
and date and time as YYMMDDhhmm

hW*Number* Setting of high voltage correction bit
Number from 0...255

5.6.5 Dead time correction

x Read dead time
Response: dead time in ns for
- Counter 1 (Rate 1)
- Counter 2 (Rate 2)
- Counter 3 (Rate 3)

5.6.6 History output

5.6.6.1 History readout

```

256 716612088 721 999 12 50 120 23
---
|      |      |      |      |      |      |      |
|      |      |      |      |      |      | \---- Temperature in °C
|      |      |      |      |      |      | \----- Measuring time in seconds

```

```

|      |      |      |      | \----- Max count rate in 0.01 cps,
|      |      |      |      | Activity in 0.01 Bq
|      |      |      |      | alpha channel
|      |      |      |      | \----- Mean count rate in 0.01 cps,
|      |      |      |      | Activity in 0.01 Bq
|      |      |      |      | alpha channel
|      |      |      |      | \----- Max count rate in 0.01 cps,
|      |      |      |      | Activity in 0.01 Bq
|      |      |      |      | alpha/beta channel
|      |      |      |      | \----- Mean count rate in 0.01 cps,
|      |      |      |      | Activity in 0.01 Bq
|      |      |      |      | alpha/beta channel
|      |      |      |      | \----- Date and Time as a decimal value
\----- History status as a decimal value
see below

```

End of History:

End

5.6.6.2 *History status*

Decimal value converted in HEX:

Bit number			
0	Net value	0:No	1:Yes
1	Operation mode	0: Ratemeter	1:Scaler
2	Beta channel: display	0: Gross	1: Beta
3	not used		
4	Background measurement	0:No	1:Yes
5	Scaler with	0: Preset Counts	1: Preset Time
6	not used		
7	Accumulated counts	0:No	1:Yes
8	Number of used nuclide		
9			
10			
11			
12	0x00 Display unit cps	0x03 Display unit dps	
13	0x01 Display unit cpm	0x04 Display unit dpm	
14	0x02 Display unit Bq	0x08 Display unit Bq/cm²	
15			

5.6.7 Event log

6656 520549251

```

-----
|          |
|          | \----- Date and time as a decimal value (see 3.2)
|          | \----- Event log as a decimal value
\-----

```

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	Not used
9	
10	
11	Sound 0: off 1:on
12	LED 0: off 1:on
13	Vibration alarm 0: off 1:on
14	Not used
15	1: Alarm threshold changed
16	1: Count rate or activity -alarm
17	Not used
18	Not used
19	Not used
20	Value > alarm threshold 1
21	Value > alarm threshold 2
22	Not used
23	Not used
24	1: Scaler or Background parameter changed
25	Not used
26	Power off
27	Power on
28	Not used
29	Not used
30	Not used
31	Not used

5.6.8 Automatic sending

```

<STX>721  0 1 0 14 REAB10 123 AB<ETX><CR><LF>
|         | | | | |         |         |
|         | | | | |         |         | \- BCC in HEX
|         | | | | |         |         | \----- Dose in µR/h
|         | | | | |         |         | \----- information RadEye AB100
|         | | | | |         |         | \----- Status as Hex-value
|         | | | | |         |         | \----- Fixed to "0"
|         | | | | |         |         | \----- Display value alpha channel
|         | | | | |         |         | \----- Fixed to "o"
|         | | | | |         |         | \----- Display value alpha/beta channel

```

Formation of a BCC (block check character):

Modulo 256 sum of <STX> up to the last character before the BCC (including), coded as hexadecimal ASCII-number (e.g. 1F).

Status

Bit number	
0	Not used
1	1: Overload
2	1: Alarm
3	Not used
4	Not used
5	1: Battery voltage low
6	not used
7	not used

5.6.9 Status information

F Reading status information
 Response Number with status information

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	Not used
9	
10	
11	not used,
12	Alarm threshold read-only 0: off 1:on
13	Flag for overload
14	Temperature display 0: off 1:on
15	Not used
16	1: Alarm
17	Not used
18	Not used
19	Not used
20	Value > alarm threshold 1 (depending on count rate or dose rate display)
21	Value > alarm threshold 2 (depending on count rate or dose rate display)
22	Not used
23	Not used
24	Not used
25	Not used
26	Not used
27	Not used
28	Not used
29	Not used
30	Not used
31	Not used

5.7 RadEye B20/G20/G20-10

Used firmware version: 1.52

5.7.1 Limit values

AR1	Reading the threshold 1 for activity alarm. Response: number in 0.01 Bq units.
AR2	Reading the threshold 2 for activity alarm. Response: number in 0.01 Bq units.
AR3	Reading the threshold 1 for count rate alarm. Response: number in 0.01 cps units.
AR4	Reading the threshold 2 for count rate alarm. Response: number in 0.01 cps units.
AR5	Reading the threshold 1 dose rate alarm. Response: number in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units
AR6	Reading the threshold 2 for count rate alarm. Response: number in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units
AR7	Reading the threshold 1 for dose alarm. Response: number in μR , μrem or $0.01\mu\text{Sv}$ units
AR8	Reading the threshold 2 for dose alarm. Response number in μR , μrem or $0.01\mu\text{Sv}$ units
ARP	Reading scaler parameter: <ul style="list-style-type: none">- Preset count. Response: Number in counts- Preset time. Response: Number in seconds
AW1Number	Setting the threshold 1 for activity alarm Number: value in 0.01 Bq units.
AW2Number	Setting the threshold 2 for activity alarm Number: value in 0.01 Bq units
AW3Number	Setting the threshold 1 for count rate alarm Number: value in 0.01 cps units.
AW4Number	Setting the threshold 2 for count rate alarm Number: value in 0.01 cps units

AW5Number Setting the threshold 1 dose rate alarm.
Number: value in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units.

AW6Number Setting the threshold 2 dose rate alarm.
Number: value in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units.

AW7Number Setting the threshold 1 for dose alarm.
Number: value in μR , μrem or $0.01\mu\text{Sv}$ units

AW8Number Setting the threshold 2 for dose alarm.
Number: value in μR , μrem or $0.01\mu\text{Sv}$ units

5.7.2 Measurement values

Z Read raw count rates with dead time correction
Response:
 - Counter 1 in cps
 - HV power in cps

z Read filtered count rate
Response: Number in 0.01 cps units

A Read display value and status
Response:
 - display value Dose rate in $1\mu\text{R/h}$, $1\mu\text{rem/h}$, $0.01\mu\text{Sv/h}$ units
 or 0.01 cps, cpm, Bq, dps, dpm, Bq/cm^2 units
 - Status (see 5.5.6.2)

5.7.3 Configuration flags

5.7.3.1 Configuration flags 1 with *kR* / *kW*

Bit number	
0	Not used
1	Not used
2	Not used
3	Alarming Sound 0: off 1: on
4	Alarming LED 0: off 1: on
5	Alarming Vibration 0: off 1: on
6	Not used
7	Single Pulse 0: off 1: on

5.7.3.2 Configuration flags 2 with *fR* / *fW*

Bit number	
0	Not used
1	Not used
2	Not used, write “1”
3	Alarm threshold read-only 0: off 1:on
4	Flag for overload (readonly)
5	Display of temperature 0: off 1:on
6	Temperature unit 0: °C 1: °F
7	Not used

5.7.3.3 *Configuration flags 3 with KR / KW*

Bit number	
0	Key tone 0:off 1:on
1	Autosend 0:off 1:on
2	Finder 0:off 1:on
3	Mean-Max 0:off 1:on
4	not used, write "0"
5	display of dose 0:off 1:on
6	not used
7	not used

5.7.3.4 *Configuration flags 4 with jR / jW*

Bit number	
0	Not used
1	Battery type 0:Alkaline 1:NiMh
2	Display rotation 0:No 1:yes
3	Scaler mode 0: Preset count 1: Preset time
4	Not used
5	Scaler, after measurement 0: Stop 1: automatic restart
6	Scaler net 0:No 1:yes
7	Ratemeter net 0:No 1:yes
8	Not used
9	Not used
10	Not used
11	Not used
12	Not used
13	Not used
14	Not used
15	Not used

5.7.3.5 *Measuring unit with uR /uW*

uR Read measuring unit and operation mode
 Response: Hex-value. See below
 uRHex Write measuring unit and operation mode
 See below

Bit number	
0	0x00 Display unit cps 0x04 Display unit dpm
1	0x01 Display unit cpm 0x05 Display unit Sv/h
2	0x02 Display unit Bq 0x06 Display unit R/h
3	0x03 Display unit dps 0x07 Display unit rem/h 0x08 Display unit Bq/cm ²
4	0: Ratemeter, 1: Scaler
5	Last used display unit (dose rate)
6	0x00 Sv/h 0x01 R/h 0x02 rem/h
7	Not used

5.7.3.6 *Menu configuration*

mR Read configuration for main menu and submenu “Settings”
Response: Hex-values. See below

mRHex Hex Write configuration for main menu and submenu “Settings”
See below

Bit number			
0	Switch off	0:hidden	1:visible
1	Background	0: hidden	1:visible
2	not used, write “0”		
3	Backlight	0: hidden	1:visible
4	Measuring unit	0: hidden	1:visible
5	Operation mode	0: hidden	1:visible
6	Scaler parameter	0: hidden	1:visible
7	Nuclide table	0: hidden	1:visible
8	Alarm count rate	0: hidden	1:visible
9	Alarm dose rate	0: hidden	1:visible
10	Alarm activity	0: hidden	1:visible
11	Alarm dose	0: hidden	1:visible
12	Clear dose	0: hidden	1:visible
13	Settings	0: hidden	1:visible
14	Alarm indication	0: hidden	1:visible
15	Show alarm	0: hidden	1:visible
16	Text Info	0: hidden	1:visible
17	not used, write “0”		
18	not used, write “0”		
19	not used, write “0”		
20	not used, write “0”		
21	not used, write “0”		
22	not used, write “0”		
23	not used, write “0”		
24	not used, write “0”		
25	not used, write “0”		
26	not used, write “0”		
27	not used, write “0”		
28	not used, write “0”		
29	not used, write “0”		
30	not used, write “0”		
31	not used, write “0”		

Bit number			
0	Batt. type	0:hidden	1:visible
1	Autosend	0: hidden	1:visible
2	Single Pulse	0: hidden	1:visible
3	Finder	0: hidden	1:visible
4	Set Date/Time	0: hidden	1:visible
5	not used, write "0"		
6	not used, write "0"		
7	not used, write "0"		
8	not used, write "0"		
9	not used, write "0"		
10	not used, write "0"		
11	not used, write "0"		
12	not used, write "0"		
13	not used, write "0"		
14	not used, write "0"		
15	not used, write "0"		

5.7.4 High voltage

HR Reading high voltage bit value
Response: Bit value 0...255.

5.7.5 Dead time correction

x Read dead time
Response: dead time in ns

5.7.6 History output

5.7.6.1 History readout ratemeter

256 716612088 721 999 120 23

```
-----
|      |      |      |      |      |
|      |      |      |      | \--- Temperature in °C
|      |      |      | \----- Measuring time in seconds
|      |      | \----- Max value
|      |      |      Count rate in 0.01 cps, cpm
|      |      |      Activity in 0.01 Bq, dps, dpm, Bq/cm²
|      |      |      Dose rate µR/h, rem/h or 0.01µSv/h
|      | \----- Mean value
|      |      Count rate in 0.01 cps, cpm
|      |      Activity in 0.01 Bq, dps, dpm, Bq/cm²
|      |      Dose rate µR/h, rem/h or 0.01µSv/h
| \----- Date and Time as a decimal value
\----- History status as a decimal value
      see below
```

5.7.6.2 History readout scaler

256 716612088 721 999 120 23

```
-----
|      |      |      |      |      |
|      |      |      |      | \--- Temperature in °C
|      |      |      | \----- Measuring time in seconds
|      |      | \----- Background in 0.01 cps
|      | \----- Mean value
|      |      Count rate in 0.01 cps, cpm
|      |      Activity in 0.01 Bq, dps, dpm, Bq/cm²
|      |      Dose rate µR/h, rem/h or 0.01µSv/h
| \----- Date and Time as a decimal value
\----- History status as a decimal value
      see below
```

End of History:

End

5.7.6.3 *History status*

Decimal value converted in HEX:

Bit number			
0	Net value	0:No	1:Yes
1	Operation mode	0: Ratemeter	1:Scaler
2	Not used		
3	not used		
4	Background measurement	0:No	1:Yes
5	Scaler with	0: Preset Counts	1: Preset Time
6	Used filter		
7	0: No filter 1: Alpha blocker, 2:H*(10), 3:Hx		
8	Number of used nuclide		
9			
10			
11			
12			
12	0x00 Display unit cps	0x04 Display unit dpm	
13	0x01 Display unit cpm	0x05 Display unit Sv/h	
14	0x02 Display unit Bq	0x06 Display unit R/h	
15	0x03 Display unit dps	0x07 Display unit rem/h	
		0x08 Display unit Bq/cm ²	

5.7.7 Nuclide table

nRNumber Reading nuclide data.
Number: consecutive number
 Response: Nuclide data (see below)

nWNumberString Write nuclide data
Number: consecutive number
String: nuclide data
 e.g. nW02Sr-90 500 555 3333

nRA Reading number of stored nuclides.
 Response: value from 0...15

nWANumber Write the number of stored nuclides.
Number: value from 0...15

nRG Reading active nuclide.
 Response: value from 0...15

nWANumber Write the number of active nuclides.
Number: value from 0...15

Nuclide data:

Sr-90 500 555 3333

```

-----
|      |      |      |
|      |      |      \----- Factor for activity calculation with gamma filter
|      |      \----- Factor for activity calculation with alpha blocker
|      \----- Factor for activity calculation without filter
\----- Nuclide name. Up to 6 characters. Do not use space (0x20)

```

Activity calculation:

$F = 1 / E_{eff}$

$A_r = C_r * F$

F: Factor for Efficiency calculation

E_{eff} Efficiency for this nuclide

A_r Activity

C_r Measured count rate

For example:

Count rate is 67.3 cps without filter and efficiency for Sr-90 is

- 29% without filter.
- 25% with alpha blocker
- 10% with gamma filter

Factor F is $1/0.29=3.448$.

- $1/0.29 = 3.448$ without filter.
- $1/0.25 = 4.0$ with alpha blocker
- $1/0.10 = 10.0$ with gamma filter

To set the parameter, sent: nW00Sr-90 345 400 1000.

The activity is $67.3*3.45=232.18$ Bq

5.7.8 Event log

6656 520549251

```

-----
|          |
|          | \----- Date and time as a decimal value (see 3.2)
|          | \----- Event log as a decimal value
\-----

```

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	Not used
9	
10	
11	Sound 0: off 1:on
12	LED 0: off 1:on
13	Vibration alarm 0: off 1:on
14	Clear dose 0: no 1:yes
15	1: Alarm threshold changed
16	1: Dose rate, count rate or activity -alarm
17	1: Alarm dose
18	Not used
19	Not used
20	Dose rate, count rate or activity > alarm threshold 1
21	Dose rate, count rate or activity > alarm threshold 2
22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24	1: Scaler or Background parameter changed
25	Not used
26	Power off
27	Power on
28	Not used
29	Not used
30	Not used
31	Not used

5.7.9 Automatic sending

```

<STX>721  2 0 0 14 B20ER 1234 AB<ETX><CR><LF>
|         | | | | |         |         |         |
|         | | | | |         |         |         \- BCC in HEX
|         | | | | |         |         \----- Dose in µR/h
|         | | | | |         \----- information RadEye *)
|         | | | | \----- Status as Hex-value
|         | | | \----- Fixed to "0"
|         | | \----- Fixed to "0"
|         | \----- Fixed to "2"
|         \----- Dose rate in µR/h

```

*) Information RadEye:

B20	RadEye B20
B20ER	RadEye B20-ER
G20	RadEye G20
G20ER	RadEye G20-ER
G2010	RadEye G20-10
G20ER1	RadEye G20-ER10

Formation of a BCC (block check character):

Modulo 256 sum of <STX> up to the last character before the BCC (including), coded as hexadecimal ASCII-number (e.g. 1F).

Status

Bit number	
0	Not used
1	1: Overload
2	1: Alarm
3	Not used
4	Not used
5	1: Battery voltage low
6	not used
7	not used

5.7.10 Status information

F Reading status information
 Response Number with status information

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	Not used
9	
10	
11	not used,
12	Alarm threshold read-only 0: off 1:on
13	Flag for overload
14	Temperature display 0: off 1:on
15	Not used
16	1: Alarm
17	Not used
18	Not used
19	Not used
20	Value > alarm threshold 1 (depending on display mode)
21	Value > alarm threshold 2 (depending on display mode)
22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24	Not used
25	Not used
26	Not used
27	Not used
28	Not used
29	Not used
30	Not used
31	Not used

5.8 RadEye DW

Used firmware version: V1.53

5.8.1 Limit values

AR5	Reading threshold 1 for dose rate alarm. Response: number in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ unit. Depending on used measuring unit e.g. 1 2 3 means 1 2 3 $\mu\text{R/h}$ resp. 1.23 $\mu\text{Sv/h}$
AR7	Reading threshold 1 for dose alarm Response: number in μSv , 100 μR or 100 μrem units. 0...10.000.000
AR8	Reading threshold 2 for dose alarm Response: see command AR7
ARK	Reading threshold 3 for dose alarm Response: see command AR7
ARL	Reading threshold 4 for dose alarm Response: see command AR7
AR0	Reading initial threshold Response: value 0 or 1
ARP	Reading test pulses Response: number from 0 to 9999
aD	Read active threshold Response: number from 0 to 3, representing threshold 1...4
AW5Number	Setting the threshold 1 for dose rate alarm. <i>Number</i> = in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units. e. g. AW3123 means 123 $\mu\text{R/h}$.
AW7Number	Setting the threshold 1 for dose alarm. <i>Number</i> = in 100 μR , 100 μrem or 1 μSv units. e. g. AW7123 means 12.3mR.
AW8Number	Setting the threshold 2 for dose alarm. <i>Number</i> = see command AW7
AWKNumber	Set threshold 3 for dose alarm Response: see command AR7

<i>AWLNumber</i>	Set threshold 4 for dose alarm Response: see command AR7
<i>AW0Number</i>	Set initial threshold Response: value 0 or 1
<i>AWPNumber</i>	Set test pulses <i>Number</i> : value from 0 to 9999

5.8.2 Measurement values

<i>Z</i>	Read raw count rates with dead time correction Response: <ul style="list-style-type: none">- Counter 1 in cps- HV power in cps
----------	--

5.8.3 Configuration flags

5.8.3.1 Configuration flags 1 with kR / kW

Bit number	
0	Fixed to “0”
1	Fixed to “0”
2	Fixed to “1”
3	Alarming Sound 0: off 1: on
4	Alarming LED 0: off 1: on
5	Alarming Vibration 0: off 1: on
6	not used, write “0”
7	Single Pulse 0: off 1: on

5.8.3.2 Configuration flags 2 with fR / fW

Bit number	
0	Measuring unit: 0: Sievert 1: Roentgen
1	Fixed to “0”
2	Fixed to “1”
3	Not used
4	Flag for overload (readonly)
5	Temp.display 0: off 1:on
6	Temperature unit 0: °C 1: °F
7	Not used

5.8.3.3 *Configuration flags 3 with KR / KW*

Bit number	
0	Key tone 0:off 1:on
1	Autosend 0:off 1:on
2	Finder 0:off 1:on
3	Mean-Max 0:off 1:on
4	not used, write "0"
5	display of dose 0:off 1:on
6	Fixed to "1"
7	not used, write "0"

5.8.3.4 *Configuration flags 4 with jR / jW*

Bit number	
0	not used, write "0"
1	Battery type 0:Alkaline 1:NiMh
2	Display rotation 0:No 1:yes
3	Show positions after decimal point (from 1Sv) 0:No 1:yes
4	not used, write "0"
5	not used, write "0"
6	not used, write "0"
7	not used, write "0"
8	not used, write "0"
9	not used, write "0"
10	not used, write "0"
11	not used, write "0"
12	not used, write "0"
13	not used, write "0"
14	not used, write "0"
15	not used, write "0"

5.8.3.5 *Menu configuration*

Bit number			
0	Switch off	0:hidden	1:visible
1	Sound	0: hidden	1:visible
2	LED	0: hidden	1:visible
3	Vibrator	0: hidden	1:visible
4	not used, write "0"		
5	not used, write "0"		
6	not used, write "0"		
7	not used, write "0"		
8	not used, write "0"		
9	Change alarm dose	0: hidden	1:visible
10	Pre dose		
11	Autosend	0: hidden	1:visible
12	Test mode		
13	not used, write "0"		
14	not used, write "0"		
15	Single Pulse	0: hidden	1:visible
16	Backlight	0: hidden	1:visible
17	Show Alarm	0: hidden	1:visible
18	Settings	0: hidden	1:visible
19	Text Info	0: hidden	1:visible
20	not used, write "0"		
21	not used, write "0"		
22	not used, write "0"		
23	not used, write "0"		
24	not used, write "0"		
25	not used, write "0"		
26	not used, write "0"		
27	not used, write "0"		
28	not used, write "0"		
29	not used, write "0"		
30	not used, write "0"		
31	not used, write "0"		

5.8.4 **High voltage**

HR Reading high voltage bit value
 Response: Bit value 0...255.

5.8.5 Dead time correction

x Read dead time
 Response: dead time in ns

& Read dead time coefficient

5.8.6 History output

5.8.6.1 History readout

History status for 1. readout and change of history cycle time

```
1073741944 520534610
-----
|           |
|           |
|           | \----- Date and Time as a decimal value
| \----- History status as a decimal value
```

following readout:

```
15 11 +19
-- -- --
|   |   |
|   |   | \--- Temperature in °C
| \----- Rate mean in µR/h, µrem/h or 0.01µSv/h units
| \----- Rate max in µR/h, µrem/h or 0.01µSv/h units
```

End of History:

End

5.8.6.2 History status

decimal value converted in HEX:

```
0x40000078
|   ----
|   |
|   | \-- History cycle time (0x78 -> 120s)
|
| \----- 4: identifier for status information
|           6: identifier for status information +
|           last power on without power off
```

5.8.7 Event log

6656 520549251

```

|      |
|      \-----
\-----

```

Date and time as a decimal value (see 3.2)

Event log as a decimal value

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	Fixed to "0"
9	Fixed to "0"
10	Fixed to "1"
11	Sound 0: off 1:on
12	LED 0: off 1:on
13	Vibration alarm 0: off 1:on
14	1: Dose cleared
15	1: Alarm threshold changed
16	1: Alarm dose rate
17	1: Alarm dose
18	Not used
19	Not used
20	Dose rate > alarm threshold 1
21	Not used
22	Dose > alarm threshold 1
23	Not used
24	Not used
25	Not used
26	Power off
27	Power on
28	Pre dose enabled
29	Pre dose disabled
30	Not used
31	Fixed to "0"

5.8.8 Automatic sending

<STX>7 2 0 0 14 REGDW 123 AB<ETX><CR><LF>

```

| | | | | | | |
| | | | | | | \- BCC in HEX
| | | | | | | \----- Dose in µR/h
| | | | | | | \----- information RadEye DW
| | | | | | | \----- Status as Hex-value
| | | | | | | \----- „0“ not used
| | | | | | | \----- „0“ not used
| | | | | | | \----- „2“ not used
| | | | | | | \----- Dose rate in µR/h

```

Formation of a BCC (block check character):

Modulo 256 sum of <STX> up to the last character before the BCC (including), coded as hexadecimal ASCII-number (e.g. 1F).

Status

Bit number	
0	Not used
1	1: Overload
2	1: Alarm dose rate
3	1: Alarm dose
4	not used
5	1: Battery voltage low
6	not used
7	not used

5.8.9 Status information

F Reading status information
 Response Number with status information

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	Measuring unit: 0: Sievert 1: Roentgen
9	Fixed to “0”
10	Fixed to “1”
11	not used,
12	Alarm threshold read-only 0: off 1:on
13	Flag for overload
14	Temperature display 0: off 1:on
15	Fixed to “0”
16	1: Alarm dose rate
17	1: Alarm dose
18	Fixed to “0”
19	Not used
20	Dose rate > alarm threshold 1
21	Not used
22	Dose > alarm threshold 1
23	Not used
24	Not used
25	Not used
26	Not used
27	Not used
28	Not used
29	Not used
30	Not used
31	Not used

5.8.10 Miscellaneous

vRs Read preset dose
 Response: number from 1 to 250,000

vWs*Number* Set preset dose
 Number: value from 1 to 250,000

DT

Read time to alarm

Response: value in minutes.

5.9 RadEye DLW

Used firmware version: V1.53

5.9.1 Limit values

AR5	Reading threshold 1 for dose rate alarm. Response: number in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ unit. Depending on used measuring unit e.g. 1 2 3 means 1 2 3 $\mu\text{R/h}$ resp. 1.23 $\mu\text{Sv/h}$
ARPNumber	Read test pulses Response: number from 0 to 9999
AW5Number	Setting the threshold 1 for dose rate alarm. Number = in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units. e. g. AW3123 means 123 $\mu\text{R/h}$.
AWPNumber	Set test pulses Number: value from 0 to 9999

5.9.2 Measurement values

Z	Read raw count rates with dead time correction Response: <ul style="list-style-type: none">- Counter 1 in cps- HV power in cps
---	--

5.9.3 Configuration flags

5.9.3.1 Configuration flags 1 with kR / kW

Bit number	
0	Fixed to “0”
1	Fixed to “0”
2	Fixed to “1”
3	Alarming Sound 0: off 1: on
4	Alarming LED 0: off 1: on
5	Alarming Vibration 0: off 1: on
6	not used, write “0”
7	Single Pulse 0: off 1: on

5.9.3.2 Configuration flags 2 with fR / fW

Bit number	
0	Measuring unit: 0: Sievert 1: Roentgen
1	Fixed to “0”
2	Fixed to “1”
3	Not used
4	Flag for overload (read-only)
5	Temp.display 0: off 1:on
6	Temperature unit 0: °C 1: °F
7	Not used

5.9.3.3 *Configuration flags 3 with KR / KW*

Bit number	
0	Key tone 0:off 1:on
1	Autosend 0:off 1:on
2	Finder 0:off 1:on
3	Mean-Max 0:off 1:on
4	not used, write "0"
5	display of dose 0:off 1:on
6	Fixed to "1"
7	not used, write "0"

5.9.3.4 *Configuration flags 4 with jR / jW*

Bit number	
0	not used, write "0"
1	Battery type 0:Alkaline 1:NiMh
2	Display rotation 0:No 1:yes
3	not used, write "0"
4	not used, write "0"
5	not used, write "0"
6	not used, write "0"
7	not used, write "0"
8	not used, write "0"
9	not used, write "0"
10	not used, write "0"
11	not used, write "0"
12	not used, write "0"
13	not used, write "0"
14	not used, write "0"
15	not used, write "0"

5.9.3.5 *Menu configuration*

Bit number			
0	Switch off	0:hidden	1:visible
1	Sound	0: hidden	1:visible
2	LED	0: hidden	1:visible
3	Vibrator	0: hidden	1:visible
4	not used, write "0"		
5	not used, write "0"		
6	not used, write "0"		
7	not used, write "0"		
8	not used, write "0"		
9	Change alarm dose	0: hidden	1:visible
10	not used, write "0"		
11	Autosend	0: hidden	1:visible
12	Test mode		
13	not used, write "0"		
14	not used, write "0"		
15	Single Pulse	0: hidden	1:visible
16	Backlight	0: hidden	1:visible
17	Show Alarm	0: hidden	1:visible
18	Settings	0: hidden	1:visible
19	Text Info	0: hidden	1:visible
20	not used, write "0"		
21	not used, write "0"		
22	not used, write "0"		
23	not used, write "0"		
24	not used, write "0"		
25	not used, write "0"		
26	not used, write "0"		
27	not used, write "0"		
28	not used, write "0"		
29	not used, write "0"		
30	not used, write "0"		
31	not used, write "0"		

5.9.4 **High voltage**

HR Reading high voltage bit value
 Response: Bit value 0...255.

5.9.5 Dead time correction

x Read dead time
 Response: dead time in ns

5.9.6 History output

5.9.6.1 History readout

History status for 1. readout and change of history cycle time

```
1073741944 520534610
-----
|           |
|           |
|           | \----- Date and Time as a decimal value
| \----- History status as a decimal value
```

following readout:

```
15 11 +19
-- -- --
|   |   |
|   |   | \--- Temperature in °C
| \----- Rate mean in µR/h, µrem/h or 0.01µSv/h units
| \----- Rate max in µR/h, µrem/h or 0.01µSv/h units
```

End of History:

End

5.9.6.2 History status

decimal value converted in HEX:

```
0x40000078
|   ----
|   |
|   | \-- History cycle time (0x78 -> 120s)
|
| \----- 4: identifier for status information
|           6: identifier for status information +
|           last power on without power off
```

5.9.7 Event log

6656 520549251

```

|      |
|      \-----
\-----

```

Date and time as a decimal value (see 3.2)

Event log as a decimal value

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	Fixed to "0"
9	Fixed to "0"
10	Fixed to "1"
11	Sound 0: off 1:on
12	LED 0: off 1:on
13	Vibration alarm 0: off 1:on
14	1: Dose cleared
15	1: Alarm threshold changed
16	1: Alarm dose rate
17	1: Alarm dose
18	Not used
19	Not used
20	Dose rate > alarm threshold 1
21	Not used
22	Not used
23	Not used
24	Not used
25	Not used
26	Power off
27	Power on
28	Not used
29	Not used
30	Not used
31	Fixed to "0"

5.9.8 Automatic sending

<STX>7 2 0 0 14 REGDLW 123 AB<ETX><CR><LF>

17	1: Alarm dose
18	Fixed to “0”
19	Not used
20	Dose rate > alarm threshold 1
21	Not used
22	Not used
23	Not used
24	Not used
25	Not used
26	Not used
27	Not used
28	Not used
29	Not used
30	Not used
31	Not used

5.10 RadEye AB100 (from V1.67)

Used firmware version: 1.75

5.10.1 Limit values

AR0	Reading parametrs for sigma threshold. Response: <ul style="list-style-type: none">- Sigma value 0 (=off), 2...9.- Min. count rate for sigma alarm in cps- Max. background value in cps- Min. background value in cps- Actual sigma alarm threshold in 0.01 cps units
AR1	Reading the threshold 1 for activity alarm β threshold. Response: number in 0.01 Bq units.
AR2	Reading the threshold 2 for activity alarm β threshold. Response: number in 0.01 Bq units.
AR3	Reading the threshold 1 for activity alarm α threshold. Response: number in 0.01 Bq units.
AR4	Reading the threshold 2 for activity alarm α threshold. Response: number in 0.01 Bq units.
AR5	Reading the threshold 1 for count rate alarm β threshold. Response: number in 0.01 cps units
AR6	Reading the threshold 2 for count rate alarm β threshold. Response: number in 0.01 cps units
AR7	Reading the threshold 1 for count rate alarm α threshold. Response: number in 0.01 cps units
AR8	Reading the threshold 2 for count rate alarm α threshold. Response: number in 0.01 cps units
AR9	Reading alarm thresholds for dose rate and dose Response: <ul style="list-style-type: none">- threshold 1 for dose rate alarm in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units- threshold 2 for dose rate alarm in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units- threshold 1 for dose alarm in μR, μrem or $0.01\mu\text{Sv}$ units- threshold 2 for dose alarm in μR, μrem or $0.01\mu\text{Sv}$ units
ARP	Reading scaler parameter:

- Preset count threshold 1. Response: Number in counts
- Preset count threshold 3. Response: Number in counts
- Preset time. Response: Number in seconds

ARN

Reading background parameter:

- Preset count threshold 1. Response: Number in counts
- Preset count threshold 3. Response: Number in counts
- Preset time. Response: Number in seconds

ARB

Reading background values:

- Background value threshold 1. Response: Value in 0.01 cps
- Background value threshold 3. Response: Value in 0.01 cps

AW0Number Number... Writing parametrs for sigma threshold.

Number: Sigma value 0 (=off), 2...9.

Number: Min. count rate for sigma alarm in cps

Number: Max. background value in cps

Number: Min. background value in cps

AW1Number

Setting the threshold 1 for activity alarm β threshold.

Number: value in 0.01 Bq units.

AW2Number

Setting the threshold 2 for activity alarm β threshold

Number: value in 0.01 Bq units

AW3Number

Setting the threshold 1 for activity alarm α threshold

Number: value in 0.01 Bq units.

AW4Number

Setting the threshold 2 for activity alarm α threshold

Number: value in 0.01 Bq units

AW5Number

Setting the threshold 1 for count rate alarm β threshold..

Number: value in 0.01 cps units.

AW6Number

Setting the threshold 2 for count rate alarm β threshold.

Number: value in 0.01 cps units.

AW7Number

Setting the threshold 1 for count rate alarm α threshold.

Number: value in 0.01 cps units

AW8Number

Setting the threshold 2 for count rate alarm α threshold.

Number: value in 0.01 cps units

AW9 Number Number... Setting alarm thresholds for dose rate and dose

Number: threshold 1 for dose rate alarm in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units

Number: threshold 2 for dose rate alarm in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units
Number: threshold 1 for dose alarm in μR , μrem or $0.01\mu\text{Sv}$ units
Number: threshold 2 for dose alarm in μR , μrem or $0.01\mu\text{Sv}$ units

AWPNumber Number.. Set scaler parameter:

Number: Preset count threshold 1 in counts (from 0 to 65000 counts)
Number: Preset count threshold 3 in counts (from 0 to 65000 counts)
Number: Preset time in seconds (from 0 to 9999 seconds)

AWNNumber Number.. Set background parameter:

Number: Preset count threshold 1 in counts (from 0 to 65000 counts)
Number: Preset count threshold 3 in counts (from 0 to 65000 counts)
Number: Preset time in seconds (from 0 to 9999 seconds)

AWBNumber Number Set background value:

Number: Background value threshold 1 in 0.01 cps (from 0 to 65000 cps)
Number: Background value threshold 3 in 0.01 cps (from 0 to 65000 cps)

5.10.2 Measurement values

Z Read raw count rates with dead time correction

Response:

- Counter 1 in cps
- Counter 2 in cps
- Counter 3 in cps
- Counter 4 in cps
- HV power index
- Probe current in 0.1 μA units

z Read filtered count rate β and α threshold

Response:

- Value β threshold 0.01 cps
- Value α threshold in 0.01 cps

A Read display value and status

Response:

- display value β threshold in 0.01 cps, cpm, Bq, dps, dpm or Bq/cm² units
- display value α threshold in 0.01 cps, cpm, Bq, dps dpm or Bq/cm² units or in μR , μrem or $0.01\mu\text{Sv}$ units in dose rate mode
- Status (see 5.5.6.2)

5.10.3 Scaler remote control (from V1.73)

CG Start scaler with parameter previously set

CS Stop scaler

CI

Scaler information

Response:

- 0: Scaler stopped, 1: Scaler active
- Scaler counting time
- accumulated counts threshold 1 or window 1
- accumulated counts threshold 3 or window 2
- accumulated counts threshold 1
- accumulated counts threshold 2
- accumulated counts threshold 3
- accumulated counts threshold 4

5.10.4 Configuration flags

5.10.4.1 Configuration flags 1 with *kR* / *kW*

Bit number	
0	Keylock 0: enable 1: disable
1	Not used
2	Not used
3	Alarming Sound 0: off 1: on
4	Alarming LED 0: off 1: on
5	Alarming Vibration 0: off 1: on
6	Not used
7	Single Pulse 0: off 1: on

5.10.4.2 Configuration flags 2 with *fR* / *fW*

Bit number	
0	Show Gamma Symbol 0: no 1:yes
1	Not used
2	Not used
3	Alarm threshold read-only 0: off 1:on
4	Flag for overload (read-only)
5	Display of temperature 0: off 1:on
6	Temperature unit 0: °C 1: °F
7	Not used

5.10.4.3 Configuration flags 3 with KR / KW

Bit number	
0	Key tone 0:off 1:on
1	Autosend 0:off 1:on
2	Finder 0:off 1:on
3	Mean-Max 0:off 1:on
4	not used, write "0"
5	not used
6	not used
7	not used

5.10.4.4 Configuration flags 4 with jR / jW

Bit number	
0	Not used
1	Battery type 0:Alkaline 1:NiMh
2	Display rotation 0:No 1:yes
3	Scaler mode 0: Preset count 1: Preset time
4	Accumulated counts 0:no 1:yes
5	Scaler, after measurement 0: Stop 1: automatic restart
6	Scaler net 0:No 1:yes
7	Ratemeter net 0:No 1:yes
8	Alpha-LED 0:No 1:yes
9	Alpha-Sound 0:No 1:yes
10	Not used
11	Alarm acknowledge (from 1.75) 0:No 1:yes
12	Set Alarm1 to 3xBackground (from 1.75) 0:no 1:yes
13	Ratemeter Tau active 0:No 1:yes
14	Not used
15	Not used

5.10.4.5 *Measuring unit with uR /uW*

uR Read measuring unit and operation mode

Response: Hex-value. See below

uRHex Write measuring unit and operation mode

See below

Bit number	
0	0x00 Display unit cps 0x04 Display unit dpm
1	0x01 Display unit cpm 0x08 Display unit Bq/cm ²
2	0x02 Display unit Bq
3	0x03 Display unit dps
4	0: Ratemeter (Tau or ADF), 1: Scaler
5	Display mode:
6	00: measurement value β threshold 01: measurement value α threshold 02: dual display 03: graphic display
7	0: two decimal places 1: one decimal place

5.10.4.6 Menu configuration

mR Read configuration for main menu and submenu “Settings”
Response: Hex-values. See below

mRHex Hex Write configuration for main menu and submenu “Settings”
See below

5.10.4.6.1 Main menu

Bit number			
0	Switch off	0:hidden	1:visible
1	Background	0: hidden	1:visible
2	Select counter tube	0: hidden	1:visible
3	Backlight	0: hidden	1:visible
4	Measuring unit	0: hidden	1:visible
5	Operation mode	0: hidden	1:visible
6	Scaler parameter	0: hidden	1:visible
7	Nuclide table	0: hidden	1:visible
8	Alarm β Alarm act. β Alarm dose rate	0: hidden	1:visible
9	Alarm α Alarm act. α Alarm dose	0: hidden	1:visible
10	not used		
11	not used		
12	Clear dose	0: hidden	1:visible
13	Settings	0: hidden	1:visible
14	Alarm indication	0: hidden	1:visible
15	Show alarm	0: hidden	1:visible
16	Text Info	0: hidden	1:visible
17	not used, write “0”		
18	not used, write “0”		
19	not used, write “0”		
20	not used, write “0”		
21	not used, write “0”		
22	not used, write “0”		
23	not used, write “0”		
24	not used, write “0”		
25	not used, write “0”		
26	not used, write “0”		
27	not used, write “0”		
28	not used, write “0”		
29	not used, write “0”		
30	not used, write “0”		
31	not used, write “0”		

5.10.4.6.2 Submenu “Settings”

Bit number			
0	Batt. type	0:hidden	1:visible
1	Autosend	0: hidden	1:visible
2	Single Pulse	0: hidden	1:visible
3	Finder	0: hidden	1:visible
4	Alpha-LED	0: hidden	1:visible
5	Alpha Sound	0: hidden	1:visible
6	Set Date/Time	0: hidden	1:visible
7	Display $\alpha\beta$ / β	0: hidden	1:visible
8	Not used		
9	Set HV	0: hidden	1:visible
10	Not used		
11	Language (from V1.75)	0: hidden	1:visible
12	Edit Tau (from V1.75)	0: hidden	1:visible
13	Contrast (from V1.75)	0: hidden	1:visible
14, 15	not used, write “0”		

5.10.4.6.3 Submenu “Alarm indication”

Bit number			
0	Sound	0:hidden	1:visible
1	LED	0: hidden	1:visible
2	Vibrator	0: hidden	1:visible
3...8	not used, write “0”		

5.10.4.6.4 Submenu “Operation mode”

Bit number			
0	Ratemeter ADF	0:hidden	1:visible
1	Scaler	0: hidden	1:visible
2	Ratemeter Tau	0: hidden	1:visible

5.10.5 High voltage

HR Reading high voltage
 Response: value in Volt.

5.10.6 History output

5.10.6.1 History readout ratemeter

```
256 716612088 721 999 12 50 120 23
-----
|      |      |      |      |      |      |      |
|      |      |      |      |      |      | \---- Temperature in °C
|      |      |      |      |      | \----- Measuring time in seconds
|      |      |      |      | \----- Max value α channel
|      |      |      |      |      Count rate in 0.01 cps
|      |      |      |      |      Activity in 0.01 Bq
|      |      |      |      |      Dose rate in 1 μR/h
|      |      |      | \----- Mean value α channel
|      |      |      |      count rate in 0.01 cps,
|      |      |      |      Activity in 0.01 Bq
|      |      |      |      Dose rate in 1 μR/h
|      |      | \----- Max value β channel
|      |      |      count rate in 0.01 cps,
|      |      |      Activity in 0.01 Bq
|      | \----- Mean value β channel
|      |      count rate in 0.01 cps,
|      |      Activity in 0.01 Bq
| \----- Date and Time as a decimal value
\----- History status as a decimal value
see below
```

5.10.6.2 History readout scaler

```
256 716612088 721 999 12 50 120 23
-----
|      |      |      |      |      |      |      |
|      |      |      |      |      |      | \---- Temperature in °C
|      |      |      |      |      | \----- Measuring time in seconds
|      |      |      |      | \----- Background value α channel
|      |      |      |      |      Count rate in 0.01 cps
|      |      |      | \----- Measuring value α channel
|      |      |      |      count rate in 0.01 cps,
|      |      |      |      Activity in 0.01 Bq
|      |      |      |      Dose rate in 1 μR/h
|      |      | \----- Background value β channel
|      |      |      count rate in 0.01 cps,
|      | \----- Measuring value β channel
|      |      count rate in 0.01 cps,
|      |      Activity in 0.01 Bq
| \----- Date and Time as a decimal value
\----- History status as a decimal value
see below
```

End of History:

End

5.10.6.3 *History status*

Decimal value converted in HEX:

Bit number			
0	Net value	0:No	1:Yes
1	Operation mode	0: Ratemeter	1:Scaler
2	Number of used probe		
3			
4			
5			
6	Background measurement	0:No	1:Yes
7	Display of	0: $\alpha\beta$	1: β
8	Number of used nuclide		
9			
10			
11			
12	0x00 Display unit cps 0x08 Display unit Bq/cm ² 0x01 Display unit cpm 0x02 Display unit Bq 0x03 Display unit dps 0x04 Display unit dpm		
13			
14			
15			
16	Probe style		
17	000: $\alpha\beta$ -Probe		
18	001: α -Probe		
	010: β -Probe		
19	Not used		
20	Accumulated counts	0:No	1:Yes
21	Not used		
22	Not used		
23	Not used		
24	Not used		
25	Not used		
26	Not used		
27	Not used		
28	Not used		
29	Not used		
30	Not used		
31	Not used		

5.10.7 Nuclide table

nRNumber Reading nuclide data.
Number: number of probe
 Response: Nuclide data (see below)

nWNumberString Write nuclide data
Number: number of probe
String: nuclide data (see below)

Nuclide data: every string contains up to 16 nuclide data of the corresponding probe

For example:

00 01 Sr-90 471 999 Am-241 330 123 usw. 0815

\--	Checksum (in HEX)
\-----	Factor for α threshold
\-----	Factor for β threshold
\-----	Name of nuclide 2 (do not use space)
\-----	Factor for α threshold
\-----	Factor for β threshold
\-----	Name of nuclide 1 (do not use space)
\-----	Active nuclide number
\-----	Number of nuclides in this string

Activity calculation:

$$F=1/E_{eff}$$

$$Ar=Cr \cdot F$$

For example:

Count rate is 67.3 cps and efficiency for Sr-90 is 29%. Factor F is $1/0.29=3.448$.

The activity is $67.3 \cdot 3.45=232.18$ Bq.

Write "0" for null efficiency. In this case the RadEye displays "### Bq".

Checksum calculation: Checksum starts from the beginning of the message to the end, with the checksum. Checksum is modulo 256.

5.10.8 Probe data

<i>zRNumber</i>	Reading probe data. <i>Number</i> : consecutive number Response: probe data (see below)
<i>zWNumberString</i>	Write probe data <i>Number</i> : consecutive number <i>String</i> : probe data
<i>zRA</i>	Read the number of stored probes. Response: value from 0...15
<i>zWANumber</i>	Read the number of stored probes. <i>Number</i> : value from 0...15
<i>zRG</i>	Reading active probe. Response: value from 0...15
<i>zWGNumber</i>	Write active probe. <i>Number</i> : value from 0...15

For example:

HP-380AB 10000 5 5 5 5 600 100 100 30 100 10 1 3 5 8 30 225 1000 1250 1234

1 2 3 4 5 6 7 8 9 A B C D E F G H I J K L

1	Probe name. Max. 11 characters. Do not use space
2	Calibration factor in (nSv/h)/cps
3	Dead time Rate 1 μ s
4	Dead time Rate 2 μ s
5	Dead time Rate 3 μ s
6	Dead time Rate 4 μ s
7	High voltage in Volt
8	Overload Rate 1 in kcps
9	Overload Rate 3 in kcps
A	Overload probe current μ A
B	Area of this probe
C	Timeout for detector error in seconds
D	Factor Rate 1 for dose rate calculation. Range: -128...+127
E	Factor Rate 2 for dose rate calculation. Range: -128...+127
F	Factor Rate 3 for dose rate calculation. Range: -128...+127

G
H
I
J
K
L

Flags

Bit number		
0	Window:	
1	00: No window. Display of Rate 1 Rate 3 01: β window. Display Rate 1 – Rate 2 and Rate 3 10: α window. Display of Rate 1 and Rate 3 – Rate 4 11: β and α window. Display of Rate 1 – Rate 2 and Rate 3 – Rate 4	
2	Pulse fade out 0: off 1: on	
3	Probe style	
4	000: $\alpha\beta$ -Probe	011: Dose rate probe
5	001: α -Probe	100: Gamma probe multi channel
	010: β -Probe	
6	Reserved, write “0”	
7	Reserved, write “0”	
8	Single Pulse 0: every Count, 1: with Divide Ratio (see Command ARC)	
9	Not used	
10	Display Rate 4 in place of Rate 3 (from V2.00)	
11..15	Not used	

Dead time correction:

$$C1 = \frac{C1_{raw}}{1 - \tau1 * C1_{raw}}$$

τ_1	Dead time
----------	-----------

C1raw: raw count rate

1- $\tau_1 \cdot C_{1\text{raw}}$ is limited to min. 0.1

Dose rate calculation:

$$DR = Cf * (\boxed{D} * Rate1 + \boxed{E} * Rate2 + \boxed{F} * Rate3 + \boxed{G} * Rate4)$$

DR: Dose rate in nSv/h

Cf: Calibration factor in nSv/h/cps

Rate 1..4 Count rate of threshold 1..4 in cps

5.10.9 Event log

6656 520549251

```

-----
|          |
|          \----- Date and time as a decimal value (see 3.2)
\----- Event log as a decimal value

```

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Probe Type Part 1
7	
8	Operation mode 0: Ratemeter 1:Scaler
9	Probe Type Part 2
10	
11	Sound 0: off 1:on
12	LED 0: off 1:on
13	Vibration alarm 0: off 1:on
14	Not used
15	1: Mode activity
16	1: count rate or activity -alarm
17	Not used
18	Not used
19	
20	Count rate or activity > alarm threshold 1
21	Count rate or activity > alarm threshold 2
22	Not used
23	Not used
24	Not used
25	Not used
26	Power off
27	Power on
28	Not used
29	Not used
30	1: Window Beta
31	Not used

5.10.10 Automatic sending

```

<STX>721  5 0 0 14 SX 1234 AB<ETX><CR><LF>
|         | | | | | | | |
|         | | | | | | | \---- BCC in HEX
|         | | | | | | | \----- Dose in µR/h
|         | | | | | | \----- information RadEye GX
|         | | | | | \----- Status as Hex-value
|         | | | \----- Fixed to "0"
|         | | \----- Fixed to "0"
|         | \----- Fixed to "5"
|         \----- Count rate 1 in 0.01cps

```

Formation of a BCC (block check character):

Modulo 256 sum of <STX> up to the last character before the BCC (including), coded as hexadecimal ASCII-number (e.g. 1F).

Status

Bit number	
0	Not used
1	1: Overload
2	1: Alarm
3	Not used
4	Not used
5	1: Battery voltage low
6	not used
7	not used

5.10.11 Status information

F Reading status information
 Response Number with status information

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	Not used
9	
10	
11	not used,
12	Alarm threshold (read-only) 0: off 1:on
13	Flag for overload
14	Temperature display 0: off 1:on
15	Not used
16	1: Alarm
17	Not used
18	Not used
19	Not used
20	Value > alarm threshold 1 (depending on display mode)
21	Value > alarm threshold 2 (depending on display mode)
22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24	Not used
25	Not used
26	Not used
27	Not used
28	Not used
29	Not used
30	Not used
31	Not used

5.10.12 Miscellaneous

ARC Read single pulse divider
 Response: number for single pulse divider

ARt	Read fixed time constant (from V1.75) Response: time in 0.1s units
<i>AWCNumber</i>	Write single pulse divider Number: value from 1 to 65535
AWt	Write fixed time constant (from V1.75) Number: value from 10 to 65535

5.11 RadEye GX

Used firmware version: 1.62

5.11.1 Limit values

AR0	Reading sigma parameter: <ul style="list-style-type: none">- Sigma value (0...9)- min. count rate for sigma alarm (cps)- max. Background (cps)- min. Background (cps)
AR1	Reading the threshold 1 for activity alarm. Response: number in 0.01 Bq units.
AR2	Reading the threshold 2 for activity alarm. Response: number in 0.01 Bq units.
AR3	Reading the threshold 1 for count rate alarm. Response: number in 0.01 cps units.
AR4	Reading the threshold 2 for count rate alarm. Response: number in 0.01 cps units.
AR5	Reading the threshold 1 dose rate alarm. Response: number in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units
AR6	Reading the threshold 2 for count rate alarm. Response: number in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units
AR7	Reading the threshold 1 for dose alarm. Response: number in μR , μrem or $0.01\mu\text{Sv}$ units
AR8	Reading the threshold 2 for dose alarm. Response number in μR , μrem or $0.01\mu\text{Sv}$ units
ARP	Reading scaler parameter: <ul style="list-style-type: none">- Preset count. Response: Number in counts- Preset time. Response: Number in seconds
ARN	Reading background parameter: <ul style="list-style-type: none">- Preset count. Response: Number in counts- Preset time. Response: Number in seconds

ARB	Reading background value: Response: Number in 0.01 cps
AW0	Set sigma parameter: - Sigma value (0...9) - min. count rate for sigma alarm - max. Background - min. Background
AW1Number	Setting the threshold 1 for activity alarm Number: value in 0.01 Bq units.
AW2Number	Setting the threshold 2 for activity alarm Number: value in 0.01 Bq units
AW3Number	Setting the threshold 1 for count rate alarm Number: value in 0.01 cps units.
AW4Number	Setting the threshold 2 for count rate alarm Number: value in 0.01 cps units
AW5Number	Setting the threshold 1 dose rate alarm. Number: value in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units.
AW6Number	Setting the threshold 2 dose rate alarm. Number: value in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units.
AW7Number	Setting the threshold 1 for dose alarm. Number: value in μR , μrem or $0.01\mu\text{Sv}$ units
AW8Number	Setting the threshold 2 for dose alarm. Number: value in μR , μrem or $0.01\mu\text{Sv}$ units
AWPNumber Number	Set scaler parameter: Number: Preset count in counts (from 0 to 9999 counts) Number: Preset time in seconds (from 0 to 9999 seconds)
AWNNumber Number	Set background parameter: Number: Preset count in counts (from 0 to 9999 counts) Number: Preset time in seconds (from 0 to 9999 seconds)
AWBNumber	Set background value: Number: Value in 0.01 cps (from 0 to 100.00 cps)

5.11.2 Measurement values

Z Read raw count rates with dead time correction

Response:

- Counter 1 in cps
- HV power index
- Counter 1 in cps without dead time
- HV power index

A Read display value and status

Response:

- display value in $1\mu\text{R/h}$, $1\mu\text{rem/h}$, $0.01\mu\text{Sv/h}$ units
or 0.01 cps, cpm, Bq, dps, dpm, Bq/cm² units
- Status (see 5.5.6.2)

5.11.3 Scaler remote control (from V1.66)

CG Start scaler with parameter previously set

CS Stop scaler

CI Scaler information

Response:

- 0: Scaler stopped, 1: Scaler active
- Scaler counting time
- accumulated counts

5.11.4 Configuration flags

5.11.4.1 Configuration flags 1 with kR / kW

Bit number	
0	Not used
1	Not used
2	Not used
3	Alarming Sound 0: off 1: on
4	Alarming LED 0: off 1: on
5	Alarming Vibration 0: off 1: on
6	Not used
7	Single Pulse 0: off 1: on

5.11.4.2 Configuration flags 2 with fR / fW

Bit number	
0	Show Gamma symbol 0: off 1: on
1	Not used

2	Not used
3	Alarm threshold read-only 0: off 1:on
4	Flag for overload (read-only)
5	Display of temperature 0: off 1:on
6	Temperature unit 0: °C 1: °F
7	Not used

5.11.4.5 *Measuring unit with uR /uW*

uR Read measuring unit and operation mode
 Response: Hex-value. See below
uRHex Write measuring unit and operation mode
 See below

Bit number	
0	0x00 Display unit cps 0x05 Display unit Sv/h
1	0x01 Display unit cpm 0x06 Display unit R/h
2	0x02 Display unit Bq 0x07 Display unit rem/h
3	0x03 Display unit dps 0x08 Display unit Bq/cm ² 0x04 Display unit dpm 0x09 Display unit Gy/h
4	0: Ratemeter, 1: Scaler
5	Not used
6	Not used
7	Not used

5.11.4.6 Menu configuration

mR	Read configuration for main menu and submenu “Settings” Response: Hex-values. See below
mRHex Hex	Write configuration for main menu and submenu “Settings” See below
mRHex HexHex	Write configuration for main menu, submenu “Settings” and submenu “Alarm indication” (from V1.64) See below

Bit number			
0	Switch off	0:hidden	1:visible
1	Background	0: hidden	1:visible
2	Select counter tube	0: hidden	1:visible
3	Backlight	0: hidden	1:visible
4	Measuring unit	0: hidden	1:visible
5	Operation mode	0: hidden	1:visible
6	Scaler parameter	0: hidden	1:visible
7	Nuclide table	0: hidden	1:visible
8	Alarm count rate	0: hidden	1:visible
9	Alarm dose rate	0: hidden	1:visible
10	Alarm activity	0: hidden	1:visible
11	Alarm dose	0: hidden	1:visible
12	Clear dose	0: hidden	1:visible
13	Settings	0: hidden	1:visible
14	Alarm indication	0: hidden	1:visible
15	Show alarm	0: hidden	1:visible
16	Text Info	0: hidden	1:visible
17	not used, write “0”		
18	not used, write “0”		
19	not used, write “0”		
20	not used, write “0”		
21	not used, write “0”		
22	not used, write “0”		
23	not used, write “0”		
24	not used, write “0”		
25	not used, write “0”		
26	not used, write “0”		
27	not used, write “0”		
28	not used, write “0”		
29	not used, write “0”		
30	not used, write “0”		
31	not used, write “0”		

Submenu “Settings”

Bit number			
0	Batt. type	0:hidden	1:visible

1	Autosend	0: hidden	1:visible
2	Single Pulse	0: hidden	1:visible
3	Finder	0: hidden	1:visible
4	Set Date/Time	0: hidden	1:visible
5	Set HV	0: hidden	1:visible
6	Bluetooth	0: hidden	1:visible
7	not used, write "0"		
8	not used, write "0"		
9	not used, write "0"		
10	not used, write "0"		
11	not used, write "0"		
12	not used, write "0"		
13	not used, write "0"		
14	not used, write "0"		
15	not used, write "0"		

Submenu "Alarm indication"

Bit number			
0	Sound	0:hidden	1:visible
1	LED	0: hidden	1:visible
2	Vibrator	0: hidden	1:visible
3	not used, write "0"		
4	not used, write "0"		
5	not used, write "0"		
6	not used, write "0"		
7	not used, write "0"		
8	not used, write "0"		

5.11.5 High voltage

HR Reading high voltage bit value
 Response: Bit value 0...255.

5.11.6 Dead time correction

x Read dead time
 Response: dead time in ns

5.11.7 History output

5.11.7.1 History readout ratemeter

256 716612088 721 999 120 23

```
---  -----  ---  ---  ---  ---
|      |      |      |      |      |
|      |      |      |      |      | \--- Temperature in °C
|      |      |      |      |      | \----- Measuring time in seconds
|      |      |      |      |      | \----- Max value
|      |      |      |      |      | Count rate in 0.01 cps, cpm
|      |      |      |      |      | Activity in 0.01 Bq, dps, dpm, Bq/cm²
|      |      |      |      |      | Dose rate µR/h, rem/h or 0.01µSv/h
|      |      |      |      |      | \----- Mean value
|      |      |      |      |      | Count rate in 0.01 cps, cpm
|      |      |      |      |      | Activity in 0.01 Bq, dps, dpm, Bq/cm²
|      |      |      |      |      | Dose rate 0.01 µR/h, rem/h or 0.1nSv/h
|      |      |      |      |      | \----- Date and Time as a decimal value
|      |      |      |      |      | \----- History status as a decimal value
|      |      |      |      |      | see below
```

5.11.7.2 History readout scaler

256 716612088 721 999 120 23

```
---  -----  ---  ---  ---  ---
|      |      |      |      |      |
|      |      |      |      |      | \--- Temperature in °C
|      |      |      |      |      | \----- Measuring time in seconds
|      |      |      |      |      | \----- Background in 0.01 cps
|      |      |      |      |      | \----- Mean value
|      |      |      |      |      | Count rate in 0.01 cps, cpm
|      |      |      |      |      | Activity in 0.01 Bq, dps, dpm, Bq/cm²
|      |      |      |      |      | Dose rate µR/h, rem/h or 0.01µSv/h
|      |      |      |      |      | \----- Date and Time as a decimal value
|      |      |      |      |      | \----- History status as a decimal value
|      |      |      |      |      | see below
```


End of History:

End

5.11.7.3 *History status*

Decimal value converted in HEX:

Bit number			
0	Net value	0:No	1:Yes
1	Operation mode	0: Ratemeter	1:Scaler
2	Number of used probe		
3			
4			
5			
6	Background measurement	0:No	1:Yes
7	Accumulated counts	0:No	1:Yes
8	Number of used nuclide		
9			
10			
11			
12	0x00 Display unit cps	0x04 Display unit dpm	
13	0x01 Display unit cpm	0x05 Display unit Sv/h	
14	0x02 Display unit Bq	0x06 Display unit R/h	
15	0x03 Display unit dps	0x07 Display unit rem/h	
		0x08 Display unit Bq/cm ²	

5.11.8 Nuclide table

nRNumber Reading nuclide data.
Number: consecutive number
 Response: Nuclide data (see below)

nWNumberString Write nuclide data
Number: consecutive number
String: nuclide data
 e.g. nW02Sr-90 321

nRA Reading number of stored nuclides.
 Response: value from 0...15

nWANumber Write the number of stored nuclides.
Number: value from 0...15

nRG Reading active nuclide.
 Response: value from 0...15

nWANumber Write the number of active nuclides.
Number: value from 0...15

Nuclide data:

For example:

Sr-90 321

| |

| \----- Factor for activity calculation with factor 100

\----- Nuclide name. Up to 6 characters. Do not use space (0x20)

Activity calculation:

$F = 1/E_{eff}$

$Ar = Cr * F$

For example:

Count rate is 67.3 cps and efficiency for Sr-90 is 29%. Factor F is $1/0.29=3.448$. To set the parameter, sent: nW00Sr-90 345.

The activity is $67.3*3.45=232.18$ Bq.

5.11.9 Probe data

<i>zRNumber</i>	Reading probe data. <i>Number</i> : consecutive number Response: probe data (see below)
<i>zWNumberString</i>	Write probe data <i>Number</i> : consecutive number <i>String</i> : probe data
<i>zRA</i>	Read the number of stored probes. Response: value from 0...15
<i>zWNumber</i>	Read the number of stored probes. <i>Number</i> : value from 0...15
<i>zRG</i>	Reading active probe. Response: value from 0...15
<i>zWGNumber</i>	Write active probe. <i>Number</i> : value from 0...15

For example:

HP-270 510 90 900 10 600 20 60 -40 00

1 2 3 4 5 6 7 8 9 10

1	Probe name. Max. 11 characters. Do not use space
2	Calibration factor in (nSv/h)/cps
3	Dead time μ s
4	High voltage in Volt
5	Overload counter in kcps
6	Overload high voltage power indicator
7	Area of this probe
8	Timeout for detector error in seconds
9	Dead time correction in 0.1 μ s
10	Flags

Flags

Bit number	
0	Not used
1	Not used
2	Not used
3	Not used
4	Not used
5	Not used
6	Not used
7	Reserved

Dead time correction:

$$C1 = \frac{C1_{raw}}{1 - \tau1 * C1_{raw} - \tau2 * C1_{raw}^2}$$

$\tau1$ Dead time

$\tau2$ Dead time correction

$C1_{raw}$: raw count rate

$\tau2 * C1_{raw}^2$ is limited to min. -0.8

$1 - \tau1 * C1_{raw} - \tau2 * C1_{raw}^2$ limited to min. 0.15

5.11.10 Event log

6656 520549251

| |
| \-----
\-----

Date and time as a decimal value (see 3.2)

Event log as a decimal value

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	Operation mode 0: Ratemeter 1:Scaler
9	Not used
10	
11	Sound 0: off 1:on
12	LED 0: off 1:on
13	Vibration alarm 0: off 1:on
14	Clear dose 0: no 1:yes
15	1: Alarm threshold changed
16	1: Dose rate, count rate or activity -alarm
17	1: Alarm dose
18	1: Overload
19	Not used
20	Dose rate, count rate or activity > alarm threshold 1
21	Dose rate, count rate or activity > alarm threshold 2
22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24	1: Scaler or Background parameter changed
25	Not used
26	Power off
27	Power on
28	Not used
29	Not used
30	Not used
31	Not used

5.11.11 Automatic sending

```

<STX>721  5 0 0 14 GX 1234 AB<ETX><CR><LF>
|         | | | | | | | |
|         | | | | | | | \---- BCC in HEX
|         | | | | | | | \----- Dose in µR/h
|         | | | | | | \----- information RadEye GX
|         | | | | | \----- Status as Hex-value
|         | | | \----- Fixed to "0"
|         | \----- Fixed to "0"
|         \----- Fixed to "5"
\----- Count rate in 0.01cps

```

Formation of a BCC (block check character):

Modulo 256 sum of <STX> up to the last character before the BCC (including), coded as hexadecimal ASCII-number (e.g. 1F).

Status

Bit number	
0	Not used
1	1: Overload
2	1: Alarm
3	Not used
4	Not used
5	1: Battery voltage low
6	not used
7	not used

5.11.12 Status information

F Reading status information
 Response Number with status information

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	Not used
9	
10	
11	not used,
12	Alarm threshold (read-only) 0: off 1:on
13	Flag for overload
14	Temperature display 0: off 1:on
15	Not used
16	1: Alarm
17	Not used
18	Not used
19	Not used
20	Value > alarm threshold 1 (depending on display mode)
21	Value > alarm threshold 2 (depending on display mode)
22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24	Not used
25	Not used
26	Not used
27	Not used
28	Not used
29	Not used
30	Not used
31	Not used

5.12 RadEye SX

Used firmware version: 2.00

5.12.1 Limit values

AR0	Reading parameters for sigma threshold. Response: <ul style="list-style-type: none">- Sigma value 0 (=off), 2...9.- Min. count rate for sigma alarm in cps- Max. background value in cps- Min. background value in cps- Actual sigma alarm threshold in 0.01 cps units
AR1	Reading the threshold 1 for activity alarm β threshold. Response: number in 0.01 Bq units.
AR2	Reading the threshold 2 for activity alarm β threshold. Response: number in 0.01 Bq units.
AR3	Reading the threshold 1 for activity alarm α threshold. Response: number in 0.01 Bq units.
AR4	Reading the threshold 2 for activity alarm α threshold. Response: number in 0.01 Bq units.
AR5	Reading the threshold 1 for count rate alarm β threshold. Response: number in 0.01 cps units
AR6	Reading the threshold 2 for count rate alarm β threshold. Response: number in 0.01 cps units
AR7	Reading the threshold 1 for count rate alarm α threshold. Response: number in 0.01 cps units
AR8	Reading the threshold 2 for count rate alarm α threshold. Response: number in 0.01 cps units
AR9	Reading alarm thresholds for dose rate and dose Response: <ul style="list-style-type: none">- threshold 1 for dose rate alarm in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units- threshold 2 for dose rate alarm in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units- threshold 1 for dose alarm in μR, μrem or $0.01\mu\text{Sv}$ units- threshold 2 for dose alarm in μR, μrem or $0.01\mu\text{Sv}$ units
ARP	Reading scaler parameter:

- Preset count threshold 1. Response: Number in counts
- Preset count threshold 3. Response: Number in counts
- Preset time. Response: Number in seconds

ARN

Reading background parameter:

- Preset count threshold 1. Response: Number in counts
- Preset count threshold 3. Response: Number in counts
- Preset time. Response: Number in seconds

ARB

Reading background values:

- Background value threshold 1. Response: Value in 0.01 cps
- Background value threshold 3. Response: Value in 0.01 cps

ARt

Reading Tau (from V2.00):

Response: Value in 0.1 s

AW0Number Number... Writing parametrs for sigma threshold.

Number: Sigma value 0 (=off), 2...9.

Number: Min. count rate for sigma alarm in cps

Number: Max. background value in cps

Number: Min. background value in cps

AW1Number

Setting the threshold 1 for activity alarm β threshold.

Number: value in 0.01 Bq units.

AW2Number

Setting the threshold 2 for activity alarm β threshold

Number: value in 0.01 Bq units

AW3Number

Setting the threshold 1 for activity alarm α threshold

Number: value in 0.01 Bq units.

AW4Number

Setting the threshold 2 for activity alarm α threshold

Number: value in 0.01 Bq units

AW5Number

Setting the threshold 1 for count rate alarm β threshold..

Number: value in 0.01 cps units.

AW6Number

Setting the threshold 2 for count rate alarm β threshold.

Number: value in 0.01 cps units.

AW7Number

Setting the threshold 1 for count rate alarm α threshold.

Number: value in 0.01 cps units

AW8Number

Setting the threshold 2 for count rate alarm α threshold.

Number: value in 0.01 cps units

AW9 *Number Number...* Setting alarm thresholds for dose rate and dose

Number: threshold 1 for dose rate alarm in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units

Number: threshold 2 for dose rate alarm in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units

Number: threshold 1 for dose alarm in μR , μrem or $0.01\mu\text{Sv}$ units

Number: threshold 2 for dose alarm in μR , μrem or $0.01\mu\text{Sv}$ units

AWP *Number Number..* Set scaler parameter:

Number: Preset count threshold 1 in counts (from 0 to 65000 counts)

Number: Preset count threshold 3 in counts (from 0 to 65000 counts)

Number: Preset time in seconds (from 0 to 9999 seconds)

AWN *Number Number..* Set background parameter:

Number: Preset count threshold 1 in counts (from 0 to 65000 counts)

Number: Preset count threshold 3 in counts (from 0 to 65000 counts)

Number: Preset time in seconds (from 0 to 9999 seconds)

AWB *Number Number* Set background value:

Number: Background value threshold 1 in 0.01 cps (from 0 to 65000 cps)

Number: Background value threshold 3 in 0.01 cps (from 0 to 65000 cps)

AWt *Number* Set Tau:

Number: Value in 0.1 s (from 10 to 600 means from 1.0 to 60.0s)

5.12.2 Measurement values

Z Read raw count rates with dead time correction

Response:

- Counter 1 in cps
- Counter 2 in cps
- Counter 3 in cps
- Counter 4 in cps
- HV power index
- Probe current in 0.1 μA units

z Read filtered count rate β and α threshold

Response:

- Value β threshold 0.01 cps
- Value α threshold in 0.01 cps

A Read display value and status
 Response:
 - display value β threshold in 0.01 cps, cpm, Bq, dps, dpm or Bq/cm² units
 - display value α threshold in 0.01 cps, cpm, Bq, dps dpm or Bq/cm² units
 or in μ R, μ rem or 0.01 μ Sv units in dose rate mode
 - Status (see 5.5.6.2)

5.12.3 Scaler remote control (from V1.73)

CG Start scaler with parameter previously set
 CS Stop scaler
 CI Scaler information
 Response:
 - 0: Scaler stopped, 1: Scaler active
 - Scaler counting time
 - accumulated counts threshold 1 or window 1
 - accumulated counts threshold 3 or window 2
 - accumulated counts threshold 1
 - accumulated counts threshold 2
 - accumulated counts threshold 3
 - accumulated counts threshold 4

5.12.4 Configuration flags

5.12.4.1 *Configuration flags 1 with kR / kW*

Bit number	
0	Disable key lock 0: no 1: yes
1	Not used
2	Show NBR Bargraph 0: off 1: on
3	Alarming Sound 0: off 1: on
4	Alarming LED 0: off 1: on
5	Alarming Vibration 0: off 1: on
6	NBR 0: off 1: on
7	Single Pulse 0: off 1: on

5.12.4.2 Configuration flags 2 with fR / fW

Bit number	
0	Show Gamma Symbol 0: no 1:yes
1	Not used
2	Not used
3	Alarm threshold read-only 0: off 1:on
4	Flag for overload (read-only)
5	Display of temperature 0: off 1:on
6	Temperature unit 0: °C 1: °F
7	Not used

5.12.4.3 Configuration flags 3 with KR / KW

Bit number	
0	Key tone 0:off 1:on
1	Autosend 0:off 1:on
2	Finder 0:off 1:on
3	Mean-Max 0:off 1:on
4	not used, write "0"
5	display of dose 0:off 1:on
6	not used
7	not used

5.12.4.4 Configuration flags 4 with jR / jW

Bit number	
0	Not used
1	Battery type 0:Alkaline 1:NiMh
2	Display rotation 0:No 1:yes
3	Scaler mode 0: Preset count 1: Preset time
4	Scaler, accumulated counts 0:No 1:yes
5	Scaler, after measurement 0: Stop 1: automatic restart
6	Scaler net 0:No 1:yes
7	Ratemeter net 0:No 1:yes
8	Alpha-LED 0:No 1:yes
9	Alpha-Sound 0:No 1:yes
10	Not used
11	Alarm, Background LCD 0:No 1:yes
12	Not used
13	Unit R/h and rem/h: lowest prefix 'm' instead of 'μ'
14	Not used
15	Not used

5.12.4.5 *Measuring unit with uR /uW*

uR Read measuring unit and operation mode
Response: Hex-value. See below

uRHex Write measuring unit and operation mode
See below

Bit number	
0	0x00 Display unit cps 0x05 Display unit Sv/h
1	0x01 Display unit cpm 0x06 Display unit R/h
2	0x02 Display unit Bq 0x07 Display unit rem/h
3	0x03 Display unit dps 0x08 Display unit Bq/cm ² 0x04 Display unit dpm 0x09 Display unit Gy/h (from V2.00)
4	0: Ratemeter, 1: Scaler
5	Display mode:
6	00: measurement value β threshold 01: measurement value α threshold 02: dual display 03: graphic display
7	Ratemeter 0: ADF 1: Tau

5.12.4.6 *Menu configuration*

mR Read configuration for main menu, submenu “Settings”, submenu “Alarm indication” and , submenu “Operation mode”
Response: Hex-values. See below

mRHex Hex Write configuration for main menu , submenu “Settings”, submenu “Alarm indication” and , submenu “Operation mode”
See below

5.12.4.6.1 Main menu

Bit number	
0	Switch off 0:hidden 1:visible
1	Background 0: hidden 1:visible
2	Select counter tube 0: hidden 1:visible
3	Backlight 0: hidden 1:visible
4	Measuring unit 0: hidden 1:visible
5	Operation mode 0: hidden 1:visible
6	Scaler parameter 0: hidden 1:visible
7	Nuclide table 0: hidden 1:visible
8	Alarm β 0: hidden 1:visible Alarm act. β Alarm dose rate

9	Alarm α Alarm act. α Alarm dose	0: hidden	1:visible
10	not used		
11	not used		
12	Clear dose	0: hidden	1:visible
13	Settings	0: hidden	1:visible
14	Alarm indication	0: hidden	1:visible
15	Show alarm	0: hidden	1:visible
16	Text Info	0: hidden	1:visible
17..31	not used, write "0"		

5.12.4.6.2 Sub menu "Settings"

Bit number			
0	Batt. type	0:hidden	1:visible
1	Autosend	0: hidden	1:visible
2	Single Pulse	0: hidden	1:visible
3	Finder	0: hidden	1:visible
4	Alpha-LED	0: hidden	1:visible
5	Alpha Sound	0: hidden	1:visible
6	Set Date/Time	0: hidden	1:visible
7	Not used, write "0"		
8	Alarm-NBR	0: hidden	1:visible
9	Set HV	0: hidden	1:visible
10	Bluetooth	0: hidden	1:visible
11	Set fixed TC (from V2.00)	0: hidden	1:visible
12	Language (from V2.00)	0: hidden	1:visible
13	Contrast (from V2.00)	0: hidden	1:visible
14...15	not used, write "0"		

5.12.4.6.3 Submenu "Alarm indication"


Bit number			
0	Sound	0:hidden	1:visible
1	LED	0: hidden	1:visible
2	Vibrator	0: hidden	1:visible
3	LCD LED	0: hidden	1:visible
4...8	not used, write "0"		

5.12.4.6.4 Submenu "Operation mode"

Bit number			
0	Ratemeter ADF	0:hidden	1:visible
1	Scaler	0: hidden	1:visible

5.12.6.2 History readout scaler

256 716612088 721 999 12 50 120 23



End of History:

End

5.12.6.3 History status

Decimal value converted in HEX:

Bit number		
0	Net value	0:No 1:Yes
1	Operation mode	0: Ratemeter 1:Scaler
2	Number of used probe	
3		
4		
5		
6	Background measurement	0:No 1:Yes
7	Display of	0: $\alpha\beta$ 1: β
8	Number of used nuclide	
9		
10		
11		
12	0x00 Display unit cps	0x05 Display unit Sv/h
13	0x01 Display unit cpm	0x06 Display unit R/h
14	0x02 Display unit Bq	0x07 Display unit rem/h

15	0x03 Display unit dps 0x04 Display unit dpm	0x08 Display unit Bq/cm ² 0x09 Display unit Gy/h (from V2.00)
16	Probe style	
17	0000: $\alpha\beta$ -Probe	0011: Dose rate probe
18	0001: α -Probe	0100: Gamma probe multi channel
19	0010: β -Probe	
	1011: GM-Probe (Dose rate)	1000: GM-Probe (Contamination)
20	Accumulated counts	0:No 1:Yes
21	Not used	
22	Not used	
23	Not used	
24	Window:	
25	000: No window. Display of Rate 1 Rate 3	
26	001: β window. Display Rate 1 – Rate 2 and Rate 3	
	010: α window. Display of Rate 1 and Rate 3 – Rate 4	
	011: β and α window. Display of Rate 1 – Rate 2 and Rate 3 – Rate 4	
	101: Display Rate 1 – Rate 2 and Rate 4	
27	Not used	
28	Not used	
29	Not used	
30	Reserved	
31	Reserved	

5.12.7 Nuclide table

nR	Reading nuclide data.
Number:	number of probe
Response:	Nuclide data (see below)

nW	NumberString	Write nuclide data
	Number:	number of probe
	String:	nuclide data (see below)

Nuclide data: every string contains up to 16 nuclide data of the corresponding probe

For example:

00 01 Sr-90 471 999 Am-241 330 123 0815

Field	Length	Start	End	Offset	Field
Checksum (in HEX)	4	0	3	0	Checksum (in HEX)
Factor for α threshold	4	4	7	4	Factor for α threshold
Factor for β threshold	4	8	11	8	Factor for β threshold
Name of nuclide 2 (do not use space)	16	12	27	12	Name of nuclide 2 (do not use space)
Factor for α threshold	4	28	31	20	Factor for α threshold
Factor for β threshold	4	32	35	24	Factor for β threshold
Name of nuclide 1 (do not use space)	16	36	51	28	Name of nuclide 1 (do not use space)
Active nuclide number	4	52	55	32	Active nuclide number
Number of nuclides in this string	4	56	59	36	Number of nuclides in this string

Activity calculation:

$$F=1/E_{eff}$$

$$\text{Ar}=\text{Cr}^*\text{F}$$

For example:

Count rate is 67.3 cps and efficiency for Sr-90 is 29%. Factor F is $1/0.29=3.448$.

The activity is $67.3 \times 3.45 = 232.18$ Bq.

Write “0” for null efficiency. In this case the RadEye displays “### Bq”.

Checksum calculation: Checksum starts from the beginning of the message to the end, with the checksum. Checksum is modulo 256.

5.12.8 Probe data

cR	Reading channel text for Dual Probes. Response: <i>C1 C3</i>
cW <i>C1 C3</i>	write channel text for Dual Probes. <i>C1</i> : Name for lower channel (e.g. R1-R2). Max. 5 Character <i>C3</i> : Name for upper channel (e.g. R4). Max. 5 Character
zR <i>Number</i>	Reading probe data. <i>Number</i> : consecutive number Response: probe data (see below)
zW <i>NumberString</i>	Write probe data <i>Number</i> : consecutive number <i>String</i> : probe data
zRA	Read the number of stored probes. Response: value from 0...15
zW <i>ANumber</i>	Read the number of stored probes. <i>Number</i> : value from 0...15
zRG	Reading active probe. Response: value from 0...15
zW <i>GNumber</i>	Write active probe. <i>Number</i> : value from 0...15

For example:

HP-380AB 10000 5 5 5 5 600 100 100 30 100 10 1 3 5 8 30 225 1000 1250 1234

1	2	3	4	5	6	7	8	9	A	B	C	D	E	F	G	H	I	J	K	L
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

1	Probe name. Max. 11 characters. Do not use space
2	Calibration factor in (nSv/h)/cps
3	Dead time Rate 1 μ s
4	Dead time Rate 2 μ s
5	Dead time Rate 3 μ s
6	Dead time Rate 4 μ s
7	High voltage in Volt
8	Overload Rate 1 in kcps
9	Overload Rate 3 in kcps
A	Overload probe current μ A

B	Area of this probe
C	Timeout for detector error in seconds
D	Factor Rate 1 for dose rate calculation. Range: -128...+127
E	Factor Rate 2 for dose rate calculation. Range: -128...+127
F	Factor Rate 3 for dose rate calculation. Range: -128...+127
G	Factor Rate 4 for dose rate calculation. Range: -128...+127
H	Threshold Rate 1 in mV. Range: 30..1220mV
I	Threshold Rate 2 in mV. Range: 30..1220mV
J	Threshold Rate 3 in mV. Range: 30..1220mV
K	Threshold Rate 4 in mV. Range: 30..1220mV
L	Flags

Flags

Bit number		
0	Window:	
1	00: No window. Display of Rate 1 Rate 3 01: β window. Display Rate 1 – Rate 2 and Rate 3 10: α window. Display of Rate 1 and Rate 3 – Rate 4 11: β and α window. Display of Rate 1 – Rate 2 and Rate 3 – Rate 4	
2	Pulse fade out 0: off 1: on	
3	Probe style	
4	0000: $\alpha\beta$ -Probe	0011: Dose rate probe
5	0001: α -Probe	0100: Gamma probe multi channel
6	0010: β -Probe	
	1011: GM-Probe (Dose rate)	1000: GM-Probe (Contamination)
7	FHZ674-NBR	
8	Single Pulse 0: every Count, 1: with Divide Ratio (see Command ARC)	
9	Not used	
10	Display Rate 4 in place of Rate 3 (from V2.00)	
11	Use long integration time 0:No 1:Yes	
12	Calibration factor 0: Hx 1: H*10	
13	Not used	
14	Probe is calibrated 0:No 1:Yes	
15	Must be "0"	

Dead time correction:

$$C1 = \frac{C1_{raw}}{1 - \tau1 * C1_{raw}}$$

$\tau1$ Dead time
C1raw: raw count rate

1- $\tau_1 * C_{1raw}$ is limited to min. 0.1

Dose rate calculation:

$$DR = Cf * (D * Rate1 + E * Rate2 + F * Rate3 + G * Rate4)$$

DR: Dose rate in nSv/h

Cf: Calibration factor in nSv/h/cps

Rate 1..4 Count rate of threshold 1..4 in cps

5.12.9 Event log

6656 520549251

	\-----	Date and time as a decimal value (see 3.2)
\-----		Event log as a decimal value

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Probe Type Part 1 *)
7	
8	Operation mode 0: Ratemeter 1:Scaler
9	Probe Type Part 2 *)
10	
11	Sound 0: off 1:on
12	LED 0: off 1:on
13	Vibration alarm 0: off 1:on
14	Clear dose 0: no 1:yes
15	1: Mode activity
16	1: Dose rate, count rate or activity -alarm
17	1: Alarm dose
18	NBR-Alarm
19	
	01: Low energy
	10: High energy
	11: Anomaly
20	Dose rate, count rate or activity > alarm threshold 1
21	Dose rate, count rate or activity > alarm threshold 2
22	Dose > alarm threshold 1 (Dose rate probe) count rate or activity R3 > alarm threshold 1 (Contamination, Dual Channel)
23	Dose > alarm threshold 1 (Dose rate probe) count rate or activity R3 > alarm threshold 1 (Contamination, Dual Channel)
24	1: Window Alpha
25	1: Safety alarm

26	Power off
27	Power on
28	1: NBR active
29	1: Mid Energy
30	1: Window Beta
31	Not used

*) Probe Type

xxxx x00x 00xx xxxx

|||| |||| |||| ||||

|||| |++---+----- Probe Type

00x00: $\alpha\beta$ -Probe	00x11: Dose rate probe
00x01: α -Probe	01x00: Gamma probe multi channel
00x10: β -Probe	
10x11: GM-Probe (Dose rate)	10x00: GM-Probe (Contamination)

5.12.10 Automatic sending

<STX>721 5 0 0 14 SX 1234 AB<ETX><CR><LF>

								\---- BCC in HEX
								\----- Dose in μ R/h-units
								\----- information RadEye SX
								\----- Status as Hex-value
								\----- Measuring unit Rate 3
								\----- Rate 3
								\----- Measuring unit Rate 1
								\----- Rate 1

Formation of a BCC (block check character):

Modulo 256 sum of <STX> up to the last character before the BCC (including), coded as hexadecimal ASCII-number (e.g. 1F).

Status

Bit number	
0	Not used
1	1: Overload
2	1: Alarm
3	Not used
4	Not used
5	1: Battery voltage low

6	not used
7	not used

Measuring unit for Rate 1 and Rate 3

5	cps
3	cpm
6	Bq
9	dps
7	dpm
8	Bq/cm ²
2	R/h
0	Sv/h
10	Rem/h
1	Gy/h

5.12.11 Status information

F Reading status information
 Response Number with status information

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	Not used
9	
10	
11	not used,
12	Alarm threshold (read-only) 0: off 1:on
13	Flag for overload
14	Temperature display 0: off 1:on
15	Not used
16	1: Alarm
17	Not used
18	Not used
19	Not used
20	Value > alarm threshold 1 (depending on display mode)
21	Value > alarm threshold 2 (depending on display mode)

22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24	Not used
25	Not used
26	Not used
27	Not used
28	Not used
29	Not used
30	Not used
31	Not used

5.12.12 Miscellaneous

ARC Read single pulse divider
 Response: number for single pulse divider

ARt Read fixed time constant (from V2.00)
 Response: time in 0.1s units

AWC*Number* Write single pulse divider
 Number: value from 1 to 65535

AWt Write fixed time constant (from V2.00)
 Number: value from 10 to 65535

5.13 RadEye PX

Used firmware version: 1.63

5.13.1 Limit values

AR1	Reading the threshold 1 for activity alarm β threshold. Response: number in 0.01 Bq units.
AR2	Reading the threshold 2 for activity alarm β threshold. Response: number in 0.01 Bq units.
AR3	Reading the threshold 1 for activity alarm α threshold. Response: number in 0.01 Bq units.
AR4	Reading the threshold 2 for activity alarm α threshold. Response: number in 0.01 Bq units.
AR5	Reading the threshold 1 for count rate alarm β threshold. Response: number in 0.01 cps units
AR6	Reading the threshold 2 for count rate alarm β threshold. Response: number in 0.01 cps units
AR7	Reading the threshold 1 for count rate alarm α threshold. Response: number in 0.01 cps units
AR8	Reading the threshold 2 for count rate alarm α threshold. Response: number in 0.01 cps units
AR9	Reading alarm thresholds for dose rate and dose Response: <ul style="list-style-type: none">- threshold 1 for dose rate alarm in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units- threshold 2 for dose rate alarm in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units- threshold 1 for dose alarm in μR, μrem or $0.01\mu\text{Sv}$ units- threshold 2 for dose alarm in μR, μrem or $0.01\mu\text{Sv}$ units
ARP	Reading scaler parameter: <ul style="list-style-type: none">- Preset count threshold 1. Response: Number in counts- Preset count threshold 3. Response: Number in counts- Preset time. Response: Number in seconds
ARN	Reading background parameter: <ul style="list-style-type: none">- Preset count threshold 1. Response: Number in counts

- Preset count threshold 3. Response: Number in counts
- Preset time. Response: Number in seconds

ARB	Reading background values: <ul style="list-style-type: none"> - Background value threshold 1. Response: Value in 0.01 cps - Background value threshold 3. Response: Value in 0.01 cps
AW1Number	Setting the threshold 1 for activity alarm β threshold. <i>Number:</i> value in 0.01 Bq units.
AW2Number	Setting the threshold 2 for activity alarm β threshold <i>Number:</i> value in 0.01 Bq units
AW3Number	Setting the threshold 1 for activity alarm α threshold <i>Number:</i> value in 0.01 Bq units.
AW4Number	Setting the threshold 2 for activity alarm α threshold <i>Number:</i> value in 0.01 Bq units
AW5Number	Setting the threshold 1 for count rate alarm β threshold.. <i>Number:</i> value in 0.01 cps units.
AW6Number	Setting the threshold 2 for count rate alarm β threshold. <i>Number:</i> value in 0.01 cps units.
AW7Number	Setting the threshold 1 for count rate alarm α threshold. <i>Number:</i> value in 0.01 cps units
AW8Number	Setting the threshold 2 for count rate alarm α threshold. <i>Number:</i> value in 0.01 cps units
AW9Number	Setting alarm thresholds for dose rate and dose <i>Number:</i> threshold 1 for dose rate alarm in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units <i>Number:</i> threshold 2 for dose rate alarm in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units <i>Number:</i> threshold 1 for dose alarm in μR , μrem or $0.01\mu\text{Sv}$ units <i>Number:</i> threshold 2 for dose alarm in μR , μrem or $0.01\mu\text{Sv}$ units
AWPNumber	<i>Number..</i> Set scaler parameter: <i>Number:</i> Preset count threshold 1 in counts (from 0 to 9999 counts) <i>Number:</i> Preset count threshold 3 in counts (from 0 to 9999 counts) <i>Number:</i> Preset time in seconds (from 0 to 9999 seconds)
AWNNumber	<i>Number..</i> Set background parameter: <i>Number:</i> Preset count threshold 1 in counts (from 0 to 9999 counts) <i>Number:</i> Preset count threshold 3 in counts (from 0 to 9999 counts)

Number: Preset time in seconds (from 0 to 9999 seconds)

AWBNumber Number Set background value:

Number: Background value threshold 1 in 0.01 cps (from 0 to 100 cps)

Number: Background value threshold 3 in 0.01 cps (from 0 to 100 cps)

5.13.2 Measurement values

Z Read raw count rates with dead time correction

Response:

- Counter 1 in cps
- Counter 2 in cps
- Counter 3 in cps
- Reserved
- HV power index

z Read filtered count rate β and α threshold

Response:

- Value β threshold 0.01 cps
- Value α threshold in 0.01 cps

A

Read display value and status

Response:

- display value β threshold in 0.01 cps, cpm, Bq, dps, dpm or Bq/cm² units
- display value α threshold in 0.01 cps, cpm, Bq, dps dpm or Bq/cm² units or in μ R, μ rem or 0.01 μ Sv units in dose rate mode
- Status (see 5.5.6.2)

5.13.3 Configuration flags

5.13.3.1 Configuration flags 1 with *kR / kW*

Bit number	
0	Not used
1	Not used
2	Not used
3	Alarming Sound 0: off 1: on
4	Alarming LED 0: off 1: on
5	Alarming Vibration 0: off 1: on
6	Not used
7	Single Pulse 0: off 1: on

5.13.3.2 Configuration flags 2 with *fR / fW*

Bit number	
0	Not used
1	Not used
2	Not used
3	Alarm threshold read-only 0: off 1:on
4	Flag for overload (read-only)
5	Display of temperature 0: off 1:on
6	Temperature unit 0: °C 1: °F
7	Not used

5.13.3.3 *Configuration flags 3 with KR / KW*

Bit number	
0	Key tone 0:off 1:on
1	Autosend 0:off 1:on
2	Finder 0:off 1:on
3	Mean-Max 0:off 1:on
4	not used, write "0"
5	display of dose 0:off 1:on
6	not used
7	not used

5.13.3.4 *Configuration flags 4 with jR / jW*

Bit number	
0	Not used
1	Battery type 0:Alkaline 1:NiMh
2	Display rotation 0:No 1:yes
3	Scaler mode 0: Preset count 1: Preset time
4	Not used
5	Scaler, after measurement 0: Stop 1: automatic restart
6	Scaler net 0:No 1:yes
7	Ratemeter net 0:No 1:yes
8	Alpha-LED 0:No 1:yes
9	Alpha-Sound 0:No 1:yes
10	Not used
11	Not used
12	Not used
13	Ratemeter Tau active 0:No 1:yes
14	Not used
15	Not used

5.13.3.5 *Measuring unit with uR /uW*

uR Read measuring unit and operation mode
 Response: Hex-value. See below
 uRHex Write measuring unit and operation mode
 See below

Bit number	
0	0x00 Display unit cps 0x05 Display unit Sv/h
1	0x01 Display unit cpm 0x06 Display unit R/h
2	0x02 Display unit Bq 0x07 Display unit rem/h
3	0x03 Display unit dps 0x08 Display unit Bq/cm ² 0x04 Display unit dpm 0x09 Display unit Gy/h
4	0: Ratemeter (Tau or ADF), 1: Scaler
5	Display mode:
6	00: measurement value β threshold 01: measurement value α threshold 02: dual display 03: graphic display
7	0: two decimal places 1: one decimal place

5.13.3.6 *Menu configuration*

mR Read configuration for main menu and submenu “Settings”
Response: Hex-values. See below

mRHex Hex Write configuration for main menu and submenu “Settings”
See below

5.13.3.6.1 Main menu

Bit number			
0	Switch off	0:hidden	1:visible
1	Background	0: hidden	1:visible
2	Select counter tube	0: hidden	1:visible
3	Backlight	0: hidden	1:visible
4	Measuring unit	0: hidden	1:visible
5	Operation mode	0: hidden	1:visible
6	Scaler parameter	0: hidden	1:visible
7	Nuclide table	0: hidden	1:visible
8	Alarm β Alarm act. β Alarm dose rate	0: hidden	1:visible
9	Alarm α Alarm act. α Alarm dose	0: hidden	1:visible
10	not used		
11	not used		
12	Clear dose	0: hidden	1:visible
13	Settings	0: hidden	1:visible
14	Alarm indication	0: hidden	1:visible
15	Show alarm	0: hidden	1:visible
16	Text Info	0: hidden	1:visible
17	not used, write “0”		
18	not used, write “0”		
19	not used, write “0”		
20	not used, write “0”		
21	not used, write “0”		
22	not used, write “0”		
23	not used, write “0”		
24	not used, write “0”		
25	not used, write “0”		
26	not used, write “0”		
27	not used, write “0”		
28	not used, write “0”		
29	not used, write “0”		
30	not used, write “0”		
31	not used, write “0”		

5.13.3.6.2 Submenu “Settings”

Bit number			
0	Batt. type	0:hidden	1:visible
1	Autosend	0: hidden	1:visible
2	Single Pulse	0: hidden	1:visible
3	Finder	0: hidden	1:visible
4	Alpha-LED	0: hidden	1:visible
5	Alpha Sound	0: hidden	1:visible
6	Set Date/Time	0: hidden	1:visible
7	Display $\alpha\beta$ / β	0: hidden	1:visible
8	Not used		
9	Set HV	0: hidden	1:visible
10	Bluetooth	0: hidden	1:visible
11	Edit Tau	0: hidden	1:visible
12	Language	0: hidden	1:visible
13	Contrast	0: hidden	1:visible
14	not used		
15	not used		

5.13.3.6.3 Submenu “Alarm indication”

Bit number			
0	Sound	0:hidden	1:visible
1	LED	0: hidden	1:visible
2	Vibrator	0: hidden	1:visible
3...8	not used, write “0”		

5.13.3.6.4 Submenu “Operation mode”

Bit number			
0	Ratometer ADF	0:hidden	1:visible
1	Scaler	0: hidden	1:visible
2	Ratometer Tau	0: hidden	1:visible

5.13.4 High voltage

HR Reading high voltage
 Response: value in Volt.

5.13.5 History output

5.13.5.1 History readout ratemeter

```
256 716612088 721 999 12 50 120 23
-----
|      |      |      |      |      |      |      |
|      |      |      |      |      |      | \---- Temperature in °C
|      |      |      |      |      |      | \----- Measuring time in seconds
|      |      |      |      |      | \----- Max value α channel
|      |      |      |      |      |      Count rate in 0.01 cps
|      |      |      |      |      |      Activity in 0.01 Bq
|      |      |      |      |      |      Dose rate in 1 μR/h
|      |      |      |      | \----- Mean value α channel
|      |      |      |      |      count rate in 0.01 cps,
|      |      |      |      |      Activity in 0.01 Bq
|      |      |      |      |      Dose rate in 1 μR/h
|      |      |      | \----- Max value β channel
|      |      |      |      count rate in 0.01 cps,
|      |      |      |      Activity in 0.01 Bq
|      |      | \----- Mean value β channel
|      |      |      count rate in 0.01 cps,
|      |      |      Activity in 0.01 Bq
|      | \----- Date and Time as a decimal value
| \----- History status as a decimal value
\----- see below
```

5.13.5.2 History readout scaler

```
256 716612088 721 999 12 50 120 23
-----
|      |      |      |      |      |      |      |
|      |      |      |      |      |      | \---- Temperature in °C
|      |      |      |      |      |      | \----- Measuring time in seconds
|      |      |      |      |      | \----- Background value α channel
|      |      |      |      |      |      Count rate in 0.01 cps
|      |      |      |      | \----- Measuring value α channel
|      |      |      |      |      count rate in 0.01 cps,
|      |      |      |      |      Activity in 0.01 Bq
|      |      |      |      |      Dose rate in 1 μR/h
|      |      |      | \----- Background value β channel
|      |      |      |      count rate in 0.01 cps,
|      |      | \----- Measuring value β channel
|      |      |      count rate in 0.01 cps,
|      |      |      Activity in 0.01 Bq
|      | \----- Date and Time as a decimal value
| \----- History status as a decimal value
\----- see below
```

End of History:

End

5.13.5.3 *History status*

Decimal value converted in HEX:

Bit number			
0	Net value	0:No	1:Yes
1	Operation mode	0: Ratemeter	1:Scaler
2	Number of used probe		
3			
4			
5			
6	Background measurement	0:No	1:Yes
7	Display of	0: $\alpha\beta$	1: β
8	Number of used nuclide		
9			
10			
11			
12	0x00 Display unit cps	0x05 Display unit Sv/h	
13	0x01 Display unit cpm	0x06 Display unit R/h	
14	0x02 Display unit Bq	0x07 Display unit rem/h	
15	0x03 Display unit dps	0x08 Display unit Bq/cm ²	
	0x04 Display unit dpm	0x09 Display unit Gy/h	
16	Probe style		
17	000: $\alpha\beta$ -Probe		011: Dose rate probe
18	001: α -Probe		100: Gamma probe multi channel
	010: β -Probe		
19	Not used		
20	Accumulated counts	0:No	1:Yes
21	Not used		
22	Not used		
23	Not used		
24	Not used		
25	Not used		
26	Not used		
27	Not used		
28	Not used		
29	Not used		
30	Not used		
31	Not used		

5.13.6 Nuclide table

nRNumber Reading nuclide data.
Number: number of probe
 Response: Nuclide data (see below)

nWNumberString Write nuclide data
Number: number of probe
String: nuclide data (see below)

Nuclide data: every string contains up to 16 nuclide data of the corresponding probe

For example:

00 01 Sr-90 471 999 Am-241 330 123 usw. 0815

\--	Checksum (in HEX)
\-----	Factor for α threshold
\-----	Factor for β threshold
\-----	Name of nuclide 2 (do not use space)
\-----	Factor for α threshold
\-----	Factor for β threshold
\-----	Name of nuclide 1 (do not use space)
\-----	Active nuclide number
\-----	Number of nuclides in this string

Activity calculation:

$$F=1/E_{eff}$$

$$Ar=Cr \cdot F$$

For example:

Count rate is 67.3 cps and efficiency for Sr-90 is 29%. Factor F is $1/0.29=3.448$.

The activity is $67.3 \cdot 3.45=232.18$ Bq.

Write "0" for null efficiency. In this case the RadEye displays "### Bq".

Checksum calculation: Checksum starts from the beginning of the message to the end, with the checksum. Checksum is modulo 256.

5.13.7 Probe data

<i>zRNumber</i>	Reading probe data. <i>Number</i> : consecutive number Response: probe data (see below)
<i>zWNumberString</i>	Write probe data <i>Number</i> : consecutive number <i>String</i> : probe data
<i>zRA</i>	Read the number of stored probes. Response: value from 0...15
<i>zWANumber</i>	Read the number of stored probes. <i>Number</i> : value from 0...15
<i>zRG</i>	Reading active probe. Response: value from 0...15
<i>zWGNumber</i>	Write active probe. <i>Number</i> : value from 0...15

For example:

HP-380AB 10000 5 5 5 5 600 100 100 30 100 10 1 3 5 8 30 225 1000 1250 1234

1 2 3 4 5 6 7 8 9 A B C D E F G H I J K L

1	Probe name. Max. 11 characters. Do not use space
2	Calibration factor in 10(nSv/h)/cps units
3	Dead time Rate 1 μ s
4	Dead time Rate 2 μ s
5	Dead time Rate 3 μ s
6	Not used
7	High voltage in Volt
8	Overload Rate 1 in kcps
9	Overload Rate 3 in kcps
A	Not used
B	Area of this probe
C	Timeout for detector error in seconds
D	Factor Rate 1 for dose rate calculation. Range: -128...+127
E	Factor Rate 2 for dose rate calculation. Range: -128...+127
F	Factor Rate 3 for dose rate calculation. Range: -128...+127

G
H
I
J
K
L

Flags

Bit number		
0	Window: 00: No window. Display of Rate 1 Rate 3 01: β window. Display Rate 1 – Rate 2 and Rate 3	
1		
2	Pulse fade out 0: off 1: on	
3	Probe style	
4	000: $\alpha\beta$ -Probe	011: Dose rate probe
5	001: α -Probe	100: Gamma probe multi channel
	010: β -Probe	
6	Reserved, write “0”	
7	Reserved, write “0”	
8	Single Pulse 0: every Count, 1: with Divide Ratio (see Command ARC)	
9..15	Not used	

Dead time correction:

$$C1 = \frac{C1_{raw}}{1 - \tau1 * C1_{raw}}$$

τ_1	Dead time
----------	-----------

C1raw: raw count rate

1- $\tau_1 \cdot C_{1\text{raw}}$ is limited to min. 0.1

Dose rate calculation:

$$DR = Cf * (\overline{D} * Rate1 + \overline{E} * Rate2 + \overline{F} * Rate3)$$

DR: Dose rate in nSv/h

Cf: Calibration factor in nSv/h/cps

Rate 1..3	Count rate of threshold 1..3 in cps
-----------	-------------------------------------

5.13.8 Event log

6656 520549251

```

|      |
|      \-----
\-----

```

Date and time as a decimal value (see 3.2)

Event log as a decimal value

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	Operation mode 0: Ratemeter 1:Scaler
9	Not used
10	
11	Sound 0: off 1:on
12	LED 0: off 1:on
13	Vibration alarm 0: off 1:on
14	Clear dose 0: no 1:yes
15	1: Alarm threshold changed
16	1: Dose rate, count rate or activity -alarm
17	1: Alarm dose
18	NBR-Alarm 01: Low energy 10: High energy 11: Anomaly
19	
20	Dose rate, count rate or activity > alarm threshold 1
21	Dose rate, count rate or activity > alarm threshold 2
22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24	1: Scaler or Background parameter changed
25	Not used
26	Power off
27	Power on
28	Not used
29	Not used
30	Not used
31	Not used

5.13.9 Automatic sending

```

<STX>721  5 0 0 14 PX 1234 AB<ETX><CR><LF>
|         | | | | | | | |
|         | | | | | | | \---- BCC in HEX
|         | | | | | | | \----- Dose in µR/h
|         | | | | | | \----- information RadEye GX
|         | | | | | \----- Status as Hex-value
|         | | | | \----- Fixed to "0"
|         | | | \----- Fixed to "0"
|         | | \----- Fixed to "5"
|         | \----- Count rate 1 in 0.01cps
| \-----

```

Formation of a BCC (block check character):

Modulo 256 sum of <STX> up to the last character before the BCC (including), coded as hexadecimal ASCII-number (e.g. 1F).

Status

Bit number	
0	Not used
1	1: Overload
2	1: Alarm
3	Not used
4	Not used
5	1: Battery voltage low
6	not used
7	not used

5.13.10 Status information

F Reading status information
 Response Number with status information

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	Not used
9	
10	
11	not used,
12	Alarm threshold (read-only) 0: off 1:on
13	Flag for overload
14	Temperature display 0: off 1:on
15	Not used
16	1: Alarm
17	Not used
18	Not used
19	Not used
20	Value > alarm threshold 1 (depending on display mode)
21	Value > alarm threshold 2 (depending on display mode)
22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24	Not used
25	Not used
26	Not used
27	Not used
28	Not used
29	Not used
30	Not used
31	Not used

5.13.11 Miscellaneous

ARC Read single pulse divider
 Response: number for single pulse divider

ARt	Read fixed time constant (from V1.69) Response: time in 0.1s units
AWC <i>Number</i>	Write single pulse divider Number: value from 1 to 65535
AWt	Write fixed time constant (from V1.69) Number: value from 10 to 65535

5.14 RadEye HEC

Used firmware version: 1.69

5.14.1 Limit values

AR1	Reading the threshold 1 for activity alarm β threshold. Response: number in 0.01 Bq units.
AR2	Reading the threshold 2 for activity alarm β threshold. Response: number in 0.01 Bq units.
AR3	Reading the threshold 1 for activity alarm α threshold. Response: number in 0.01 Bq units.
AR4	Reading the threshold 2 for activity alarm α threshold. Response: number in 0.01 Bq units.
AR5	Reading the threshold 1 for count rate alarm β threshold. Response: number in 0.01 cps units
AR6	Reading the threshold 2 for count rate alarm β threshold. Response: number in 0.01 cps units
AR7	Reading the threshold 1 for count rate alarm α threshold. Response: number in 0.01 cps units
AR8	Reading the threshold 2 for count rate alarm α threshold. Response: number in 0.01 cps units
AW1Number	Setting the threshold 1 for activity alarm β threshold. Number: value in 0.01 Bq units.
AW2Number	Setting the threshold 2 for activity alarm β threshold Number: value in 0.01 Bq units
AW3Number	Setting the threshold 1 for activity alarm α threshold Number: value in 0.01 Bq units.
AW4Number	Setting the threshold 2 for activity alarm α threshold Number: value in 0.01 Bq units
AW5Number	Setting the threshold 1 for count rate alarm β threshold.. Number: value in 0.01 cps units.

<i>AW6Number</i>	Setting the threshold 2 for count rate alarm β threshold. <i>Number</i> : value in 0.01 cps units.
<i>AW7Number</i>	Setting the threshold 1 for count rate alarm α threshold. <i>Number</i> : value in 0.01 cps units
<i>AW8Number</i>	Setting the threshold 2 for count rate alarm α threshold. <i>Number</i> : value in 0.01 cps units

5.14.2 Measurement values

Z	Read raw count rates with dead time correction Response: <ul style="list-style-type: none"> - Counter 1 in cps - Counter 2 in cps - Counter 3 in cps - Counter 4 in cps - HV power index - Probe current in 0.1 μA units
z	Read filtered count rate β and α threshold Response: <ul style="list-style-type: none"> - Value β threshold 0.01 cps - Value α threshold in 0.01 cps

5.14.3 Scaler remote control (from V1.69)

CG	Start scaler with parameter previously set
CS	Stop scaler
CI	Scaler information Response: <ul style="list-style-type: none"> - 0: Scaler stopped, 1: Scaler active - Scaler counting time - accumulated counts threshold 1 or window 1 - accumulated counts threshold 3 or window 2 - accumulated counts threshold 1 - accumulated counts threshold 2 - accumulated counts threshold 3 - accumulated counts threshold 4

A

Read display value and status

Response:

- display value β threshold in 0.01 cps, cpm, Bq, dps, dpm or Bq/cm² units
- display value α threshold in 0.01 cps, cpm, Bq, dps dpm or Bq/cm² units
- Status (see 5.5.6.2)

5.14.4 Configuration flags

5.14.4.1 Configuration flags 1 with *kR* / *kW*

Bit number	
0	Not used
1	Not used
2	Not used
3	Alarming Sound 0: off 1: on
4	Alarming LED 0: off 1: on
5	Alarming Vibration 0: off 1: on
6	Not used
7	Single Pulse 0: off 1: on

5.14.4.2 Configuration flags 2 with *fR* / *fW*

Bit number	
0	Not used
1	Not used
2	Not used
3	Alarm threshold read-only 0: off 1:on
4	Flag for overload (read-only)
5	Display of temperature 0: off 1:on
6	Temperature unit 0: °C 1: °F
7	Not used

5.14.4.3 Configuration flags 3 with KR / KW

Bit number			
0	Key tone	0:off	1:on
1	Autosend	0:off	1:on
2	Finder	0:off	1:on
3	Mean-Max	0:off	1:on
4	not used, write "0"		
5	not used, write "0"		
6	Preset Count is read only	0:off	1:on
7	Beep at the end of a scaler measurement	0:off	1:on

5.14.4.4 Configuration flags 4 with jR / jW

Bit number			
0	Not used		
1	not used, write “1”		
2	not used, write “0”		
3	Scaler mode	0: Preset count 1: Preset time	
4	Not used		
5	Scaler, after measurement	0: Stop	1: automatic restart
6	Scaler net	0:No	1:yes
7	Ratemeter net	0:No	1:yes
8	Alpha-LED	0:No	1:yes
9	Alpha-Sound	0:No	1:yes
10	Not used		
11	Not used		
12	Not used		
13	Not used		
14	Not used		
15	Not used		

5.14.4.5 *Measuring unit with uR /uW*

uR Read measuring unit and operation mode
 Response: Hex-value. See below
 uRHex Write measuring unit and operation mode
 See below

Bit number	
0	0x00 Display unit cps 0x04 Display unit dpm
1	0x01 Display unit cpm 0x08 Display unit Bq/cm ²
2	0x02 Display unit Bq
3	0x03 Display unit dps
4	0: Ratemeter, 1: Scaler
5	Display mode:
6	00: measurement value β threshold 01: measurement value α threshold 02: dual display
7	0: two decimal places 1: one decimal place

5.14.4.6 Menu configuration

mR Read configuration for main menu and submenu “Settings”
 Response: Hex-values. See below

mRHex Hex Write configuration for main menu and submenu “Settings”
 See below

Bit number			
0	Switch off	0:hidden	1:visible
1	Background	0: hidden	1:visible
2	Select counter tube	0: hidden	1:visible
3	Backlight	0: hidden	1:visible
4	Measuring unit	0: hidden	1:visible
5	Operation mode	0: hidden	1:visible
6	Scaler parameter	0: hidden	1:visible
7	Nuclide table	0: hidden	1:visible
8	Alarm β Alarm act. β	0: hidden	1:visible
9	Alarm α Alarm act. α	0: hidden	1:visible
10	not used		
11	not used		
12	not used		
13	Settings	0: hidden	1:visible
14	Alarm indication	0: hidden	1:visible
15	Show alarm	0: hidden	1:visible
16	Text Info	0: hidden	1:visible
17	not used, write “0”		
18	not used, write “0”		
19	not used, write “0”		
20	not used, write “0”		
21	not used, write “0”		
22	not used, write “0”		
23	not used, write “0”		
24	not used, write “0”		
25	not used, write “0”		
26	not used, write “0”		
27	not used, write “0”		
28	not used, write “0”		
29	not used, write “0”		
30	not used, write “0”		
31	not used, write “0”		

Submenu “Settings”

Bit number			
0	Batt. type	0:hidden	1:visible
1	Autosend	0: hidden	1:visible
2	Single Pulse	0: hidden	1:visible
3	Finder	0: hidden	1:visible
4	Alpha-LED	0: hidden	1:visible
5	Alpha Sound	0: hidden	1:visible
6	Set Date/Time	0: hidden	1:visible
7	Display α/β	0: hidden	1:visible
8	Source test		
9	Set HV	0: hidden	1:visible
10	not used, write “0”		
11	Language (<i>from V1.69</i>)	0: hidden	1:visible
12	Contrast (<i>from V1.69</i>)	0: hidden	1:visible
13..15	not used, write “0”		

5.14.5 High voltage

HR Reading high voltage
 Response: value in Volt.

5.14.6 Remote control (from V1.68)

CG Start a scaler measurement
 CS Stop a scaler measurement
 CI Scaler measurement information:

- Scaler measurement 0: stopped 1: active
- expired time
- accumulated counts alpha/beta or beta channel
- accumulated counts alpha channel (Rate 3)
- accumulated counts counter 1
- accumulated counts counter 2
- accumulated counts counter 3
- accumulated counts counter 4

5.14.7 History output

5.14.7.1 History readout ratemeter

```
256 716612088 721 999 12 50 120 23
-----
|      |      |      |      |      |      |      |
|      |      |      |      |      |      | \---- Temperature in °C
|      |      |      |      |      |      | \----- Measuring time in seconds
|      |      |      |      |      | \----- Max value α channel
|      |      |      |      |      |      Count rate in 0.01 cps
|      |      |      |      |      |      Activity in 0.01 Bq
|      |      |      |      | \----- Mean value α channel
|      |      |      |      |      count rate in 0.01 cps,
|      |      |      |      |      Activity in 0.01 Bq
|      |      |      | \----- Max value β channel
|      |      |      |      count rate in 0.01 cps,
|      |      |      |      Activity in 0.01 Bq
|      |      | \----- Mean value β channel
|      |      |      count rate in 0.01 cps,
|      |      |      Activity in 0.01 Bq
|      | \----- Date and Time as a decimal value
| \----- History status as a decimal value
| see below
```

5.14.7.2 History readout scaler

```
256 716612088 721 999 12 50 120 23
-----
|      |      |      |      |      |      |      |
|      |      |      |      |      |      | \---- Temperature in °C
|      |      |      |      |      |      | \----- Measuring time in seconds
|      |      |      |      |      | \----- Background value α channel
|      |      |      |      |      |      Count rate in 0.01 cps
|      |      |      |      | \----- Measuring value α channel
|      |      |      |      |      count rate in 0.01 cps,
|      |      |      |      |      Activity in 0.01 Bq
|      |      |      | \----- Background value β channel
|      |      |      |      count rate in 0.01 cps,
|      |      | \----- Measuring value β channel
|      |      |      count rate in 0.01 cps,
|      |      |      Activity in 0.01 Bq
|      | \----- Date and Time as a decimal value
| \----- History status as a decimal value
| see below
```

End of History:

End

5.14.7.3 *History status*

Decimal value converted in HEX:

Bit number			
0	Net value	0:No	1:Yes
1	Accumulated counts	0:No	1:Yes
2	Not used		
3	Number of used parameter set		
4			
5			
6	Background measurement	0:No	1:Yes
7	Beta Window	0: $\alpha\beta$	1: β
8	Number of used nuclide		
9			
10			
11			
12	0x00 Display unit cps	0x04 Display unit dpm	
13	0x01 Display unit cpm	0x08 Display unit Bq/cm ²	
14	0x02 Display unit Bq		
15	0x03 Display unit dps		

5.14.8 Nuclide table

nRNumber Reading nuclide data.
Number: number of probe
 Response: Nuclide data (see below)

nWNumberString Write nuclide data
Number: number of probe
String: nuclide data (see below)

Nuclide data: every string contains up to 16 nuclide data of the corresponding parameter set

For example:

00 01 Sr-90 471 999 Am-241 330 123 usw. 0815

\--	Checksum (in HEX)
\-----	Factor for α threshold
\-----	Factor for β threshold
\-----	Name of nuclide 2 (do not use space)
\-----	Factor for α threshold
\-----	Factor for β threshold
\-----	Name of nuclide 1 (do not use space)
\-----	Active nuclide number
\-----	Number of nuclides in this string

Activity calculation:

$$F=1/E_{eff}$$

$$Ar=Cr \cdot F$$

For example:

Count rate is 67.3 cps and efficiency for Sr-90 is 29%. Factor F is $1/0.29=3.448$.

The activity is $67.3 \cdot 3.45=232.18$ Bq.

Write "0" for null efficiency. In this case the RadEye displays "### Bq".

Checksum calculation: Checksum starts from the beginning of the message to the end, with the checksum. Checksum is modulo 256.

5.14.9 Parameter data

zRNumber Reading parameter data.
Number: consecutive number
 Response: parameter data (see below)

zWNumberString Write parameter data
Number: consecutive number
String: parameter data

zRA Read the number of stored parameters.
 Response: value from 0...7

zWANumber Read the number of stored parameters.
Number: value from 0...7

zRG Reading active parameter.
 Response: value from 0...7

zWGNumber Write active parameter.
Number: value from 0...7

For example:

HEC-Def 10000 5 5 5 5 600 100 100 30 100 10 1 3 5 8 30 225 1000 1250 1234

1 2 3 4 5 6 7 8 9 A B C D E F G H I J K L

1	Probe name. Max. 11 characters. Do not use space
2	Not used
3	Dead time Rate 1 μ s
4	Dead time Rate 2 μ s
5	Dead time Rate 3 μ s
6	Dead time Rate 4 μ s
7	High voltage in Volt
8	Overload Rate 1 in kcps
9	Overload Rate 3 in kcps
A	Overload probe current μ A
B	Reference area
C	Timeout for detector error in seconds
D	Not used
E	Not used
F	Not used

G	Not used
H	Threshold Rate 1 in mV. Range: 30..1220mV
I	Threshold Rate 2 in mV. Range: 30..1220mV
J	Threshold Rate 3 in mV. Range: 30..1220mV
K	Not used, write “1220”mV
L	Flags

Flags

Bit number	
0	Window:
1	00: No window. Display of Rate 1 Rate 3 01: β window. Display Rate 1 – Rate 2 and Rate 3
2	Pulse fade out 0: off 1: on
3	Not used
4	
5	
6	Reserved, write “0”
7	Reserved, write “0”
8	Single Pulse 0: every Count, 1: with Divide Ratio (see Command ARC)
9..15	Not used

Dead time correction:

$$C1 = \frac{C1_{raw}}{1 - \tau1 * C1_{raw}}$$

C1raw: raw count rate

τ1 Dead time

1- τ1*C1raw is limited to min. 0.1

5.14.10 Event log

6656 520549251

```

|      |
|      \-----
\-----

```

Date and time as a decimal value (see 3.2)

Event log as a decimal value

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	Operation mode 0: Ratemeter 1:Scaler
9	Not used
10	
11	Sound 0: off 1:on
12	LED 0: off 1:on
13	Vibration alarm 0: off 1:on
14	Clear dose 0: no 1:yes
15	1: Alarm threshold changed
16	1: count rate or activity -alarm
17	Not used
18	Not used
19	
20	count rate or activity > alarm threshold 1
21	count rate or activity > alarm threshold 2
22	Not used
23	Not used
24	1: Scaler or Background parameter changed
25	Not used
26	Power off
27	Power on
28	Not used
29	Not used
30	Not used
31	Not used

5.14.11 Automatic sending

```

<STX>721  5 0 0 14 HEC 1234 AB<ETX><CR><LF>
|         | | | | | | | |
|         | | | | | | | \---- BCC in HEX
|         | | | | | | | \----- Dose in µR/h
|         | | | | | | | \----- information RadEye GX
|         | | | | | | | \----- Status as Hex-value
|         | | | | | | | \----- Fixed to "0"
|         | | | | | | | \----- Fixed to "0"
|         | | | | | | | \----- Fixed to "5"
|         | | | | | | | \----- Count rate 1 in 0.01cps
\-----

```

Formation of a BCC (block check character):

Modulo 256 sum of <STX> up to the last character before the BCC (including), coded as hexadecimal ASCII-number (e.g. 1F).

Status

Bit number	
0	Not used
1	1: Overload
2	1: Alarm
3	Not used
4	Not used
5	1: Battery voltage low
6	not used
7	not used

5.14.12 Status information

F Reading status information
 Response Number with status information

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	Not used
9	
10	
11	not used,
12	Alarm threshold (read-only) 0: off 1:on
13	Flag for overload
14	Temperature display 0: off 1:on
15	Not used
16	1: Alarm
17	Not used
18	Not used
19	Not used
20	Value > alarm threshold 1 (depending on display mode)
21	Value > alarm threshold 2 (depending on display mode)
22	Not used
23	Not used
24	Not used
25	Not used
26	Not used
27	Not used
28	Not used
29	Not used
30	Not used
31	Not used

5.14.13 Checksource

IRA	Read number of check sources. Response: number from 0 to 15
IRR	Read repeat time and timestamp of last complete test. Response: <ul style="list-style-type: none">- repeat time in hh from 0 to 65535 (0=off)- date and time of last complete test
IR0..15	Read check source information <ul style="list-style-type: none">- Name of source- Serial of source- Activity of Source (Bq)- Day of activity- Month of activity- Year of activity- Half life (with factor 1000)- Beta efficiency (with factor 100)- Alpha efficiency (with factor 100)- Crossover- Flags- Timestamp of last Test
	<ul style="list-style-type: none">- Beta efficiency (with factor 100)

5.15 RadEye GFex/GF-10ex/Gex/G-10ex

Used firmware version: V1.53

5.15.1 Limit values

AR5	Reading the threshold 1 for dose rate alarm. Response: number in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ unit. Depending on used measuring unit e.g. 1 2 3 means 1 2 3 $\mu\text{R/h}$ resp. 1.23 $\mu\text{Sv/h}$
AR6	Reading the threshold 2 for dose rate alarm. Response: see command AR5
AR7	Reading the threshold 1 for the dose alarm Response: number in μSv , $100\mu\text{R}$ or $100\mu\text{rem}$ units
AR8	Reading the threshold 2 for the dose alarm see command AR7
AW5Number	Setting the threshold 1 for dose rate alarm. Number = in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units. e. g. AW3123 means 123 $\mu\text{R/h}$.
AW6Number	Setting the threshold 2 for dose rate alarm. see command AW3
AW7Number	Setting the threshold 1 for dose alarm. Number = in $100\mu\text{R}$, $100\mu\text{rem}$ or $1\mu\text{Sv}$ units. e. g. AW7123 means 12.3mR.
AW8Number	Setting the threshold 2 for dose alarm. Number = see command AW7

5.15.2 Measurement values

Z	Read raw count rates with dead time correction Response: <ul style="list-style-type: none">- Counter 1 in cps- TTP value- HV power in cps
---	--

5.15.3 Configuration flags

5.15.3.1 Configuration flags 1 with *kR* / *kW*

Bit number	
0	Not used
1	Not used
2	Not used
3	Alarming Sound 0: off 1: on
4	Alarming LED 0: off 1: on
5	Alarming Vibration 0: off 1: on
6	not used, write "0"
7	Single Pulse 0: off 1: on

5.15.3.2 Configuration flags 2 with *fR* / *fW*

Bit number	
0	Disable keylock: 0: no 1:yes
1	
2	Fixed to "1"
3	Alarm threshold read-only 0: off 1:on
4	Flag for overload (readonly)
5	Temp.display 0: off 1:on
6	Temperature unit 0: °C 1: °F
7	Don not silence audio alarm 0: off 1:on

5.15.3.3 Configuration flags 3 with *KR* / *KW*

Bit number	
0	Key tone 0:off 1:on
1	Autosend 0:off 1:on
2	Finder 0:off 1:on
3	Mean-Max 0:off 1:on
4	not used, write "0"
5	display of dose 0:off 1:on
6	Fixed to "1"
7	Not used

5.15.3.4 Configuration flags 4 with *jR* / *jW*

Bit number	
0	not used, write "0"

1	Battery type	0:Alkaline	1:NiMh
2	Display rotation	0:No	1:yes
3	Scaler mode	0: preset count	1: preset time
4	Scaler “Auto restart”	0:off	1:on
5	Net Scaler	0:off	1:on
6	Net Ratemeter	0:off	1:on
7	not used, write “0”		
8	not used, write “0”		
9	not used, write “0”		
10	not used, write “0”		
11	not used, write “0”		
12	not used, write “0”		
13	Display mR/h instead μ R/h	0:No	1:yes
14	not used, write “0”		
15	not used, write “0”		

5.15.3.5 *Measuring unit with uR /uW*

uR Read measuring unit and operation mode

Response: Hex-value. See below

uRHex Write measuring unit and operation mode

See below

Bit number	
0	0x05 Display unit Sv/h 0x06 Display unit R/h 0x07 Display unit rem/h
1	
2	
3	
4	0: Ratemeter, 1: Scaler
5	Last used display unit (dose rate) 0x00 Sv/h 0x01 R/h 0x02 rem/h
6	
7	
7	Not used

5.15.3.6 *Menu configuration*

Bit number			
0	Switch off	0:hidden	1:visible
1	Background	0: hidden	1:visible
2	not used, write "0"		
3	Backlight	0: hidden	1:visible
4	Measuring unit		
5	Operation mode		
6	Scaler parameter		
7	not used, write "0"		
8	not used, write "0"		
9	Alarm Dose Rate	0: hidden	1:visible
10	not used, write "0"		
11	Alarm Dose	0: hidden	1:visible
12	Clear Dose	0: hidden	1:visible
13	Settings	0: hidden	1:visible
14	Alarm indication	0: hidden	1:visible
15	Show alarm	0: hidden	1:visible
16	Text info	0: hidden	1:visible
17	not used, write "0"		
18	not used, write "0"		
19	not used, write "0"		
20	not used, write "0"		
21	not used, write "0"		
22	not used, write "0"		
23	not used, write "0"		
24	not used, write "0"		
25	not used, write "0"		
26	not used, write "0"		
27	not used, write "0"		
28	not used, write "0"		
29	not used, write "0"		
30	not used, write "0"		
31	not used, write "0"		

Submenu Settings:

Bit number			
0	Batt. type	0:hidden	1:visible
1	Autosend	0: hidden	1:visible
2	Single Pulse	0: hidden	1:visible
3	Finder	0: hidden	1:visible
4	Set Date/Time	0: hidden	1:visible
5	Source test	0: hidden	1:visible
6	not used, write "0"		


```

15 11 +19
-- -- ---
|   |   |
|   |   \--- Temperature in °C
|   \----- Rate mean in µR/h, µrem/h or 0.01µSv/h units
\----- Rate max in µR/h, µrem/h or 0.01µSv/h units

```

End of History:

End

5.15.6.2 *History status*

decimal value converted in HEX:

0x40000078

```

|   ----
|   |
|   \-- History cycle time (0x78 -> 120s)
|
\----- 4: identifier for status information
          6: identifier for status information +
            last power on without power off

```


5.15.7 Event log

6656 520549251

```

-----
|          |
|          | \----- Date and time as a decimal value (see 3.2)
|          | \----- Event log as a decimal value
\-----

```

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	Fixed to "0"
9	Fixed to "0"
10	Fixed to "1"
11	Sound 0: off 1:on
12	LED 0: off 1:on
13	Vibration alarm 0: off 1:on
14	1: Dose cleared
15	1: Alarm threshold changed
16	1: Alarm dose rate
17	1: Alarm dose
18	Not used
19	Not used
20	Dose rate > alarm threshold 1
21	Dose rate > alarm threshold 2
22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24	Not used
25	Not used
26	Power off
27	Power on
28	Not used
29	Not used
30	Not used
31	Fixed to "0"

5.15.8 Automatic sending

<STX>7 2 0 0 14 FH41B2 123 AB<ETX><CR><LF>

16	1: Alarm dose rate
17	1: Alarm dose
18	Fixed to "0"
19	Not used
20	Dose rate > alarm threshold 1
21	Dose rate > alarm threshold 2
22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24	Low energy alarm
25	High energy alarm
26	Not used
27	Not used
28	Not used
29	Not used
30	Not used
31	Not used

5.16 RadEye G/G-10/B20/G20/G20-10/GF/GF-10 (from V2.00)

Used firmware version: V2.06

5.16.1 Limit values

AR1	Reading the threshold 1 for activity alarm. Response: number in 0.01 Bq units.
AR2	Reading the threshold 2 for activity alarm. Response: number in 0.01 Bq units.
AR3	Reading the threshold 1 for count rate alarm. Response: number in 0.01 cps units.
AR4	Reading the threshold 2 for count rate alarm. Response: number in 0.01 cps units.
AR5	Reading the threshold 1 dose rate alarm. Response: number in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units
AR6	Reading the threshold 2 for dose rate alarm. Response: number in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units
AR7	Reading the threshold 1 for dose alarm. Response: number in μR , μrem or $0.01\mu\text{Sv}$ units
AR8	Reading the threshold 2 for dose alarm. Response number in μR , μrem or $0.01\mu\text{Sv}$ units
ARP	Reading scaler parameter: <ul style="list-style-type: none">- Preset count. Response: Number in counts- Preset time. Response: Number in seconds
ARN	Reading background parameter: <ul style="list-style-type: none">- Preset count. Response: Number in counts- Preset time. Response: Number in seconds
ARB	Reading background value: Response: Value in 0.01cps
AW1Number	Setting the threshold 1 for activity alarm Number: value in 0.01 Bq units.
AW2Number	Setting the threshold 2 for activity alarm Number: value in 0.01 Bq units

<i>AW3Number</i>	Setting the threshold 1 for count rate alarm <i>Number</i> : value in 0.01 cps units.
<i>AW4Number</i>	Setting the threshold 2 for count rate alarm <i>Number</i> : value in 0.01 cps units
<i>AW5Number</i>	Setting the threshold 1 dose rate alarm. <i>Number</i> : value in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units.
<i>AW6Number</i>	Setting the threshold 2 dose rate alarm. <i>Number</i> : value in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units.
<i>AW7Number</i>	Setting the threshold 1 for dose alarm. <i>Number</i> : value in μR , μrem or $0.01\mu\text{Sv}$ units
<i>AW8Number</i>	Setting the threshold 2 for dose alarm. <i>Number</i> : value in μR , μrem or $0.01\mu\text{Sv}$ units
<i>AWPNumber Number</i>	Set scaler parameter: <i>Number</i> : Preset count in counts (from 0 to 9999 counts) <i>Number</i> : Preset time in seconds (from 0 to 9999 seconds)
<i>AWNNumber Number</i>	Set background parameter: <i>Number</i> : Preset count in counts (from 0 to 9999 counts) <i>Number</i> : Preset time in seconds (from 0 to 9999 seconds)
<i>AWBNumber Number</i>	Set background value: <i>Number</i> : Value in 0.01 cps (from 0 to 100 cps)

5.16.2 Measurement values

Z	Read raw count rates with dead time correction Response: - Counter 1 in cps - HV power in cps
z	Read filtered count rate Response: Number in 0.01 cps units
A	Read display value and status Response: - display value Dose rate in $1\mu\text{R/h}$, $1\mu\text{rem/h}$, $0.01\mu\text{Sv/h}$ units or 0.01 cps, cpm, Bq, dps, dpm, Bq/cm ² units - Status (see 5.16.6.3)

5.16.3 Configuration flags

5.16.3.1 Configuration flags 1 with *kR* / *kW*

Bit number	
0	Not used
1	Not used
2	Not used
3	Alarming Sound 0: off 1: on
4	Alarming LED 0: off 1: on
5	Alarming Vibration 0: off 1: on
6	Not used
7	Single Pulse 0: off 1: on

5.16.3.2 Configuration flags 2 with *fR* / *fW*

Bit number	
0	Not used
1	Not used
2	Not used, write “1”
3	Alarm threshold read-only 0: off 1: on
4	Flag for overload (readonly)
5	Display of temperature 0: off 1: on
6	Temperature unit 0: °C 1: °F
7	Not used

5.16.3.3 Configuration flags 3 with KR / KW

Bit number	
0	Key tone 0:off 1:on
1	Autosend 0:off 1:on
2	Finder 0:off 1:on
3	Mean-Max 0:off 1:on
4	not used, write "0"
5	display of dose 0:off 1:on
6	not used
7	not used

5.16.3.4 Configuration flags 4 with jR / jW

Bit number	
0	Not used
1	Battery type 0:Alkaline 1:NiMh
2	Display rotation 0:No 1:yes
3	Scaler mode 0: Preset count 1: Preset time
4	Not used
5	Scaler, after measurement 0: Stop 1: automatic restart
6	Scaler net 0:No 1:yes
7	Ratemeter net 0:No 1:yes
8	Not used
9	Not used
10	Not used
11	Not used
12	Not used
13	Not used
14	Not used
15	Not used

5.16.3.5 *Measuring unit with uR /uW*

uR Read measuring unit and operation mode
Response: Hex-value. See below
uRHex Write measuring unit and operation mode
See below

Bit number	
0	0x00 Display unit cps 0x04 Display unit dpm
1	0x01 Display unit cpm 0x05 Display unit Sv/h
2	0x02 Display unit Bq 0x06 Display unit R/h
3	0x03 Display unit dps 0x07 Display unit rem/h 0x08 Display unit Bq/cm ²
4	0: Ratemeter, 1: Scaler
5	Last used display unit (dose rate)
6	0x00 Sv/h 0x01 R/h 0x02 rem/h
7	Not used

5.16.3.6 *Menu configuration*

mR Read configuration for

- main menu,
- submenu “Settings”
- submenu “Alarm indication”

Response: Hex-values. See below

mRHex Hex Write configuration for

- main menu,
- submenu “Settings”
- submenu “Alarm indication”

See below

Bit number			
0	Switch off	0:hidden	1:visible
1	Background	0: hidden	1:visible
2	not used, write “0”		
3	Backlight	0: hidden	1:visible
4	Measuring unit	0: hidden	1:visible
5	Operation mode	0: hidden	1:visible
6	Scaler parameter	0: hidden	1:visible
7	Nuclide table	0: hidden	1:visible
8	Alarm count rate	0: hidden	1:visible
9	Alarm dose rate	0: hidden	1:visible
10	Alarm activity	0: hidden	1:visible
11	Alarm dose	0: hidden	1:visible
12	Clear dose	0: hidden	1:visible
13	Settings	0: hidden	1:visible
14	Alarm indication	0: hidden	1:visible
15	Show alarm	0: hidden	1:visible
16	Text Info	0: hidden	1:visible
17	not used, write “0”		
18	not used, write “0”		
19	not used, write “0”		
20	not used, write “0”		
21	not used, write “0”		
22	not used, write “0”		
23	not used, write “0”		
24	not used, write “0”		
25	not used, write “0”		
26	not used, write “0”		
27	not used, write “0”		
28	not used, write “0”		
29	not used, write “0”		

30	not used, write “0”
31	not used, write “0”

Bit number			
0	Batt. type	0:hidden	1:visible
1	Autosend	0: hidden	1:visible
2	Single Pulse	0: hidden	1:visible
3	Finder	0: hidden	1:visible
4	Set Date/Time	0: hidden	1:visible
5	not used, write “0”		
6	not used, write “0”		
7	not used, write “0”		
8	not used, write “0”		
9	not used, write “0”		
10	not used, write “0”		
11	not used, write “0”		
12	not used, write “0”		
13	not used, write “0”		
14	not used, write “0”		
15	not used, write “0”		

Bit number			
0	Sound	0:hidden	1:visible
1	LED	0: hidden	1:visible
2	Vibrator	0: hidden	1:visible
3	Finder	0: hidden	1:visible
4	not used, write “0”		
5	not used, write “0”		
6	not used, write “0”		
7	not used, write “0”		
8	not used, write “0”		

5.16.4 High voltage

HR Reading high voltage bit value
 Response: Bit value 0...255.

5.16.5 Dead time correction

x Read dead time
 Response: dead time in ns

5.16.6 History output

5.16.6.1 History readout ratemeter

```
256 716612088 721 999 120 23
---  -
|    |          |    |    |    |
|    |          |    |    |    \--- Temperature in °C
|    |          |    |    \----- Measuring time in seconds
|    |          |    \----- Max value
|    |          |          Count rate in 0.01 cps, cpm
|    |          |          Activity in 0.01 Bq, dps, dpm, Bq/cm²
|    |          |          Dose rate µR/h, rem/h or 0.01µSv/h
|    |          \----- Mean value
|    |          |          Count rate in 0.01 cps, cpm
|    |          |          Activity in 0.01 Bq, dps, dpm, Bq/cm²
|    |          |          Dose rate µR/h, rem/h or 0.01µSv/h
|    \----- Date and Time as a decimal value
\----- History status as a decimal value
           see below
```

5.16.6.2 History readout scaler

```
256 716612088 721 999 120 23
---  -
|    |          |    |    |    |
|    |          |    |    |    \--- Temperature in °C
|    |          |    |    \----- Measuring time in seconds
|    |          |    \----- Background in 0.01 cps
|    |          \----- Mean value
|    |          |          Count rate in 0.01 cps, cpm
|    |          |          Activity in 0.01 Bq, dps, dpm, Bq/cm²
|    |          |          Dose rate µR/h, rem/h or 0.01µSv/h
|    \----- Date and Time as a decimal value
\----- History status as a decimal value
           see below
```

End of History:

End

5.16.6.3 *History status*

Decimal value converted in HEX:

Bit number			
0	Net value	0:No	1:Yes
1	Operation mode	0: Ratemeter	1:Scaler
2	Reserved. For internal use only.		
3	Accumulated counts	0:No	1:Yes
4	Background measurement	0:No	1:Yes
5	Scaler with	0: Preset Counts	1: Preset Time
6	Used filter (B20 and B20-ER only)		
7	0: No filter 1: Alpha blocker, 2:H*(10), 3:Hx		
8	Number of used nuclide		
9			
10			
11			
12	0x00 Display unit cps	0x04 Display unit dpm	
13	0x01 Display unit cpm	0x05 Display unit Sv/h	
14	0x02 Display unit Bq	0x06 Display unit R/h	
15	0x03 Display unit dps	0x07 Display unit rem/h	
		0x08 Display unit Bq/cm ²	

5.16.7 Nuclide table

nRNumber Reading nuclide data.
Number: consecutive number
 Response: Nuclide data (see below)

nWNumberString Write nuclide data
Number: consecutive number
String: nuclide data
 e.g. nW02Sr-90 500 555 3333

nRA Reading number of stored nuclides.
 Response: value from 0...15

nWANumber Write the number of stored nuclides.
Number: value from 0...15

nRG Reading active nuclide.
 Response: value from 0...15

nWANumber Write the number of active nuclides.
Number: value from 0...15

Nuclide data:

Sr-90 500 555 3333

```

-----
|      |      |      |
|      |      |      \----- Factor for activity calculation with gamma filter
|      |      \----- Factor for activity calculation with alpha blocker
|      \----- Factor for activity calculation without filter
\----- Nuclide name. Up to 6 characters. Do not use space (0x20)
```

Activity calculation:

$F = 1/E_{eff}$

$Ar = Cr * F$

F: Factor for Efficiency calculation

E_{eff} Efficiency for this nuclide

Ar Activity

Cr Measured count rate

For example:

Count rate is 67.3 cps without filter and efficiency for Sr-90 is

- 29% without filter.
- 25% with alpha blocker
- 10% with gamma filter

Factor F is $1/0.29=3.448$.

- $1/0.29 = 3.448$ without filter.
- $1/0.25 = 4.0$ with alpha blocker
- $1/0.10 = 10.0$ with gamma filter

To set the parameter, sent: `nW00Sr-90 345 400 1000`.

The activity is $67.3 \times 3.45 = 232.18$ Bq

5.16.8 Event log

6656 520549251

```

-----
|          |
|          | \----- Date and time as a decimal value (see 3.2)
|          | \----- Event log as a decimal value
\-----

```

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	Not used
9	
10	
11	Sound 0: off 1:on
12	LED 0: off 1:on
13	Vibration alarm 0: off 1:on
14	Clear dose 0: no 1:yes
15	1: Alarm threshold changed
16	1: Dose rate, count rate or activity -alarm
17	1: Alarm dose
18	Not used
19	Not used
20	Dose rate, count rate or activity > alarm threshold 1
21	Dose rate, count rate or activity > alarm threshold 2
22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24	1: Scaler or Background parameter changed
25	Not used
26	Power off
27	Power on
28	Not used
29	Not used
30	Not used
31	Not used

5.16.9 Automatic sending

```

<STX>721  2 0 0 14 B20ER 1234 AB<ETX><CR><LF>
|         | | | | |         |         |
|         | | | | |         |         | \- BCC in HEX
|         | | | | |         |         | \----- Dose in µR/h
|         | | | | |         |         | \----- information RadEye *)
|         | | | | |         |         | \----- Status as Hex-value
|         | | | | |         |         | \----- Fixed to "0"
|         | | | | |         |         | \----- Fixed to "0"
|         | | | | |         |         | \----- Measurement unit **)
|         | | | | |         |         | \----- Measurement value

```

*) Information RadEye:

B20	RadEye B20
B20ER	RadEye B20-ER
G20	RadEye G20
G20ER	RadEye G20-ER
G2010	RadEye G20-10
G20ER1	RadEye G20-ER10
GF	RadEye GF
GF10	RadEye GF-10
FH41B2	RadEye G or RadEye G-10

**) Measurement unit

5	Dimension cps
3	Dimension cpm
7	Dimension dpm
8	Dimension Bq/cm ²
6	Dimension Bq
9	Dimension dps
0	Dimension Sv/h
10	Dimension rem/h
2	Dimension R/h

Formation of a BCC (block check character):

Modulo 256 sum of <STX> up to the last character before the BCC (including), coded as hexadecimal ASCII-number (e.g. 1F).

Status

Bit number	
0	Not used
1	1: Overload
2	1: Alarm
3	Not used
4	Not used
5	1: Battery voltage low
6	not used
7	not used

5.16.10 Status information

F Reading status information
 Response Number with status information

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	Not used
9	
10	
11	not used,
12	Alarm threshold read-only 0: off 1:on
13	Flag for overload
14	Temperature display 0: off 1:on
15	Not used
16	1: Alarm
17	Not used
18	Not used
19	Not used
20	Value > alarm threshold 1 (depending on display mode)
21	Value > alarm threshold 2 (depending on display mode)
22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24	Not used
25	Not used
26	Not used
27	Not used
28	Not used
29	Not used
30	Not used
31	Not used

5.17 RadEye G/G-10/B20/G20/G20-10/GF/GF-10 (from V3.00)

Used firmware version: V3.06

5.17.1 Limit values

AR1	Reading the threshold 1 for activity alarm. Response: number in 0.01 Bq units.
AR2	Reading the threshold 2 for activity alarm. Response: number in 0.01 Bq units.
AR3	Reading the threshold 1 for count rate alarm. Response: number in 0.01 cps units.
AR4	Reading the threshold 2 for count rate alarm. Response: number in 0.01 cps units.
AR5	Reading the threshold 1 dose rate alarm. Response: number in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units
AR6	Reading the threshold 2 for dose rate alarm. Response: number in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units
AR7	Reading the threshold 1 for dose alarm. Response: number in μR , μrem or $0.01\mu\text{Sv}$ units
AR8	Reading the threshold 2 for dose alarm. Response number in μR , μrem or $0.01\mu\text{Sv}$ units
ARP	Reading scaler parameter: <ul style="list-style-type: none">- Preset count. Response: Number in counts- Preset time. Response: Number in seconds
ARN	Reading background parameter: <ul style="list-style-type: none">- Preset count. Response: Number in counts- Preset time. Response: Number in seconds
ARB	Reading background value: Response: Value in 0.01cps
ARt	Reading Tau (from V3.06): Response: Value in 0.1 s
AW1Number	Setting the threshold 1 for activity alarm Number: value in 0.01 Bq units.
AW2Number	Setting the threshold 2 for activity alarm Number: value in 0.01 Bq units

<i>AW3Number</i>	Setting the threshold 1 for count rate alarm <i>Number</i> : value in 0.01 cps units.
<i>AW4Number</i>	Setting the threshold 2 for count rate alarm <i>Number</i> : value in 0.01 cps units
<i>AW5Number</i>	Setting the threshold 1 dose rate alarm. <i>Number</i> : value in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units.
<i>AW6Number</i>	Setting the threshold 2 dose rate alarm. <i>Number</i> : value in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units.
<i>AW7Number</i>	Setting the threshold 1 for dose alarm. <i>Number</i> : value in μR , μrem or $0.01\mu\text{Sv}$ units
<i>AW8Number</i>	Setting the threshold 2 for dose alarm. <i>Number</i> : value in μR , μrem or $0.01\mu\text{Sv}$ units
<i>AWPNumber Number</i>	Set scaler parameter: <i>Number</i> : Preset count in counts (from 0 to 9999 counts) <i>Number</i> : Preset time in seconds (from 0 to 9999 seconds)
<i>AWNNumber Number</i>	Set background parameter: <i>Number</i> : Preset count in counts (from 0 to 9999 counts) <i>Number</i> : Preset time in seconds (from 0 to 9999 seconds)
<i>AWBNumber Number</i>	Set background value: <i>Number</i> : Value in 0.01 cps (from 0 to 100 cps)
<i>AWtNumber</i>	Set Tau: <i>Number</i> : Value in 0.1 s (from 10 to 600 means from 1.0 to 60.0s)

5.17.2 Measurement values

<i>Z</i>	Read raw count rates with dead time correction Response: - Counter 1 in cps - TTP-Value - HV power in cps
<i>z</i>	Read filtered count rate Response: Number in 0.01 cps units

- A Read display value and status
 Response:
 - display value Dose rate in 1μR/h, 1μrem/h, 0.01μSv/h units
 or 0.01 cps, cpm, Bq, dps, dpm, Bq/cm² units
 - Status (see 5.17.6.3)

5.17.3 Configuration flags

5.17.3.1 *Configuration flags 1 with kR / kW*

Bit number			
0	Not used		
1	Live graph (from V3.06)	0: disable	1: enable
2	Not used		
3	Alarming Sound	0: off	1: on
4	Alarming LED	0: off	1: on
5	Alarming Vibration	0: off	1: on
6	Not used		
7	Single Pulse	0: off	1: on

5.17.3.2 *Configuration flags 2 with fR / fW*

Bit number			
0	Keylock(from V3.06)	0: enable, 1:disable	
1	Not used		
2	Not used, write “1”		
3	Alarm threshold read-only	0: off	1:on
4	Flag for overload (readonly)		
5	Display of temperature	0: off	1:on
6	Temperature unit	0: °C	1: °F
7	not used		

5.17.3.3 *Configuration flags 3 with KR / KW*

Bit number			
0	Key tone	0:off	1:on
1	Autosend	0:off	1:on
2	Finder	0:off	1:on
3	Mean-Max	0:off	1:on
4	not used, write “0”		
5	display of dose	0:off	1:on
6	not used		

7	not used
---	----------

5.17.3.4 *Configuration flags 4 with jR / jW*

Bit number	
0	Not used
1	Battery type 0:Alkaline 1:NiMh
2	Display rotation 0:No 1:yes
3	Scaler mode 0: Preset count 1: Preset time
4	Accumulated counts 0:No 1:yes
5	Scaler, after measurement 0: Stop 1: automatic restart
6	Scaler net 0:No 1:yes
7	Ratemeter net 0:No 1:yes
8	Not used
9	Not used
10	Not used
11	Alarm LCD LED 0:No 1:yes
12	BTCombo: Remote power off 0:No 1:yes
13	Display mR/h instead of μ R/h 0:No 1:yes
14	Not used
15	Not used

5.17.3.5 *Measuring unit with uR /uW*

uR Read measuring unit and operation mode
 Response: Hex-value. See below
 uRHex Write measuring unit and operation mode
 See below

Bit number	
0	0x00 Display unit cps 0x05 Display unit Sv/h
1	0x01 Display unit cpm 0x06 Display unit R/h
2	0x02 Display unit Bq 0x07 Display unit rem/h
3	0x03 Display unit dps 0x08 Display unit Bq/cm ² 0x04 Display unit dpm 0x09 Display unit Gy/h
4	0: Ratemeter, 1: Scaler
5	Last used display unit (dose rate)
6	00 Gy/h 10 R/h 01 Sv/h 11 rem/h
7	Ratemeter Mode: 1:ADF 0:Tau (only B20/B20-ER)

5.17.3.6 *Menu configuration*

mR Read configuration for

- main menu,
- submenu “Settings”
- submenu “Alarm indication”

Response: Hex-values. See below

mRHex Hex Write configuration for

- main menu,
- submenu “Settings”
- submenu “Alarm indication”

See below

5.17.3.6.1 Main menu RadEye G/G10

Bit number			
0	Switch off	0:hidden	1:visible
1	Sound	0: hidden	1:visible
2	LED	0: hidden	1:visible
3	Vibrator	0: hidden	1:visible
4	not used, write “0”		
5	not used, write “0”		
6	not used, write “0”		
7	not used, write “0”		
8	Alarm dose rate	0: hidden	1:visible
9	Alarm dose	0: hidden	1:visible
10	not used, write “0”		
11	Autosend	0: hidden	1:visible
12	not used, write “0”		
13	Clear dose	0: hidden	1:visible
14	Finder	0: hidden	1:visible
15	Single Pulse	0: hidden	1:visible
16	Backlight	0: hidden	1:visible
17	Show Alarm	0: hidden	1:visible
18	Settings	0: hidden	1:visible
19	Text Info	0: hidden	1:visible
20	Bluetooth	0: hidden	1:visible
21..31	not used, write “0”		

5.17.3.6.2 Main menu RadEye B20../G20../GF..

Bit number			
0	Switch off	0:hidden	1:visible
1	Background	0: hidden	1:visible
2	not used, write "0"		
3	Backlight	0: hidden	1:visible
4	Measuring unit	0: hidden	1:visible
5	Operation mode	0: hidden	1:visible
6	Scaler parameter	0: hidden	1:visible
7	Nuclide table	0: hidden	1:visible
8	Alarm count rate	0: hidden	1:visible
9	Alarm dose rate	0: hidden	1:visible
10	Alarm activity	0: hidden	1:visible
11	Alarm dose	0: hidden	1:visible
12	Clear dose	0: hidden	1:visible
13	Settings	0: hidden	1:visible
14	Alarm indication	0: hidden	1:visible
15	Show alarm	0: hidden	1:visible
16	Text Info	0: hidden	1:visible
17..31	not used, write "0"		

5.17.3.6.3 Submenu "Settings" RadEye B20../G20../GF..

Bit number			
0	Batt. type	0:hidden	1:visible
1	Autosend	0: hidden	1:visible
2	Single Pulse	0: hidden	1:visible
3	Finder	0: hidden	1:visible
4	Set Date/Time	0: hidden	1:visible
5	Source test	0: hidden	1:visible
6	Language	0: hidden	1:visible
7	Contrast	0: hidden	1:visible
8	Edit Tau	0: hidden	1:visible
8..15	not used, write "0"		

5.17.3.6.4 Submenu "Settings" RadEye G/G10

Bit number			
0	Batt. type	0:hidden	1:visible
1	Set Date/Time	0: hidden	1:visible
2	Source check	0: hidden	1:visible
3	Language	0: hidden	1:visible

4	Contrast	0: hidden	1:visible
5..15	not used, write "0"		

5.17.3.6.5 Submenu "Alarm Indication" RadEye B20../G20../GF..

Bit number			
0	Sound	0:hidden	1:visible
1	LED	0: hidden	1:visible
2	Vibrator	0: hidden	1:visible
3	not used, write "0"		
4	not used, write "0"		
5	not used, write "0"		
6	not used, write "0"		
7	not used, write "0"		
8	not used, write "0"		

5.17.4 High voltage

HR Reading high voltage bit value

Response: Bit value 0...255.

5.17.5 Dead time correction

x Read dead time

Response: dead time in ns

5.17.6 History output

5.17.6.1 History readout ratemeter

256 716612088 721 999 120 23

```

---  ---  ---  ---  ---  ---
|      |      |      |      |      |
|      |      |      |      | \--- Temperature in °C
|      |      |      |      | \----- Measuring time in seconds
|      |      |      | \----- Max value
|      |      |      |      Count rate in 0.01 cps, cpm
|      |      |      |      Activity in 0.01 Bq, dps, dpm, Bq/cm²
|      |      |      |      Dose rate µR/h, rem/h or 0.01µSv/h
|      |      | \----- Mean value
|      |      |      Count rate in 0.01 cps, cpm
|      |      |      Activity in 0.01 Bq, dps, dpm, Bq/cm²
|      |      |      Dose rate µR/h, rem/h or 0.01µSv/h
|      | \----- Date and Time as a decimal value
| \----- History status as a decimal value
| see below

```

5.17.6.2 History readout scaler

256 716612088 721 999 120 23

```

-----
|      |      |      |      |      |
|      |      |      |      | \--- Temperature in °C
|      |      |      | \----- Measuring time in seconds
|      |      | \----- Background in 0.01 cps
|      | \----- Mean value
|      |      Count rate in 0.01 cps, cpm
|      |      Activity in 0.01 Bq, dps, dpm, Bq/cm²
|      |      Dose rate µR/h, rem/h or 0.01µSv/h
| \----- Date and Time as a decimal value
\----- History status as a decimal value
      see below

```

End of History:

End

5.17.6.3 History status

Decimal value converted in HEX:

Bit number		
0	Net value	0:No 1:Yes
1	Operation mode	0: Ratemeter 1:Scaler
2	Reserved. For internal use only.	
3	Accumulated counts	0:No 1:Yes
4	Background measurement	0:No 1:Yes
5	Scaler with	0: Preset Counts 1: Preset Time
6	Used filter (B20 and B20-ER only)	
7	0: No filter 1: Alpha blocker, 2:H*(10), 3:Hx	
8	Number of used nuclide	
9		
10		
11		
12	0x00 Display unit cps	0x05 Display unit Sv/h
13	0x01 Display unit cpm	0x06 Display unit R/h
14	0x02 Display unit Bq	0x07 Display unit rem/h
15	0x03 Display unit dps	0x08 Display unit Bq/cm²
	0x04 Display unit dpm	0x09 Display unit Gy/h

5.17.7 Nuclide table

nRNumber Reading nuclide data.
Number: consecutive number
 Response: Nuclide data (see below)

nWNumberString Write nuclide data
Number: consecutive number
String: nuclide data
 e.g. nW02Sr-90 500 555 3333

nRA Reading number of stored nuclides.
 Response: value from 0...15

nWANumber Write the number of stored nuclides.
Number: value from 0...15

nRG Reading active nuclide.
 Response: value from 0...15

nWANumber Write the number of active nuclides.
Number: value from 0...15

Nuclide data:

Sr-90 500 555 3333

```

-----
|      |      |      |
|      |      |      \----- Factor for activity calculation with gamma filter
|      |      \----- Factor for activity calculation with alpha blocker
|      \----- Factor for activity calculation without filter
\----- Nuclide name. Up to 6 characters. Do not use space (0x20)

```

Activity calculation:

$F = 1/E_{eff}$

$A_r = C_r * F$

F: Factor for Efficiency calculation

E_{eff} Efficiency for this nuclide

A_r Activity

C_r Measured count rate

For example:

Count rate is 67.3 cps without filter and efficiency for Sr-90 is

- 29% without filter.
- 25% with alpha blocker

- 10% with gamma filter

Factor F is $1/0.29=3.448$.

- $1/0.29 = 3.448$ without filter.
- $1/0.25 = 4.0$ with alpha blocker
- $1/0.10 = 10.0$ with gamma filter

To set the parameter, sent: nW00Sr-90 345 400 1000.

The activity is $67.3 \times 3.45 = 232.18$ Bq

5.17.8 Event log

6656 520549251

```

-----
|          |
|          | \----- Date and time as a decimal value (see 3.2)
|          | \----- Event log as a decimal value
\-----

```

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	Not used
9	
10	
11	Sound 0: off 1:on
12	LED 0: off 1:on
13	Vibration alarm 0: off 1:on
14	Clear dose 0: no 1:yes
15	1: Alarm threshold changed
16	1: Dose rate, count rate or activity -alarm
17	1: Alarm dose
18	Not used
19	Not used
20	Dose rate, count rate or activity > alarm threshold 1
21	Dose rate, count rate or activity > alarm threshold 2
22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24	1: Scaler or Background parameter changed
25	Not used
26	Power off
27	Power on
28	Not used
29	Not used
30	Not used
31	Not used

5.17.9 Automatic sending

```

<STX>721  2 0 0 14 B20ER 1234 AB<ETX><CR><LF>
|         | | | | |         |         |         |
|         | | | | |         |         |         \- BCC in HEX
|         | | | | |         |         \----- Dose in µR/h
|         | | | | |         \----- information RadEye *)
|         | | | | \----- Status as Hex-value
|         | | | \----- Fixed to "0"
|         | | \----- Fixed to "0"
|         | \----- Measurement unit **)
|         \----- Measurement value

```

*) Information RadEye:

B20	RadEye B20
B20ER	RadEye B20-ER
G20	RadEye G20
G20ER	RadEye G20-ER
G2010	RadEye G20-10
G20ER1	RadEye G20-ER10
GF	RadEye GF
GF10	RadEye GF-10
FH41B2	RadEye G or RadEye G-10

**) Measurement unit

5	Dimension cps
3	Dimension cpm
7	Dimension dpm
8	Dimension Bq/cm ²
6	Dimension Bq
9	Dimension dps
0	Dimension Sv/h
10	Dimension rem/h
2	Dimension R/h
1	Dimension Gy/h

Formation of a BCC (block check character):

Modulo 256 sum of <STX> up to the last character before the BCC (including), coded as hexadecimal ASCII-number (e.g. 1F).

Status

Bit number	
0	Not used
1	1: Overload
2	1: Alarm
3	Not used
4	Not used
5	1: Battery voltage low
6	not used
7	not used

5.17.10 Status information

F Reading status information
 Response Number with status information

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	Not used
9	
10	
11	not used,
12	Alarm threshold read-only 0: off 1:on
13	Flag for overload
14	Temperature display 0: off 1:on
15	Not used
16	1: Alarm
17	Not used
18	Not used
19	Not used
20	Value > alarm threshold 1 (depending on display mode)
21	Value > alarm threshold 2 (depending on display mode)
22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24..31	Not used

5.18 RadEye GN/GN+

Used firmware version: 3.02

5.18.1 Limit values

AR0	Reading parameters for sigma threshold. Response: <ul style="list-style-type: none">- Sigma value gamma 0 (=off), 2...9.- '0'- Sigma value neutron 0 (=off), 2...9- '5'
AR1	Reading the threshold 1 and 2 for alarm count rate gamma. Response: number in 0.01 cps
AR2	Reading the threshold 1 and 2 for alarm count rate neutron. Response: number in 0.01 cps.
AR3	Reading the threshold 1 and 2 for alarm dose rate gamma. Response: number in R/h or 0.01μSv/h units.
AR4	Reading the threshold 1 and 2 for alarm dose gamma. Response: number in R or 0.01μSv units.
AR5	Reading the threshold 1 and 2 for alarm dose rate neutron. Response: number in R/h or 0.01μSv/h units.
AR6	Reading the threshold 1 and 2 for alarm dose neutron. Response: number in R or 0.01μSv units.
ART	Read alarm timeout Response: number in seconds
SR3	Read NBR alarm threshold <ul style="list-style-type: none">• LOW Energy• MID Energy• HIGH Energy• NBR Alarm min. count rate

<i>AW0Number Number</i>	Set parameters for sigma alarm. <ul style="list-style-type: none"> - Sigma value gamma 0 (=off), 2...9. - '0' - Sigma value neutron 0 (=off), 2...9 - '5'
<i>AW1Number Number</i>	Set the threshold 1 and 2 for alarm count rate gamma. <i>Number</i> : number in 0.01 cps
<i>AW2Number Number</i>	Set the threshold 1 and 2 for alarm count rate neutron. <i>Number</i> : number in 0.01 cps.
<i>AW3Number Number</i>	Set the threshold 1 and 2 for alarm dose rate gamma. <i>Number</i> : number in R/h or 0.01μSv/h units.
<i>AW4Number Number</i>	Set the threshold 1 and 2 for alarm dose gamma. <i>Number</i> : number in R or 0.01μSv units.
<i>AW5Number Number</i>	Set the threshold 1 and 2 for alarm dose rate neutron. <i>Number</i> : number in R/h or 0.01μSv/h units.
<i>AW6Number Number</i>	Set the threshold 1 and 2 for alarm dose neutron. <i>Number</i> : number in R or 0.01μSv units.
<i>SW3 Number... Number</i>	Set NBR alarm threshold <ul style="list-style-type: none"> • LOW Energy • MID Energy • HIGH Energy • NBR Alarm min. count rate

5.18.2 Measurement values

<i>Z</i>	Read raw count rates with dead time correction Response: <ul style="list-style-type: none"> - Counter 1 in cps - Counter 2 in cps - Counter 3 in cps - Counter 4 in cps - Counter 5 in cps - Counter 6 in cps - HV power index
<i>z</i>	Read filtered count rates Response: <ul style="list-style-type: none"> - Value gamma 0.01 cps

- Value neutron in 0.01 cps

A Read display value and status

Response:

- display value gamma in 0.01 cps, R/h, rem/h or 0.01μSv/h
- display value neutron in 0.01 cps, , R/h, rem/h or 0.01μSv/h
- Status (see 5.18.5.2)

5.18.3 Configuration flags

5.18.3.1 Configuration flags 1 with *kR* / *kW*

Bit number	
0	Not used
1	Not used
2	Alarm LCD-LED 0: off 1: on
3	Alarming Sound 0: off 1: on
4	Alarming LED 0: off 1: on
5	Alarming Vibration 0: off 1: on
6	NBR 0: off 1: on
7	Single Pulse 0: off 1: on

5.18.3.2 Configuration flags 2 with *fR* / *fW*

Bit number	
0	Keylock 0: disable 1:enable
1	Not used
2	Not used
3	Alarm threshold read-only 0: off 1:on
4	Flag for overload (read-only)
5	Display of temperature 0: off 1:on
6	Temperature unit 0: °C 1: °F
7	Not used

5.18.3.3 Configuration flags 3 with *KR* / *KW*

Bit number	
0	Key tone 0:off 1:on
1	Autosend 0:off 1:on
2	Finder 0:off 1:on
3	Mean-Max 0:off 1:on
4	not used, write "0"
5	display of dose 0:off 1:on
6	not used

7	not used
---	----------

5.18.3.4 Configuration flags 4 with jR / jW

Bit number	
0	Safety alarm 0: disable 1:enable
1	Battery type 0:Alkaline 1:NiMh
2	Display rotation 0:No 1:yes
3	Single pulse click divider 0: disable 1:enable
4	Not used
5	Not used
6	Not used
7	Not used
8	Not used
9	Not used
10	Not used
11	Not used
12	Not used
13	Not used
14	Not used
15	Not used

5.18.3.5 Measuring unit with uR /uW

uR Read measuring unit and operation mode

Response: Hex-value. See below

uRHex Write measuring unit and operation mode

See below

Gamma

Bit number	
0	0x00 Display unit cps 0x05 Display unit Sv/h 0x01 Display unit cpm 0x06 Display unit R/h 0x07 Display unit rem/h
1	
2	
3	
4	Not used
5	Dose rate: 01: Sv/h 10: R/h 11: rem/h
6	
7	
7	Cross mode

Neutron

Bit number	
0	0x00 Display unit cps 0x05 Display unit Sv/h
1	0x01 Display unit cpm 0x06 Display unit R/h
2	0x07 Display unit rem/h
3	
4	Not used
5	Display mode:
6	00: Gamma 01: Neutron 02: dual display
7	Cross mode

5.18.3.6 Menu configuration

mR Read configuration for main menu, submenu “Settings” and “Alarm indication”
 Response: Hex-values. See below
 mRHex... Hex Write configuration for main menu, submenu “Settings” and “Alarm indication”
 See below

Bit number	
0	Switch off 0:hidden 1:visible
1	Background 0: hidden 1:visible
2	Backlight 0: hidden 1:visible
3	Measuring unit Gamma 0: hidden 1:visible
4	Measuring unit Neutron 0: hidden 1:visible
5	Operation mode 0: hidden 1:visible
6	Scaler parameter 0: hidden 1:visible
7	Set DR factor 0: hidden 1:visible
8	Edit alarms 0: hidden 1:visible
9	Clear dose Gamma 0: hidden 1:visible
10	Clear dose Neutron 0: hidden 1:visible
11	Settings 0: hidden 1:visible
12	Alarm indication 0: hidden 1:visible
13	Show alarm 0: hidden 1:visible
14	Text Info 0: hidden 1:visible
15	Training 0: hidden 1:visible
16..31	not used, write “0”

Sub menu “Settings”

Bit number	
0	Batt. type 0:hidden 1:visible
1	Autosend 0: hidden 1:visible
2	Single Pulse 0: hidden 1:visible
3	Finder 0: hidden 1:visible

End of History:
End

Decimal value converted in HEX:

```
0101 0101      0101 0101      0101 0101
||||| |||||    ||||| |||||    ||||| |||||
||||| |||||    ||||| |||||    ||||| |||\-- Background sub, 0:No, 1:Yes
||||| |||||    ||||| |||||    ||||| |||\--- 1: Scaler
||||| |||||    ||||| |||||    |||++ ++----
||||| |||||    ||||| |||||    |\----- 1:Background measurement
||||| |||||    ||||| +++++----- Measuring unit Gamma:
|||||          0:cps,1:cpm,5:Sv/h, 6;R/h, 7:rem/h
|||||          +++++----- Measuring unit Neutronen:
|||||          0:cps,1:cpm,5:Sv/h, 6;R/h, 7:rem/h
||||| |||||
||||| |||||    ||||| +++++----- Used dose rate factor Neutron (0..15)
||||| |||||
||||| |||||    |||\----- 1: Accumulated counts
||||| |||||    |||\----- Training mode, 0:No, 1:Yes
```

Bit number			
0	Net value	0:No	1:Yes
1	Operation mode	0: Ratemeter	1:Scaler
2	Not used		
3			
4			
5			
6	Background measurement	0:No	1:Yes
7	dose rate unit gamma: 5: Sv/h, 6: R/h, 7rem/h		
8			
9			
10			
12	Measuring unit neutron 0: cps, 1: cpm, 5: Sv/h, 6: R/h, 7rem/h		
13			
14			
15			
16..19	Not used		
20	Accumulated counts	0:No	1:Yes
21	Training mode	0:No	1:Yes
22..31	Not used		

5.18.6 Nuclide table

nRNumber Reading nuclide data.
Number: number of nuclide
 Response: Nuclide data (see below)

nWNumberString Write nuclide data
Number: number of nuclide
String: nuclide data (see below)

Nuclide data: every string contains nuclide data

For example:

00 Sr-90 471

```

-----
|      | |
|      | \----- Factor in (nSv/h)/cps
|      \----- Name of nuclide
\----- Nuclide number
  
```

5.18.7 Event log

6656 520549251

```

-----
|      |
|      \----- Date and time as a decimal value (see 3.2)
\----- Event log as a decimal value
  
```

```

----- 0 0 0 0 0 0 0 0 = Alarmflags
      | | | | | | | |
      | | | | | | | \- Alarmschwelle 1 Gamma
      | | | | | | | \--- Alarmschwelle 2 Gamma
      | | | | +-+----- 01: Niederenergetisch
      | | | |      10: Hochenergetisch
      | | | |      11: Mittelenergetisch
      | | | |      00: NORM
      | | | \ ----- Alarmschwelle 1 Neutron
      | | \ ----- Alarmschwelle 2 Neutron
  
```

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error

5	EEPROM checksum error
6	Not used
7	Not used
8	Not used
9	
10	LCD-LED 0: off 1:on
11	Sound 0: off 1:on
12	LED 0: off 1:on
13	Vibration alarm 0: off 1:on
14	NBR on if device started 0: off 1:on
15	Single pulse 0: off 1:on
16	Dose rate or count rate gamma > alarm threshold 1
17	Dose rate or count rate gamma > alarm threshold 2
18	NBR-Alarm
19	00: NORM 10: High energy 01: Low energy 11: Mid energy
20	Dose rate or count rate neutron > alarm threshold 1
21	Dose rate or count rate neutron > alarm threshold 2
22	Dose alarm gamma
23	Dose alarm neutron
24	Training mode 0: off 1:on
25	Dose cleared 0: no 1:yes
26	Power off
27	Power on
28	Not used
29	Not used
30	Not used
31	Not used

5.18.8 Automatic sending

<STX>9 560 1 544 00 GN AB<ETX><CR><LF>

							\-----	BCC in HEX
							\-----	Information RadEye type
							\-----	Status as Hex-value
							\-----	Dimension neutron, see 5.18.3.5
							\-----	Rate neutron in 0.01cps, µR/h or 0.01µSv/h
							\-----	Dimension gamma, see 5.18.3.5
							\-----	Rate gamma in cps, µR/h or 0.01µSv/h

Formation of a BCC (block check character):

Modulo 256 sum of <STX> up to the last character before the BCC (including), coded as hexadecimal ASCII-number (e.g. 1F).

Status

Bit number	
0	Not used
1	1: Overload
2	1: Alarm gamma
3	1: Alarm neutron
4	1: Alarm dose gamma
5	1: Alarm fast neutron
6	1: NBR alarm
7	1: Safety alarm

Type

Type	RadEye
GN	RadEye GN
GN+	RadEye GN+

5.18.9 Status information

F Reading status information
 Response Number with status information, see 5.18.5.2

5.18.10 Miscellaneous

ARC Read single pulse divider
 Response: number for single pulse divider

AWCNumber Write single pulse divider
 Number: value from 1 to 65535

5.19 RadEye PRD-CD

Used firmware version: 3.10

$$C = \text{ADF}(\text{Rate1} - \text{Rate2}) - \text{ADF}(\text{Rate2}) * f$$

$$R = \text{ADF}(\text{Max}(\text{Rate1}), a * \text{Rate4}) - \text{Min}(\text{Rate1} - \text{Rate2}; 1.5 * C(\text{Cal}))$$

$$\text{Dimension less Value} = \frac{C}{C(\text{cal})} * \text{Ref}$$

C(cal): Cal measurement on left block

5.19.1 Limit values

AR0	Reading computed sigma alarm threshold 1. Response: number in cps.
AR1	Reading alarm threshold 1 and 2 for contraband Response: number in 0.01 cps.
AR2	Reading alarm threshold 1 and 2 for Rate 2 Response: number in 0.01 cps.
AR3	Reading the threshold 1 for gross count rate alarm. Response: number in 0.01 cps.
AR4	Reading the threshold 2 for gross count rate alarm. Response: number in 0.01 cps.
AR5	Reading the threshold 1 for dose rate alarm. Response: number in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ unit. Depending on used measuring unit e.g. 1 2 3 means 1 2 3 $\mu\text{R/h}$ resp. 1.23 $\mu\text{Sv/h}$
AR6	Reading the threshold 2 for dose rate alarm. Response: see command AR5
AR7	Reading the threshold 1 for the dose alarm Response: number in μSv , $100\mu\text{R}$ or $100\mu\text{rem}$ units
AR8	Reading the threshold 2 for the dose alarm see command AR7
AR9	Reading sigma value Response: Number from 2..9

ARM	Reading minimum count rate for sigma alarm Response: Number in cps
ARN	Reading: <ul style="list-style-type: none"> - NBR alarm threshold low energy Response: Number in 0.01% units - Background preset count Response: Number in counts - Background preset time Response: Number in seconds
ARP	Reading scaler parameter: <ul style="list-style-type: none"> - Preset count. Response: Number in counts - Preset time. Response: Number in seconds - Scaler wait. time Response: Number in seconds
ARB	Reading scaler background values: <ul style="list-style-type: none"> - Background value Response: Number in 0.01 cps units - Used time for background value Response: Number in seconds
ARK	Reading contraband parameter <ul style="list-style-type: none"> - a with Factor 100 (65 means 0,65) - f with Factor 1000 (1300 means 1.3) - C(cal) with Factor 1000 (150000 means 150) - Calibration Precision in 0.1%-units (16 means 1.6%) - Calibration Reference <i>Ref</i> - Timeout Contrabandmode in s. From 0s(=off) to 64800s (=18h)
ARH	NBR alarm threshold high energy Response: Number in 0.01% units
SR3	Reading minimum count rate for NBR alarm Response: Number in cps
AW1 <i>Number Number</i>	Set the alarm threshold 1 and 2 for contraband. <i>Number</i> = Number in 0.01 cps.
AW2 <i>Number Number</i>	Set the alarm threshold 1 and 2 for Rate 2. <i>Number</i> = Number in 0.01 cps.
AW3 <i>Number</i>	Setting the threshold 1 for gross count rate alarm. <i>Number</i> = Number in 0.01 cps.
AW4 <i>Number</i>	Setting the threshold 2 for gross count rate alarm. <i>Number</i> = see command AW3.
AW5 <i>Number</i>	Setting the threshold 1 for dose rate alarm. <i>Number</i> = in $\mu\text{R/h}$, $\mu\text{rem/h}$ or $0.01\mu\text{Sv/h}$ units.

	e. g. AW3123 means 123µR/h.
AW6Number	Setting the threshold 2 for dose rate alarm. see command AW3
AW7Number	Setting the threshold 1 for dose alarm. Number = in 100µR, 100µrem or 1µSv units. e. g. AW7123 means 12.3mR.
AW8Number	Setting the threshold 2 for dose alarm. Number = see command AW7
AW9Number	Setting the sigma value. Number = form 2 to 9
AWMNumber	Setting minimum count rate for sigma alarm Number= value from 0 to 255
AWNNumber Number Number	Setting NBR alarm threshold level low energy Number=value in 0.01% units. From 1.00% to 2.55% Setting background preset count Number=counts. From 0 to 9999 Setting background preset time Number=seconds. From 0 to 9999s
AWHNumber	Setting NBR alarm threshold level high energy Number=value in 0.01% units. From 0% to 1.00%
AWBNumber Number	Setting background value Number=0.01 cps units. From 0 to 10000 (100cps) Setting used time for background value Number=seconds. From 0 to 9999s
AWPNumber Number	Setting scaler preset count Number=counts. From 0 to 9999 Setting scaler preset time Number=seconds. From 0 to 9999s
AWKNumber Number	Writing contraband parameter <ul style="list-style-type: none"> - a with Factor 100 (65 means 0,65) - f with Factor 1000 (1300 means 1.3) - C(cal) with Factor 1000 (150000 means 150) - Calibration Precision in 0.1%-units (16 means 1.6%) - Calibration Reference <i>Ref</i> - Timeout Contrabandmode in s. From 0s(=off) to 64800s (=18h)

SW3 Number Setting minimum count rate for NBR alarm
Number= value from 0 to 255

5.19.2 Measurement values

Z Read count rates with dead time correction
Response:
 - Counter 1 in cps
 - Counter 2 in cps
 - Counter 3 in cps
 - Counter 4 in cps
 - PMT current index

z Read filtered count rate 1
Response: Number in cps

5.19.3 Configuration flags

5.19.3.1 Configuration flags 1 with *kR* / *kW*

Bit number	
0	not used, write “0”
1	
2	
3	Alarming Sound 0: off 1: on
4	Alarming LED 0: off 1: on
5	Alarming Vibration 0: off 1: on
6	NBR 0: off 1: on
7	Single Pulse 0: off 1: on

5.19.3.2 Configuration flags 2 with *fR* / *fW*

Bit number	
0	not used, write “0”
1	
2	not used, write “0”
3	Alarm threshold read-only 0: off 1: on
4	Flag for overload (readonly)
5	Temp.display 0: off 1: on
6	Temperature unit 0: °C 1: °F
7	In cps-mode display of dose rate 0: no 1: yes

5.19.3.3 Configuration flags 3 with KR / KW

Bit number	
0	Key tone 0:off 1:on
1	Autosend 0:off 1:on
2	Finder 0:off 1:on
3	Mean-Max 0:off 1:on
4	not used, write "0"
5	display of dose 0:off 1:on
6	Use click divider for single pulse 0:off 1:on
7	not used, write "0"

5.19.3.4 Configuration flags 4 with jR / jW

Bit number	
0	Safety-Alarm 0:off 1:on
1	Battery type 0:Alkaline 1:NiMh
2	Display rotation 0:No 1:yes
3	Scaler mode 0: preset count 1: preset time
4	Scaler "Auto restart" 0:off 1:on
5	Net Scaler 0:off 1:on
6	Net Ratemeter 0:off 1:on
7	not used, write "0"
8	not used, write "0"
9	not used, write "0"
10	Used Ba-133 activity: 0: 10 μ Ci 1: 5 μ Ci
11	Alarm LCD-LED 0:off 1:on
12	not used, write "0"
13	not used, write "0"
14	not used, write "0"
15	not used, write "0"

5.19.3.5 *Measuring unit with uR /uW*

uR Read measuring unit and operation mode
 Response: Hex-value. See below
uRHex Write measuring unit and operation mode
 See below

Bit number	
0	0x00 Display unit cps
1	0x05 Display unit Sv/h
2	0x06 Display unit R/h
3	0x07 Display unit rem/h 0x08 Display unit Bq
4	0: Ratemeter, 1: Scaler
5	0x01 Display unit Sv/h
6	0x02 Display unit R/h 0x03 Display unit rem/h
7	1: Contraband mode

5.19.3.6 *Menu configuration*

mR Read menu configuration for

- Main menu
- Submenu “Settings”
- Submenu “Alarm indication”
- Submenu “Set alarm”

Response: Hex-values. See below

mWHex Hex Hex Write menu configuration for

- Main menu
- Submenu “Settings”
- Submenu “Set alarm”

See below

5.19.3.6.1 Main menu

Bit number			
0	Switch off	0:hidden	1:visible
1	Background	0: hidden	1:visible
2	Not used, write “0”		
3	Backlight	0: hidden	1:visible
4	Measuring unit	0: hidden	1:visible
5	Operation mode	0: hidden	1:visible
6	Scaler parameter	0: hidden	1:visible
7	Not used, write “0”		
8	Set alarm	0: hidden	1:visible
9	Not used, write “0”		
10	Not used, write “0”		
11	Not used, write “0”		
12	Not used, write “0”		
13	Settings	0: hidden	1:visible
14	Alarm indication	0: hidden	1:visible
15	Show alarm	0: hidden	1:visible
16	Text info	0: hidden	1:visible
17	Bluetooth	0: hidden	1:visible
18..31	not used, write “0”		

5.19.3.6.2 Submenu “Settings”

Bit number			
0	Batt. type	0:hidden	1:visible
1	Autosend	0: hidden	1:visible
2	Single Pulse	0: hidden	1:visible
3	Finder	0: hidden	1:visible
4	Set Date/Time	0: hidden	1:visible
5	Lu-Test	0: hidden	1:visible
6	Language	0: hidden	1:visible
7	Spill-over	0: hidden	1:visible
8	Rated Alarm	0: hidden	1:visible
9	Contrast	0: hidden	1:visible
10..15	not used, write “0”		

5.19.3.6.3 Submenu “Alarm indication”

Bit number			
0	Sound	0:hidden	1:visible
1	LED	0: hidden	1:visible
2	Vibrator	0: hidden	1:visible
3	LCD-LED	0: hidden	1:visible
4...8	not used, write “0”		

5.19.3.6.4 Submenu “Set alarm”

Bit number			
0	Count Rate	0:hidden	1:visible
1	Not used, write “0”		
2	Contraband	0: hidden	1:visible
3	Dose rate	0: hidden	1:visible
4	Dose	0: hidden	1:visible
5	Alarm NBR	0: hidden	1:visible
6...8	not used, write “0”		

5.19.4 Spill-over

ARK Reading spill-over
Response: value 0...255 (0.1 % units, default 80 means 8.0%)

AWKNumber Set spill-over
Number from 0...255

5.19.5 High voltage

HR	Reading high voltage bit value Response: Bit value 0...255.
hR	Reading high voltage correction bit date and time of last successful Luthetium check Response: Number 0...255 with offset 128 and date and time as YYMMDDhhmm
hWNumber	Setting of high voltage correction bit Number from 0...255

5.19.6 Dead time correction

x	Read dead time Response: dead time in ns for <ul style="list-style-type: none">- Counter 1 (Rate 1)- Counter 2 (Rate 2)- Counter 3 (Rate 3)- Counter 4 (Rate 4)
---	---

5.19.7.1 History readout

1536 716612088 1239 1600 30 5 120 23 4

End of History:
End

decimal value converted in HEX:

Bit number			
0	Ratemeter/Scaler net value	0: no	1: yes
1	Operation mode Scaler	0: no	1: yes
2	Not used, read as '0'		
3	Not used, read as '0'		
4	Background measurement	0: no	1: yes
5	Preset time	0: no	1: yes
6	Not used, read as '0'		
7	Contraband mode	0: no	1: yes
8..11	5:Sv/h, 6:R/h, 7:rem/h, 8: Contamination (Bq)		

12..15	Not used, read as '0'
--------	-----------------------

5.19.8 Event log

6656 520549251

	\-----	Date and time as a decimal value (see 3.2)
\-----		Event log as a decimal value

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	0x02: Display count rate 0x04: Display dose rate
9	
10	
11	Sound 0: off 1:on
12	LED 0: off 1:on
13	Vibration alarm 0: off 1:on
14	1: Dose cleared
15	1: Alarm threshold changed
16	1: count rate, dose rate, contraband or Rate 2-alarm
17	1: Dose alarm
18	1: Safety alarm
19	Not used
20	Value > alarm threshold 1 (count rate, dose rate, contraband or Rate 2)
21	Value > alarm threshold 2 (depending on count rate, dose rate or level display)
22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24	Low energy alarm
25	High energy alarm
26	Power off
27	Power on
28	NBR-alarm
29	Not used
30	Not used
31	Not used

5.19.9 Calibration log

520549251 175590 0 0001

```

-----
|           |           |           |           |
|           |           |           | \-- Status
|           |           | \----- Rate "R" with factor 1000
|           | \----- Rate "C" with factor 1000
| \----- Date and time as a decimal value (see 3.2)
\-----

```

Status

Bit number	
0	Calibration value saved
1..15	Not used

CI Initialize calibration log readout. Return: Number of saved entries

CC Clear calibration log

C+ Read next entry

5.19.10 Automatic sending

```

<STX>7 2 0 0 14 PRD-CD 123 AB<ETX><CR><LF>
| | | | | | | |
| | | | | | | \- BCC in HEX
| | | | | | \----- Dose in µR/h
| | | | | \----- information RadEye PRD-CD
| | | | \----- Status as Hex-value
| | | \----- „0" not used
| | \----- „0" not used
| \----- „2" not used
\----- Dose rate in µR/h

```

Formation of a BCC (block check character):

Modulo 256 sum of <STX> up to the last character before the BCC (including), coded as hexadecimal ASCII-number (e.g. 1F).

Status

Bit number	
0	Not used
1	1: Overload
2	1: Count rate, dose rate or level-alarm
3	1: Dose alarm
4	1: NBR alarm
5	1: Battery voltage low
6	not used
7	not used

5.19.11 Status information

F Reading status information
 Response Number with status information

Bit number	
0	HV-Error
1	Detector error
2	Low Battery voltage
3	Not used
4	Watchdog error
5	EEPROM checksum error
6	Not used
7	Not used
8	0x02: Display count rate 0x04: Display dose rate
9	
10	
11	not used,
12	Alarm threshold read-only 0: off 1:on
13	Flag for overload
14	Temperature display 0: off 1:on
15	In cps-mode display of dose rate 0: no 1: yes
16	1: Count rate, dose rate or level-alarm
17	1: Dose alarm
18	1: Safety alarm
19	Not used
20	Value > alarm threshold 1 (depending on count rate, dose rate or level display)
21	Value > alarm threshold 2 (depending on count rate, dose rate or level display)
22	Dose > alarm threshold 1
23	Dose > alarm threshold 2
24	Low energy alarm
25	High energy alarm
26	Not used
27	Not used
28	Not used
29	Not used
30	Not used
31	Fixed to "0"

Thermo Fisher Scientific Messtechnik GmbH
Frauenauracher Straße 96
D-91056 Erlangen, Germany

+49 9131 998-0
+49 9131 998 475 fax

Post address:
P.O. Box 16 60
D-91051 Erlangen, Germany

www.thermofisher.com
info.rmp.erlangen@thermofisher.com